




DUOSWIFT®

The DuoSwift combination squeegee brush features a tapered nylon brush on the leading end and a series of squeegees on the trailing end to effectively remove debris from flexible endoscopes while adhering to industry guidelines.

ABOUT:

- Bristled brush to help loosen and remove bioburden with friction
- Squeegee to pull liquid and remaining residue through the endoscope channel
- Unique design to help maximize cleaning touches within compatible endoscope channels



REF	DESCRIPTION	BRISTLE \varnothing (mm)	SQUEEGEE \varnothing (mm)	QTY
BX00711618	DuoSwift combination squeegee brush	5-7	2,5-5,2	50 pcs
BX00711619	DuoSwift combination squeegee brush with valve/control head cleaning brush	5-7 / 10	2,5-5,2 (squeegee)	50 pcs

