



BAND CLAMPS

Instrument jaws need to be in an open position, and kept that way, during cleaning and sterilization. Key Surgical® Band Clamps are a helpful tool for this, when an instrument stringer is not an option.

“Ratcheted instruments should be unlatched. Racks, pins, stringers, or other specifically designed devices can be used to hold the instruments in the unlatched position.”

-ANSI/AAMI ST79:2017, Section 8.2

ABOUT:

- Helps keep joints and jaws of various instruments open for more effective/efficient cleaning
- Silicone-covered stainless steel
- Compatible with steam sterilization
- Not made with natural rubber latex
- Available in 5 sizes
- Color-coded
- 1/pkg

PRODUCT NO.	WIDTH	COLOR
IS-P0101	1.5"	white
IS-P0105	2"	yellow
IS-P0102	2.25"	blue
IS-P0104	2.75"	red
IS-P0103	3.25"	green

