



BOWIE-DICK TEST PACKS

Designed for daily use to evaluate the efficacy of air removal and steam penetration in pre-vacuum steam sterilizers. Certified to meet AAMI standards. Available in either disposable or reusable type test packs.

THE LANTOR CUBE®

The industry's first and only FDA approved reusable Bowie-Dick Test Pack. The metal clamp (sold separately) holds the test pack in place.



- Reduces waste
- Great for inventory control
- Comes with 25 consecutively numbered test sheets and is reusable for 25 cycles

Number	Description	Quantity
CS-BD001	Lantor Cube Metal Clamp	1
CS-BD025	Lantor Cube Bowie-Dick Test Pack	25 tests/box



AIRVIEW®

Disposable Bowie-Dick Test Pack.

- Instructions for use printed on the box for ease of use
- Easy to interpret indicator test sheet for clear pass/fail results
- External process indicator to show pack has been through the sterilizer



Number	Description	Quantity
CS-BD030	Airview Bowie-Dick Test Pack	30 packs/case



FOR MORE INFORMATION CALL 800.541.7995 OR VISIT KEYSURGICAL.COM

