


8101 WALLACE ROAD
EDEN PRAIRIE, MN 55344 USA


TEL: 800.541.7995 952.914.9789
FAX: 952.914.9866

info@keysurgical.com keysurgical.com



STERILE PROCESSING
EDITION 31







FOUR EASY WAYS TO PLACE AN ORDER WITH KEY SURGICAL

It's easy to place an order with us – we're available to take your order by telephone, fax, e-mail or EDI. If you would like a standard form to fax an order to us, simply visit keysurgical.com for a form that is easy to download.

- 1 EDI via GHX
- 2 E-MAIL order@keysurgical.com
- 3 TELEPHONE 800.541.7995
- 4 FAX 952.914.9866

TERMS OF PURCHASE

No minimum purchase required (customized products excluded). All shipments are FOB Eden Prairie, Minnesota, USA. Our No-Risk Return Policy, allows you to try our products for 30 days and return them if you are not completely satisfied. For full terms of purchase and our No-Risk Return Policy, please see Terms of Purchase on our website. All product returns are subject to a restocking fee. Approved credit terms are Net 30 from date of invoice.

We also accept Visa, Mastercard, and American Express. Although our prices are subject to change without notice, we take pride in offering competitive pricing.

- 1 TELEPHONE 800.541.7995
- 2 E-MAIL info@keysurgical.com
- 3 WEBSITE keysurgical.com
- 4 FAX 952.914.9866
- 5 MAIL 8101 Wallace Road
Eden Prairie, MN 55344

FOR MORE INFORMATION

If you need more information, feel free to give us a call. Too busy to talk? You can always reach us by e-mail, fax or even snail mail if it suits you. Our knowledgeable customer service department is waiting to answer your calls on product questions. We have the answer or we'll find it for you. Or, if you prefer to get your information digitally, our website is designed to be a complete reference tool for you. You can always e-mail us with your questions; we respond to all e-mail inquiries within 24 hours. Give us a try – we look forward to hearing from you.

THE KEY SURGICAL STORY

Key Surgical was founded in 1988 by three entrepreneurial women with years of experience in both the nursing profession and the medical supply industry. Motivated by the growing need for a single-source supplier that would provide excellent customer service, the partners began building an extensive product line based on market need. The goal was to serve the varied needs of hospitals, surgical ambulatory, and outpatient centers as completely and as quickly as possible.

Over 30 years later, the same customer-driven principals and patient-focused values upon which Key Surgical was founded continue to guide everything we do. With over 3,000 of the highest quality products available to customers worldwide, Key Surgical is poised to continue to be a compassionate industry leader for years to come.

OUR MISSION

Our mission at Key Surgical is three-fold. First, to hire and retain good people who understand that nothing happens until a customer decides to buy something from us. Second, to grow by providing the quality products, value and service that our customers want, when they want them. And finally, to be good citizens by providing a rewarding and fun place for employees to work and sharing the benefits of our growth with the community in which we operate.

A NOTE ABOUT QUALITY

There are no shades of gray in the medical marketplace when it comes to quality. Key Surgical is committed to manufacturing and distributing the highest quality medical products, and to maintaining the highest level of compliance to FDA and ISO regulations. Nobody can claim to be perfect, but we strive to continuously improve our processes and systems so that we will always be closer to perfect tomorrow than we are today.

LETTER TO OUR CUSTOMERS

Thank you for being our customer and for allowing us to be your trusted partner for products and solutions for your sterile processing department. Customers vote with their purchasing decisions, and we know that to earn your vote, we need to continue to set the standard for quality, customer service, and value.

Our customer base of hospitals, GPOs, surgery centers, and distributors continues to expand to over 10,000 locations in the United States and in 70 other countries around the world. We're committed to continuously making substantial investments in our products, quality system, regulatory and compliance efforts, and operations with the goal of being the easiest company to do business with. Our goal, as it was over 30 years ago when Key Surgical was founded, is to supply you with the products you need, deliver them when you want them, and support you with personalized customer service. Although seemingly simple, focusing on delivering an exceptional experience with us every step of the way is, and will continue to be, the foundation of why we do what we do.

Our mission has three simple parts and we hold ourselves to be true to it every day. First, we will hire and retain good people who understand that our business exists because customers decide to buy something from us. Second, we will grow by providing the quality products, value, and service that our customers want, when they want them. Finally, we will be good citizens by providing a rewarding and fun place for employees to work as well as sharing the benefits of our growth with the community in which we operate. Combine our mission with our steadfast focus on improving patient safety and positive patient outcomes through products and education and you have what makes up Key Surgical.

Our entire team works hard every day to be recognized by you as a company of dedication, honesty, and integrity. In striving toward that recognition, we promise to keep our feet planted firmly on the ground and to make sure we brighten both your day and ours in the process. If we can ever be of assistance to you please call or email us directly. Businesses grow because they do the right things for their customers and their employees. We intend to keep doing those things and to keep growing, and we truly appreciate your business and thank you for continued support.

Sincerely,



Scot Milchman
Chief Executive Officer



Brian O'Connell
President and Chief Operating Officer








TABLE OF CONTENTS

-  000 – Instrument Cleaning
-  058 – Instrument Protection
-  078 – Tags
-  106 – Identification Tape
-  124 – Personal Protective Equipment
-  136 – Packaging and Preparation
-  166 – Wire and Pin Management
-  186 – Instrument Tracking
-  192 – Index



PLEASE NOTE: Unless specified, all products are not made with natural rubber latex.
 PROP 65 Warning Disclaimer: Please review item information at [keysurgical.com](https://www.keysurgical.com) for the most up-to-date Proposition 65 warning information.

INSTRUMENT CLEANING

The Key Surgical Instrument Cleaning line is designed to provide you with the tools needed to clean your surgical instruments. In addition to our wide variety of cleaning brushes, we also offer enzymatic solutions to help improve the decontamination process.

With bristles available in nylon, brass, polypropylene or stainless steel, we also provide an incredible variety of lengths and diameters to fit your needs.

For organization of your cleaning brushes, our brush organizer is perfect for hanging in the department. Our stainless steel organizer is easy to keep clean and helps your team choose the right brush in a more efficient manner.

With our vast assortment of brushes - from channel cleaning to endoscope cleaning to sterilizer cleaning - we can help you match the right brush to the right cleaning task at hand.

Dr. Weigert neodisher® PreStop

Dr. Weigert neodisher® MultiZym

EfferZyme™

Channel Cleaning Brushes

Pipe Cleaners

Toothbrush-Style Cleaning Brushes

General Instrument Cleaning Brushes

Endoscope Cleaning Brushes

Brush Assortments

Sterilizer Cleaning Brushes

Brush Organizer

French Catheter Scale

Please note: Unless specified, all products are not made with natural rubber latex.



Instrument Cleaning – PreStop, MultiZym, and EfferZyme



DR. WEIGERT NEODISHER® PRESTOP

Dr. Weigert PreStop is intended for use as a corrosion inhibitor for spraying onto surgical instruments in operating rooms directly after use. Promotes the practice of point-of-use cleaning in the O.R. Protects stainless steel instruments and keeps them moist until further reprocessing can be done. Non-aerosol spray.

Number	Details	Quantity
A	DI-442647-12	25.4 fluid oz 12/case

DR. WEIGERT NEODISHER® MULTIZYM

Dr. Weigert MultiZym is a multi-enzymatic detergent with protease, lipase, and amylase. It is intended for use for the manual cleaning of thermostable and thermolabile instruments, including MIS and micro instruments, dental instruments, flexible endoscopes, anesthetic equipment, containers, and other medical utensils in immersion and ultrasonic baths.

Number	Details	Quantity
B	DI-442735	1.32 gal 1/pkg


EFFERZYME™ - EFFERVESCENT CLEANING TABLETS

EfferZyme cleaning tablets are made from biodegradable, non-chlorine ingredients that contain fast acting protease enzymes. These neutral pH tablets enhance the removal of blood and other soils that may be resistant to chemical detergents, while still being delicate enough on instruments. May be used as an ultrasonic cleaning solution or an instrument presoak.

Number	Quantity
C	IC-2152 52 tablets/box



Instrument Cleaning – Channel Cleaning Brushes




CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. FR indicates equivalent French catheter size.

Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A BR-06-039 BR-06-039-50	.039"	1.00mm	3	.25"	6"	2 50
B BR-08-039 BR-08-039-50 BR-18-039 BR-18-039-50	.039"	1.00mm	3	.50"	8" 18"	2 50 2 50
C BR-06-062 BR-06-062-50 BR-18-062 BR-18-062-50	.062"	1.58mm	5	.50"	6" 18"	2 50 2 50
D BR-08-070 BR-08-070-50 BR-18-070 BR-18-070-50 BR-24-070 BR-24-070-50	.070"	1.78mm	5	.50"	8" 18" 24"	2 50 2 50 2 50
E BR-06-079 BR-06-079-50	.079"	2.00mm	6	.50"	6"	2 50
F BR-08-079 BR-08-079-50 BR-12-079 BR-12-079-50	.079"	2.00mm	6	.75"	8" 12"	2 50 2 50
G BR-24-080 BR-24-080-50	.080"	2.03mm	6	1"	24"	2 50
H BR-08-087 BR-08-087-50 BR-18-087 BR-18-087-50	.087"	2.20mm	7	.50"	8" 18"	2 50 2 50
I BR-08-092 BR-08-092-50	.092"	2.34mm	7	.50"	8"	2 50
J BR-24-098 BR-24-098-50	.098"	2.50mm	8	.50"	24"	2 50
K BR-08-105 BR-08-105-50 BR-18-105 BR-18-105-50	.105"	2.67mm	8	1"	8" 18"	2 50 2 50
L BR-08-118 BR-08-118-50	.118"	3.00mm	9	1.50"	8"	2 50
M BR-12-118 BR-12-118-50	.118"	3.00mm	9	2.75"	12"	2 50

BRUSH TIP DETAILS SHOWN ACTUAL SIZE






CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. FR indicates equivalent French catheter size.

	Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A	BR-18-125	.125"	3.18mm	10	2"	18"	2
	BR-18-125-50						50
B	BR-08-131	.131"	3.33mm	10	.75"	8"	2
	BR-08-131-50						50
C	BR-08-144	.144"	3.67mm	11	.75"	8"	2
	BR-08-144-50						50
D	BR-08-158	.158"	4.00mm	12	.75"	8"	2
	BR-08-158-50						50
	BR-12-158	.158"	4.00mm	12	.75"	12"	2
	BR-12-158-50						50
E	BR-18-158	.158"	4.00mm	12	.75"	18"	2
	BR-18-158-50						50
	BR-12-197	.197"	5.00mm	15	3"	12"	2
	BR-12-197-50						50
	BR-16-197	.197"	5.00mm	15	3"	16"	2
	BR-16-197-50						50
F	BR-18-197	.197"	5.00mm	15	3"	18"	2
	BR-18-197-50						50
	BR-24-197	.197"	5.00mm	15	3"	24"	2
	BR-24-197-50						50
G	BR-24-200	.200"	5.08mm	15	2"	24"	2
	BR-24-200-50						50
G	BR-12-236	.236"	6.00mm	18	2"	12"	2
	BR-12-236-50						50
	BR-16-236	.236"	6.00mm	18	2"	16"	2
	BR-16-236-50						50
	BR-24-236	.236"	6.00mm	18	2"	24"	2
	BR-24-236-50						50

BRUSH TIP DETAILS SHOWN ACTUAL SIZE



CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. FR indicates equivalent French catheter size.

Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A BR-06-250	.250"	6.35mm	19	2"	6"	2
BR-06-250-50						50
BR-08-250	.250"	6.35mm	19	2"	8"	2
BR-08-250-50						50
BR-18-250	.250"	6.35mm	19	2"	18"	2
BR-18-250-50						50
B BR-16-276	.276"	7.00mm	21	3"	16"	2
BR-16-276-50						50
BR-24-276	.276"	7.00mm	21	3"	24"	2
BR-24-276-50						50
C BR-24-312	.312"	8.00mm	24	3"	24"	2
BR-24-312-50						50
D BR-06-375	.375"	9.52mm	29	1.5"	6"	2
BR-06-375-50						50
E BR-24-375	.375"	9.52mm	29	3"	24"	2
BR-24-375-50						50
F BR-16-393	.393"	10.00mm	30	3"	16"	2
BR-16-393-50						50
BR-18-393	.393"	10.00mm	30	3"	18"	2
BR-18-393-50						50
BR-24-393	.393"	10.00mm	30	3"	24"	2
BR-24-393-50						50
G BR-24-425	.425"	10.80mm	32	3"	24"	2
BR-24-425-50						50

BRUSH TIP DETAILS SHOWN ACTUAL SIZE



CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. FR indicates equivalent French catheter size.

Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A BR-08-440	.440"	11.18mm	34	3"	8"	2
BR-08-440-50						50
BR-24-440	.440"	11.18mm	34	3"	24"	2
BR-24-440-50						50
B BR-16-443	.443"	11.25mm	34	3"	16"	2
BR-16-443-50						50
BR-24-443	.443"	11.25mm	34	3"	24"	2
BR-24-443-50						50
C BR-16-472	.472"	12.00mm	N/A	3"	16"	2
BR-16-472-50						50
D BR-06-500	.500"	12.70mm	N/A	2"	6"	2
BR-06-500-50						50
E BR-16-591	.591"	15.00mm	N/A	2"	16"	2
BR-16-591-50						50
BR-18-591	.591"	15.00mm	N/A	2"	18"	2
BR-18-591-50						50

A LITTLE BIT ABOUT...

Brian O'Connell

Job Title: **President and Chief Operating Officer**

Do you collect anything? If so, what? And why? **I collect a logo golf ball from every different course that I have the opportunity to play. I've been doing this for over 25 years.** What is your favorite golf course in the U.S.? What about international? **In the US, it's Pacific Dunes on the Oregon Coast. Internationally, it's a toss-up between two Scottish courses; The Old Course in St. Andrews or Kingsbarns.** What cheers you up? **Talking with my kids and really good comedy. Sometimes they're one in the same.** What are your top 3 favorite movies? **Glengarry Glen Ross, Caddyshack, and The Empire Strikes Back.** What is your favorite meal/food? **A good pizza is always at the top of the list, but sushi can give it some strong competition.** What's your favorite recreational activity? **Golf and Boating; can't pick a winner between the two.** Who inspires you? **Great scientific thinkers such as Einstein and Hawking.** What is something you've always wanted to try but have been too afraid to? **Skydiving.**





CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end.

Number	Diameter inches	Diameter mm	Bristle Area Length	Length	Quantity/pkg
A BR-08-630 BR-08-630-50	.630"	16.00mm	2"	8"	2 50
B BR-16-787 BR-16-787-50	.787"	20.00mm	2"	16"	2 50
C BR-18-880 BR-18-880-50	.880"	22.35mm	3.50"	18"	2 50
D BR-24-880 BR-24-880-50	.880"	22.35mm	4.25"	24"	2 50

FUN FACT:

The three women who started Key Surgical in 1988 had years of experience in both the nursing profession and the medical supply industry. They were motivated by the growing need for a single-source supplier that would provide excellent products and equally excellent customer service. One of the founders is Gail Milchman, the mother of our CEO, Scot Milchman.

BRUSH TIP DETAILS SHOWN ACTUAL SIZE






CHANNEL CLEANING BRUSHES


Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end (pull through style excluded). FR indicates equivalent French catheter size.

Number	Details	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A BR-3044	nylon	.020"	.50mm	2	.25"	3"	10
BR-3045	nylon	.020"	.50mm	2	.25"	8"	10
B BR-3047	nylon	.035"	.90mm	3	.5"	4"	10
C BR-4120	nylon, tapered	.250" to .600"	6.35mm to 15.24mm	N/A	2"	18"	2
BR-4120-50							50
D BR-4140	nylon, centered, pull through	.500"	12.70mm	N/A	3"	24"	6



(PULL THROUGH)

BRUSH TIP DETAILS SHOWN ACTUAL SIZE



CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. FR indicates equivalent French catheter size.

Number	Description	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/ pkg
A	BR-8110 curved bristle tip	.394"	10mm	29	1.75"	4.6"	3
	BR-8110-50						50
B	BR-8115 curved bristle tip	.551"	13.99mm	NA	3.625"	8"	3
	BR-8115-50						50
C	BR-8210 tapered fan bristle tip	.187" to .500"	4.75mm to 12.70mm	NA	3"	7"	3
	BR-8210-50						50
D	BR-8310 straight bristle tip	.500"	12.70mm	NA	2.5"	12"	3
	BR-8310-50						50
E	BR-8315 straight bristle tip	.500"	12.70mm	NA	3"	12"	3
	BR-8315-50						50
F	BR-8410 fan bristle tip	.500"	12.70mm	NA	3"	6"	3
	BR-8410-50						50

DID YOU KNOW?

The Key Surgical Community Giving Committee seeks out volunteer opportunities for all Key Surgical employees. We partner with local organizations that provide services and assistance to low-income families, people with disabilities, children, the homeless, and those suffering from terminal illness.


BRUSH TIP DETAILS SHOWN ACTUAL SIZE




BONE REAMER BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a rigid stainless steel handle with loop end. FR indicates equivalent French catheter size.

Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A BR-31-118	.118"	3.00mm	9	3"	31"	3
B BR-31-197	.197"	5.00mm	15	3"	31"	3
C BR-31-276	.276"	7.00mm	21	3"	31"	3



BRUSH TIP DETAILS SHOWN ACTUAL SIZE




LUMEN GUARD™ CLEANING BRUSHES

Channel cleaning brushes featuring rigid, antimicrobial nylon bristles and a protective tip on a twisted stainless steel handle with loop end. Inseparable, protective tip created using patented Z-tip technology. FR indicates equivalent French catheter size.

Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A BZ-10-050 BZ-10-050-50	.050"	1.27mm	4	1"	10"	2 50
B BZ-12-062 BZ-12-062-50	.062"	1.58mm	5	1"	12"	2 50
C BZ-18-062 BZ-18-062-50	.062"	1.58mm	5	2"	18"	2 50
D BZ-18-079 BZ-18-079-50	.079"	2.00mm	6	2"	18"	2 50
E BZ-12-087 BZ-12-087-50	.087"	2.20mm	7	1"	12"	2 50
F BZ-12-098 BZ-12-098-50	.098"	2.50mm	8	1"	12"	2 50
G BZ-18-105 BZ-18-105-50	.105"	2.67mm	8	2"	18"	2 50
H BZ-08-125 BZ-08-125-50	.125"	3.18mm	10	.75"	8"	2 50
I BZ-12-125 BZ-12-125-50	.125"	3.18mm	10	1"	12"	2 50
J BZ-18-125 BZ-18-125-50 BZ-24-125 BZ-24-125-50	.125"	3.18mm	10	2"	18" 24"	2 50 2 50
K BZ-18-131 BZ-18-131-50	.131"	3.33mm	10	2"	18"	2 50
L BZ-18-144 BZ-18-144-50	.144"	3.37mm	11	2"	10"	2 50
M BZ-12-158 BZ-12-158-50	.158"	4.00mm	12	1"	12"	2 50
N BZ-18-158 BZ-18-158-50	.158"	4.00mm	12	2"	18"	2 50

BRUSH TIP DETAILS SHOWN ACTUAL SIZE




LUMEN GUARD™ CLEANING BRUSHES (CONTINUED)

Channel cleaning brushes featuring rigid, antimicrobial nylon bristles and a protective tip on a twisted stainless steel handle with loop end. Inseparable, protective tip created using patented Z-tip technology. FR indicates equivalent French catheter size.

Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A BZ-18-197 BZ-18-197-50	.197"	5.00mm	15	2"	18"	2 50
B BZ-18-236 BZ-18-236-50	.236"	6.00mm	18	2"	18"	2 50
C BZ-16-250 BZ-16-250-50 BZ-18-250 BZ-18-250-50	.250"	6.35mm	19	2"	16"	2 50 2 50
D BZ-18-312 BZ-18-312-50	.312"	8.00mm	24	2"	18"	2 50
E BZ-18-354 BZ-18-354-50	.354"	9.00mm	27	2"	18"	2 50
F BZ-18-393 BZ-18-393-50	.393"	10.00mm	30	2"	18"	2 50
G BZ-18-433 BZ-18-433-50	.433"	11.00mm	33	2"	18"	2 50
H BZ-18-472 BZ-18-472-50	.472"	12.00mm	N/A	2"	18"	2 50
I BZ-18-1000 BZ-18-1000-50	1"	25.40mm	N/A	2"	18"	2 50
J BZ-18-1543 BZ-18-1543-50	.157" to .433"	4mm to 11mm	N/A	2"	18"	2 50
K BZ-06-2550 BZ-06-2550-50	.250" to .500"	6.35mm to 12.7mm	N/A	.75"	6"	2 50

BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Channel Cleaning Brushes




FAN TIP BRUSHES

Channel cleaning brushes featuring a fan tip to assist with cleaning closed-end channels. Antimicrobial nylon bristles on a twisted stainless steel handle with loop end. FR indicates equivalent French catheter size.

Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A	.197"	5.00mm	15	2"	12"	2
					50	
	.197"	5.00mm	15	2"	18"	2
					50	
					24"	2
					50	
B	.250"	6.35mm	19	2"	12"	2
					50	
	.250"	6.35mm	19	3"	18"	2
					50	
					24"	2
					50	
C	.276"	7.00mm	21	2"	12"	2
					50	
	.276"	7.00mm	21	3"	18"	2
					50	
					24"	2
					50	
D	.375"	9.52mm	29	2"	12"	2
					50	
	.375"	9.52mm	29	3"	18"	2
					50	
					24"	2
					50	
F	.440"	11.18mm	34	2"	12"	2
					50	
G	.440"	11.18mm	34	3"	18"	2
					50	
					24"	2
					50	
H	.500"	12.70mm	N/A	2"	12"	2
					50	
I	.500"	12.70mm	N/A	3"	18"	2
					50	
					24"	2
					50	

BRUSH TIP DETAILS SHOWN ACTUAL SIZE




ACRYLIC TIP BRUSHES

Channel cleaning brushes featuring a protective acrylic tip for an extra level of protection during cleaning. Antimicrobial nylon bristles on a flexible handle. No loop end, may be used as a pull through brush. FR indicates French catheter size.

Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg	
A	AT-12-062	.062"	1.58mm	5	1"	12"	2
	AT-12-062-50						50
	AT-18-062	.062"	1.58mm	5	1"	18"	2
	AT-18-062-50						50
	AT-24-062	.062"	1.58mm	5	1"	24"	2
	AT-24-062-50						50
B	AT-12-087	.087"	2.20mm	7	1"	12"	2
	AT-12-087-50						50
	AT-18-087	.087"	2.20mm	7	1"	18"	2
	AT-18-087-50						50
	AT-24-087	.087"	2.20mm	7	1"	24"	2
	AT-24-087-50						50
C	AT-12-105	.105"	2.67mm	8	1"	12"	2
	AT-12-105-50						50
	AT-18-105	.105"	2.67mm	8	1"	18"	2
	AT-18-105-50						50
	AT-24-105	.105"	2.67mm	8	1"	24"	2
	AT-24-105-50						50
D	AT-12-118	.118"	3.00mm	9	1"	12"	2
	AT-12-118-50						50
	AT-18-118	.118"	3.00mm	9	1"	18"	2
	AT-18-118-50						50
	AT-24-118	.118"	3.00mm	9	1"	24"	2
	AT-24-118-50						50
E	AT-12-125	.125"	3.18mm	10	1"	12"	2
	AT-12-125-50						50
	AT-18-125	.125"	3.18mm	10	1"	18"	2
	AT-18-125-50						50
	AT-24-125	.125"	3.18mm	10	1"	24"	2
	AT-24-125-50						50
F	AT-12-131	.131"	3.33mm	10	1"	12"	2
	AT-12-131-50						50
	AT-18-131	.131"	3.33mm	10	1"	18"	2
	AT-18-131-50						50
	AT-24-131	.131"	3.33mm	10	1"	24"	2
	AT-24-131-50						50

BRUSH TIP DETAILS SHOWN ACTUAL SIZE




ACRYLIC TIP BRUSHES (CONTINUED)

Channel cleaning brushes featuring a protective acrylic tip for an extra level of protection during cleaning. Antimicrobial nylon bristles on a flexible handle. No loop end, may be used as a pull through brush. FR indicates French catheter size.

Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A AT-12-158	.158"	4.00mm	12	1.5"	12"	2
AT-12-158-50						50
AT-18-158	.158"	4.00mm	12	1.5"	18"	2
AT-18-158-50						50
AT-24-158	.158"	4.00mm	12	1.5"	24"	2
AT-24-158-50						50
B AT-12-197	.197"	5.00mm	15	2"	12"	2
AT-12-197-50						50
AT-18-197	.197"	5.00mm	15	2"	18"	2
AT-18-197-50						50
AT-24-197	.197"	5.00mm	15	2"	24"	2
AT-24-197-50						50
C AT-12-250	.250"	6.35mm	19	2"	12"	2
AT-12-250-50						50
AT-18-250	.250"	6.35mm	19	2"	18"	2
AT-18-250-50						50
AT-24-250	.250"	6.35mm	19	2"	24"	2
AT-24-250-50						50
D AT-12-276	.276"	7.00mm	21	2"	12"	2
AT-12-276-50						50
AT-18-276	.276"	7.00mm	21	2"	18"	2
AT-18-276-50						50
AT-24-276	.276"	7.00mm	21	2"	24"	2
AT-24-276-50						50
E AT-12-375	.375"	9.52mm	29	3"	12"	2
AT-12-375-50						50
AT-18-375	.375"	9.52mm	29	3"	18"	2
AT-18-375-50						50
AT-24-375	.375"	9.52mm	29	3"	24"	2
AT-24-375-50						50
F AT-12-440	.440"	11.18mm	34	3"	12"	2
AT-12-440-50						50
AT-18-440	.440"	11.18mm	34	3"	18"	2
AT-18-440-50						50
AT-24-440	.440"	11.18mm	34	3"	24"	2
AT-24-440-50						50

BRUSH TIP DETAILS SHOWN ACTUAL SIZE




POLYPROPYLENE BRUSHES

Channel cleaning brushes featuring rigid polypropylene bristles on a twisted stainless steel handle with loop end. FR indicates equivalent French catheter size.

	Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A	PP-16-039	.039"	1.00mm	3	1"	16"	2
	PP-16-039-50						50
	PP-22-039	.039"	1.00mm	3	1"	22"	2
	PP-22-039-50						50
B	PP-16-062	.062"	1.58mm	5	1"	16"	2
	PP-16-062-50						50
	PP-22-062	.062"	1.58mm	5	1"	22"	2
	PP-22-062-50						50
C	PP-16-087	.087"	2.20mm	7	1"	16"	2
	PP-16-087-50						50
	PP-22-087	.087"	2.20mm	7	1"	22"	2
	PP-22-087-50						50
D	PP-16-105	.105"	2.67mm	8	2"	16"	2
	PP-16-105-50						50
	PP-22-105	.105"	2.67mm	8	2"	22"	2
	PP-22-105-50						50
E	PP-16-125	.125"	3.18mm	10	2"	16"	2
	PP-16-125-50						50
	PP-22-125	.125"	3.18mm	10	2"	22"	2
	PP-22-125-50						50
F	PP-16-158	.158"	4.00mm	12	2"	16"	2
	PP-16-158-50						50
	PP-22-158	.158"	4.00mm	12	2"	22"	2
	PP-22-158-50						50
G	PP-18-197	.197"	5.00mm	15	2"	18"	2
	PP-18-197-50						50
	PP-22-197	.197"	5.00mm	15	2"	22"	2
	PP-22-197-50						50

BRUSH TIP DETAILS SHOWN ACTUAL SIZE



POLYPROPYLENE BRUSHES (CONTINUED)

Channel cleaning brushes featuring rigid polypropylene bristles on a twisted stainless steel handle with loop end (pull through styles excluded). FR indicates equivalent French catheter size.

Number	Diameter inches	Diameter mm	FR	Bristle Area Length	Length	Quantity/pkg
A PP-18-250 PP-18-250-50	.250"	6.35mm	19	2"	18"	2 50
PP-24-250 PP-24-250-50	.250"	6.35mm	19	2"	24"	2 50
B PP-18-393 PP-18-393-50	.393"	10.00mm	30	2"	18"	2 50
C PP-18-591 PP-18-591-50	.591"	15.00mm	N/A	2"	18"	2 50
D PP-18-750 PP-18-750-50	.750"	19.05mm	N/A	6"	18"	2 50
E PT-28-118 PT-28-118-50 (pull through)	.118"	3.00mm	9	2"	28"	2 50
F PT-28-197 PT-28-197-50 (pull through)	.197"	5.00mm	15	2"	28"	2 50





PIPE CLEANERS - 50' ROLLS

Available in either cotton or polyester. Rolls allow you to cut to the size you need, minimizing waste. Single use.

Number	Material	Diameter inches	Diameter mm	Quantity	
A	PC-100C	cotton	.125"	3.18mm	1 roll/pkg
B	PC-200C	cotton	.250"	6.35mm	1 roll/pkg
C	PC-300P	polyester	.125"	3.18mm	1 roll/pkg
D	PC-400P	polyester	.250"	6.35mm	1 roll/pkg

PIPE CLEANERS - 12" LENGTHS

E	PC-112C	cotton	.125"	3.18mm	100/pkg
F	PC-212C	cotton	.250"	6.35mm	100/pkg
G	PC-312P	polyester	.125"	3.18mm	100/pkg
H	PC-412P	polyester	.250"	6.35mm	100/pkg

DOUBLE CLEANER

Soft cotton and nylon cleaner. Single use.

Number	Diameter inches	Diameter mm	Length	Quantity	
I	PC-503	.118"	3mm	6"	100/pkg

FLEECE CLEANER

Fleece cleaner for extra delicate items. Single use.

J	PC-512	.500"	12.7mm	12"	50/pkg
---	--------	-------	--------	-----	--------

FLEX CLEANER

Polyester cleaner with nylon bristles in stem. Single use.

K	PC-515	.625"	15.87mm	12"	20/pkg
---	--------	-------	---------	-----	--------

FLEX ULTRA FINE CLEANER

Polyester cleaner with nylon bristles in stem. Rolls allow you to cut to the size you need, minimizing waste. Single use.

L	PC-606	.236"	6mm	6' roll	1 roll/pkg
---	--------	-------	-----	---------	------------

DIAMETERS SHOWN ACTUAL SIZE

TOOTHBRUSH-STYLE CLEANING BRUSHES

Antimicrobial nylon, brass, or stainless bristles on a sturdy, plastic handle. Designed for external cleaning.

Number	Details	Overall Length	Quantity/pkg
A	N-2121 N-2121-50	7"	2 50
B	S-2121 S-2121-50	7"	2 50
C	N-2100 N-2100-50	7"	3 50
D	B-2100 B-2100-50	7"	3 50
E	S-2100 S-2100-50	7"	3 50
F	S-2300 S-2300-50	7"	3 50
G	S-2131 S-2131-50	7"	3 50
H	N-2221 N-2221-50	6.89"	2 50
I	N-2200	9.25"	2
DOUBLE HEADED TOOTHBRUSH-STYLE			
J	D-2100BN D-2100BN-50	7"	3 50
K	D-2100SN D-2100SN-50	7"	3 50



Instrument Cleaning – General Instrument Cleaning Brushes



GENERAL INSTRUMENT CLEANING BRUSHES

Antimicrobial nylon or stainless steel bristles on a sturdy, plastic handle. Designed for external cleaning.

	Number	Details	Overall Length	Quantity
A	N-3000	nylon	7"	1/pkg
	N-3000-3			3/pkg
B	N-4000	nylon	7.25"	1/pkg
C	S-4000	stainless steel	7.25"	1/pkg

ACETABULAR REAMER BRUSH

Designed for cleaning acetabular bone reamer heads. Rigid, antimicrobial nylon bristles. Double-ended. Round 'pigtail' style brush on one end, straight precision style brush on the other.

	Number	Details	Overall Length	Quantity
D	N-2000	double-ended	7.725"	1/pkg
	N-2000-3			3 /pkg

GENERAL INSTRUMENT CLEANING BRUSH - STAINLESS STEEL BRISTLES

Stainless steel bristles on a twisted stainless steel handle with loop end.

	Number	Details	Diameter inches	Diameter mm	Bristle Area Length	Length	Quantity
E	BR-2300	stainless steel	.438"	11.13mm	2"	6"	2/pkg

GENERAL INSTRUMENT CLEANING BRUSH - STAINLESS STEEL BRISTLES

Stainless steel bristles on a rigid stainless steel handle. No loop. Designed to reach small, challenging spaces.

	Number	Details	Overall Length	Quantity
F	BR-04-003	stainless steel	4"	1/pkg
	BR-04-003-50			50/pkg

DOUBLE-ENDED BRUSH


Soak-proof cleaning brush. Nylon bristles on twisted stainless steel handle with rosebud tip.

	Number	Brush Diameter	Bristle Area Length	Overall Length	Quantity
G	BR-5450	.500" & 1"	2" & 3.25"	12"	1/pkg


FINGER LOOP BRUSH

Nylon bristles on twisted stainless steel handle with rosebud tip. Rigid wire with a finger loop end.


	Number	Brush Diameter	Bristle Area Length	Overall Length	Quantity
H	BR-5455	1"	3"	10.5"	1/pkg



A




B



E



C



F





ENDOSCOPE BRUSHES

A variety of single, double, and triple bristled endoscope brushes as well as valve brushes and kits.

Channel Brushes

Number	Details	Brush Diameter	Bristle Area Length	Length	Quantity
A	EB10100 double-ended	10mm & 5mm	30mm & 20mm	2500mm	50/box
B	EB10103 double-ended	6mm	20mm	2500mm	50/box
C	EB10103.3 triple-cleaning	6mm	20mm	2500mm	50/box
D	EB10102 double-ended	3mm	20mm	2300mm	50/box
E	EB10104 double-ended with handle	10mm & 5mm	40mm & 20mm	2300mm	25/box
F	EB10105 double-ended with pipe cleaner	5mm & 5mm	125mm & 20mm	2300mm	50/box
G	EB10106 channel brush with handle	2mm	10mm	1200mm	50/box
H	EB10136 double-ended	3mm	20mm	1200mm	50/box






A

ENDOSCOPE BRUSHES

A variety of endoscope brushes conveniently packaged together with a valve brush.

Brush Kits


Number	Details	Quantity
A EB10103.3C	includes triple-cleaning channel EB10103.3 and valve brush EB10101	50/pkg
B EB10107C	includes double-ended channel EB10103 and valve brush EB10101	50/pkg
C EB10107G	includes double-ended channel EB10102 and valve brush EB10101	50/pkg



B

DID YOU KNOW?






Key Surgical is a global leading partner to hospitals and surgical/procedural facilities, supporting the needs in endoscopy, operating room, and instrument reprocessing. Our comprehensive portfolio is designed with a focus on patient safety; allowing our customers to efficiently source the high quality products they need with the best customer service and overall experience with every interaction. With our headquarters in the US, we also have locations in the UK and Germany and our products are distributed in over 70 countries around the world. Throughout our continued growth we have remained true to our mission; to be the easiest company in our industry with which to do business.



C

BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Endoscope Brushes



ENDOSCOPE BRUSHES

Nylon bristles attached to flexible, colored plastic tubing. Single-ended. Single use.


Number	Brush Diameter	Bristle Area Length	Length	Color	Quantity/pkg
A BR-075-18 BR-075-18-50	1.8mm	2cm	75cm	blue	1 50
B BR-075-24 BR-075-24-50	2.4mm	2cm	75cm	yellow	1 50
C BR-075-32 BR-075-32-50	3.2mm	2cm	75cm	green	1 50
D BR-075-42 BR-075-42-50	4.2mm	2cm	75cm	pink	1 50
E BR-100-18 BR-100-18-50	1.8mm	2cm	100cm	blue	1 50
F BR-100-24 BR-100-24-50	2.4mm	2cm	100cm	yellow	1 50
G BR-100-32 BR-100-32-50	3.2mm	2cm	100cm	green	1 50
H BR-100-42 BR-100-42-50	4.2mm	2cm	100cm	pink	1 50
I BR-150-18 BR-150-18-50	1.8mm	2cm	150cm	blue	1 50
J BR-150-24 BR-150-24-50	2.4mm	2cm	150cm	yellow	1 50
K BR-150-32 BR-150-32-50	3.2mm	2cm	150cm	green	1 50
L BR-150-42 BR-150-42-50	4.2mm	2cm	150cm	pink	1 50
M BR-200-18 BR-200-18-50	1.8mm	2cm	200cm	blue	1 50
N BR-200-24 BR-200-24-50	2.4mm	2cm	200cm	yellow	1 50
O BR-200-32 BR-200-32-50	3.2mm	2cm	200cm	green	1 50
P BR-200-42 BR-200-42-50	4.2mm	2cm	200cm	pink	1 50

ENDOSCOPE BRUSHES

Tapered cleaning brushes. Nylon bristles attached to flexible plastic tubing. Single-ended. Single use.

Number	Brush Diameter	Bristle Area Length	Length	Quantity/pkg
Q BR-3308 BR-3308-50	tapered 3mm to 5mm	1cm	240cm	1 50
R BR-3309 BR-3309-50	tapered 3mm to 5mm	1cm	27cm	1 50
S BR-3310 BR-3310-50	tapered 3mm to 5mm	2cm	240cm	1 50

BRUSH TIP DETAILS SHOWN ACTUAL SIZE



DOUBLE-ENDED ENDOSCOPE BRUSHES (TAPERED, DOUBLE-ENDED)

Nylon bristles attached to flexible plastic tubing. Double-ended. Single use.

	Number	Brush Diameter	Bristle Area Length	Length	Quantity
A	BR-3320	tapered 3mm to 5mm	1cm	240cm	1/pkg
	BR-3320-50				50/pkg
B	BR-3330	tapered 3mm to 5mm	2cm	240cm	1/pkg
	BR-3330-50				50/pkg

DOUBLE-ENDED ENDOSCOPE BRUSHES (NON-TAPERED, DOUBLE-ENDED)

Nylon bristles attached to flexible plastic tubing. 3mm brush on one end, 6mm brush on the other. Single use.

C	BR-3340	3mm & 6mm	2cm	240cm	1/pkg
	BR-3340-50				50/pkg

ENDOSCOPE BRUSHES (TAPERED, SINGLE-ENDED)

Nylon bristles attached to flexible plastic tubing. Single-ended. Single use.

	Number	Brush Diameter	Bristle Area Length	Length	Quantity
D	BR-3711-50	tapered 4mm-6mm	1cm	240cm	50/pkg
	BR-3711-100				100/pkg (10 bags of 10 brushes)
E	BR-3799-50	Brush Kit, BR-3711 with BR-3201		240cm	50/pkg


PULL-THROUGH ENDOSCOPE BRUSHES

Delicate pull-through cleaning brushes. Nylon bristles attached to twisted stainless steel wire. Single-ended. Single use.

F	BR-3030	1.1mm	1cm	65cm	1/pkg
	BR-3030-50				50/pkg
	BR-3035	1.1mm	1cm	230cm	1/pkg
	BR-3035-50				50/pkg
G	BR-3040	2.0mm	1cm	120cm	1/pkg
	BR-3040-50				50/pkg
H	BR-3410	3.0mm	1cm	120cm	1/pkg
	BR-3410-50				50/pkg
I	BR-3420	3.0mm	1cm	160cm	1/pkg
	BR-3420-50				50/pkg

BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Endoscope Brushes



VALVE BRUSHES

A variety of short brushes designed for cleaning valve ports and suction wells on an endoscope. Single use.

Number	Details	Brush Diameter	Bristle Area Length	Quantity
A BR-3201 BR-3201-50	suction well brush	11mm	4.5cm	3/pkg 50/pkg
B BR-3202 BR-3202-50	double-ended	4.3mm & 11mm	1cm & 4.5cm	3/pkg 50/pkg
C EB10101	double-ended	10mm & 5mm	30mm & 20mm	100/box
D BR-3307	finger loop style	11mm	4.5cm	1/pkg

ALBARRAN BRUSH

Delicate brush designed for cleaning elevator channel on a duodenoscope. Single use.

Number	Details	Brush Diameter	Bristle Area Length	Quantity
E EB10108	Albarran Brush	2.3mm	13mm	100/box

PREMIUM ENDOSCOPE BRUSHES

Nylon bristles on twisted stainless steel wire or spring coil for hard-to-turn corners. Loop handle at one end and stainless steel beaded tip at bristle end. Reusable.

Number	Details	Brush Diameter	Bristle Area Length	Length	Quantity
F BR-3300	spring coil	tapered 4.0mm to 6.0mm	1cm	246cm	1/pkg
G BR-3305	spring coil	tapered 3.0mm to 5.0mm	1cm	128cm	1/pkg
H BR-3701	twisted wire	2.4mm	6mm	95cm	1/pkg
I BR-3702	spring coil	6.0mm	5mm	95cm	1/pkg
J BR-3703	twisted wire	2.4mm	6mm	165cm	1/pkg




A LITTLE BIT ABOUT..

Scot Milchman

Job Title: **Chief Executive Officer**

What is something you've always wanted to try but been too afraid to? **Paragliding**. What 3 historical characters would you invite to dinner? **Lincoln, Einstein, & Socrates**. What are some causes you care about? **Causes that help support and improve the lives of those with disabilities and underprivileged children**. What cheers you up? **Just about everything**. If you could be any age for a week, what age would it be and why? **90**. **Why not experience something I haven't before?** What would the title of your autobiography be? **Life's Too Short for Terrible Wine**.

BRUSH TIP DETAILS SHOWN ACTUAL SIZE



ENDOSCOPE BRUSH ASSORTMENTS

Number	Includes one of each of the following brushes:	Quantity
A BR-075-99	BR-075-18, BR-075-24, BR-075-32, BR-075-42	4/pkg
B BR-100-99	BR-100-18, BR-100-24, BR-100-32, BR-100-42	4/pkg
C BR-150-99	BR-150-18, BR-150-24, BR-150-32, BR-150-42	4/pkg
D BR-200-99	BR-200-18, BR-200-24, BR-200-32, BR-200-42	4/pkg

Laparoscopic Brush Assortment


Number	Includes one package of each of the following brushes:
E BR-4699	BR-12-118, BR-16-197, BR-16-393, BR-3202, N-2121


Endoscopic Brush Assortment

Number	Includes one package of each of the following brushes:
F BR-4799	BR-12-079, BR-12-118, BR-12-197, BR-16-197, BR-16-393, BR-24-197, BR-24-393, BR-3202, N-2121


Scope Brush Assortment

Number	Includes one package of each of the following brushes:
G BR-4899	BR-3202, BR-3307, BR-3308, BR-3340







A



B




C



D



E



F

BRUSH ASSORTMENTS

Original Brush Assortment

- Number Includes one package of each of the following brushes:
- A BR-4199 BR-18-062, BR-18-087, BR-18-105, BR-18-125, BR-18-250, FT-18-500, N-2100, N-2121, S-2100

Small Diameter Channel Brush Assortment

- Number Includes one package of each of the following brushes:
- B BR-4299 BR-08-158, BR-12-079, BR-18-062, BR-18-070, BR-18-087, BR-18-105, BR-18-125, BR-24-197, BR-24-200

Large Diameter Channel Brush Assortment

- Number Includes one package of each of the following brushes:
- C BR-4399 BR-06-375, BR-08-630, BR-16-276, BR-16-393, BR-18-250, BR-18-591, BR-18-880, BR-24-312, BR-24-440

Popular Channel Brush Assortment

- Number Includes one package of each of the following brushes:
- D BR-4499 BR-06-375, BR-08-158, BR-18-062, BR-18-070, BR-18-087, BR-18-105, BR-18-125, BR-18-250, BR-24-200

Micro-Ophthalmic Brush Assortment

- Number Includes one package of each of the following brushes:
- E BR-4599 BR-18-039, BR-18-062, BR-18-070, BR-18-087, BR-18-105, BR-18-125, BR-18-158, N-2100

Suction Tube Brush Assortment - French Sizes 3-12

- Number Includes one package of each of the following brushes:
- F BR-4999 BR-08-092, BR-08-131, BR-08-144, BR-12-079, BR-12-118, BR-12-158, BR-18-039, BR-18-062, BR-18-070, BR-18-105



A ACTUAL LENGTH IS 31"



B ACTUAL LENGTH IS 36"



C

D



E ACTUAL LENGTH IS 28"



F ACTUAL LENGTH IS 24"

STERILIZER CLEANING BRUSH

Brush allows for greater reach and cleaning ability. Bristles are at a forward angle to loosen soil without marring surface.

Number	Details	Dimensions	Quantity
A BR-5000	double-surface brush, polyester bristles	6" x 5" brush head, 31" overall length	1/pkg

STERILIZER CLEANING BRUSH HANDLE

White plastic threaded handle for use with BR-5400.

Number	Details	Overall Length	Quantity
B BR-5200	threaded handle, plastic	36"	1/pkg

STERILIZER CLEANING BRUSH - SCRUB HEAD

Solid white plastic block eliminates bacteria build-up and won't crack, chip, or splinter. Flare-head design assures thorough cleaning action from any angle. Handle sold separately.

Number	Details	Dimensions	Quantity
C BR-5400	flare-head brush, nylon bristles	5.5" x 9.5" brush head	1/pkg

PISTOL-GRIP BRUSH

Soak-proof design with pistol-grip handle. Crimped nylon bristles on a block handle that resists chipping and cracking.

Number	Bristle Length	Bristle Area	Overall Length	Quantity
D BR-4000	2"	4" x 4"	8"	1/pkg

STERILIZER DRAIN CLEANING BRUSH - WITH PLASTIC HANDLE

White nylon bristles, stainless steel wire and soak-proof plastic handle.

Number	Details	Dimensions	Quantity
E BR-5500	long-handled brush, nylon bristles	1.37" diameter, 28" overall length	1/pkg

STERILIZER DRAIN CLEANING BRUSH - WITH STAINLESS STEEL HANDLE

Designed for cleaning out drains, lines, and openings. Stainless steel wire with stiff, crimped nylon bristles. Loop handle.

Number	Details	Dimensions	Quantity
F BR-5600	drain/line brush, nylon bristles	2" diameter, 24" overall length	1/pkg




STERILIZER CLEANING BRUSH ASSORTMENT

Assortment includes 1 of each of the following items.

Number	Details
BR-6099	BR-5000, BR-5200, BR-5400, BR-5500, and BR-5600




Instrument Cleaning – Brush Organizer and French Catheter Scale


A


TWO BARS INCLUDED IN BR-7000.


B


SHOWN IN USE. BRUSHES NOT INCLUDED.

C


BRUSH ORGANIZER

Stainless steel organizer with a built-in drip tray. Includes 12 permanent pegs and 2 additional pegged bars (4 pegs per bar). For easy organization and storage display of cleaning brushes.

Number	Details	Quantity
A BR-7000	16" x 35" (1.33' x 2.91') wall-mounted organizer	1/pkg

BRUSH ORGANIZER WITH BRUSH ASSORTMENT

B BR-7100	wall-mounted organizer with Original Brush Assortment	
------------------	---	--

Assortment includes 1 package of the following brushes:

BR-18-062	BR-18-087	BR-18-105
BR-18-125	BR-18-250	FT-18-500
N-2100	N-2121	S-2100

PEGGED BAR FOR BRUSH ORGANIZER


Customize your organizer to fit your needs. Stainless steel bar with 4 pegs. Adhesive-backed for universal placement. Bar measures 1" wide by 5" long.

Number	Details	Quantity
C BR-7050	pegged bar with adhesive backing	1/pkg

FRENCH CATHETER SCALE

Convenient two-sided french catheter scale with built-in 6" ruler.

Number	Details	Quantity
D MJ-370	plastic card with die-cut holes	4/pkg



INSTRUMENT PROTECTION

Instrument Tip Caps are designed to extend the life of your instruments. Our tip caps protect your valuable instruments from unnecessary wear or damage, helping reduce repair costs while giving you the peace of mind that your instrument investment is being protected. In addition to protecting and extending the life of your valuable instruments, tip caps also help to protect sterilization packs from puncture and protect employees from sharp instruments.

Available in a variety of shapes and sizes to fit all your instrument protection needs, we have the perfect tip cap for each of your various instruments – from skin hooks to scissors to osteotomes – and they are available with or without vents, depending on your preference.

Our tip caps are radiopaque, single use, and compatible with steam sterilization.

Instrument Sleeves are designed to protect various delicate instruments using an easy-to-fold design. Sleeves are single use, and compatible with steam and gas plasma sterilization.

Duraholder® Instrument Protection System (IPS) is available in several configurations to fit your instrument protection needs within a closed container or set. Duraholders help you organize and prepare your various instruments for steam sterilization in one convenient pouch.

Round Caps

Flat Caps

Osteotome Caps

Tapered Caps

Skin Hook Caps

Airflow Caps™

Self-Retaining Retractor Caps

Twin Caps

Minimally Invasive Surgical (MIS) Instrument Caps

Instrument Sleeves

Duraholder Instrument Protection System (IPS)

Please note: Unless specified, all products are not made with natural rubber latex.



Instrument Protection – Round Caps



PRODUCTS SHOWN ACTUAL SIZE

ROUND CAPS

Radiopaque. Single use. Compatible with steam sterilization. 100 caps per package.

Number	Color	Details	Diameter inches	Length inches	Diameter mm	Length mm
A	3-01-01	white	.062"	.75"	1.6mm	19.1mm
B	3-02-01	white tint	.062"	.75"	1.6mm	19.1mm
C	3-03-01	clear	.062"	.75"	1.6mm	19.1mm
D	3-04-01	clear with vents	.062"	.75"	1.6mm	19.1mm
E	3-08-01	white tint with vents	.062"	.75"	1.6mm	19.1mm
F	3-15-01	white with vents	.062"	.75"	1.6mm	19.1mm
G	3-01-02	blue	.078"	.75"	2.0mm	19.1mm
H	3-02-02	blue tint	.078"	.75"	2.0mm	19.1mm
I	3-03-02	clear	.078"	.75"	2.0mm	19.1mm
J	3-04-02	clear with vents	.078"	.75"	2.0mm	19.1mm
K	3-08-02	blue tint with vents	.078"	.75"	2.0mm	19.1mm
L	3-15-02	blue with vents	.078"	.75"	2.0mm	19.1mm
M	3-01-03	green	.109"	.75"	2.8mm	19.1mm
N	3-02-03	green tint	.109"	.75"	2.8mm	19.1mm
O	3-03-03	clear	.109"	.75"	2.8mm	19.1mm
P	3-04-03	clear with vents	.109"	.75"	2.8mm	19.1mm
Q	3-08-03	green tint with vents	.109"	.75"	2.8mm	19.1mm
R	3-15-03	green with vents	.109"	.75"	2.8mm	19.1mm
S	3-01-04	red	.125"	1"	3.2mm	25.4mm
T	3-02-04	red tint	.125"	1"	3.2mm	25.4mm
U	3-03-04	clear	.125"	1"	3.2mm	25.4mm
V	3-04-04	clear with vents	.125"	1"	3.2mm	25.4mm
W	3-08-04	red tint with vents	.125"	1"	3.2mm	25.4mm
X	3-15-04	red with vents	.125"	1"	3.2mm	25.4mm
Y	3-01-05	yellow	.187"	1"	4.7mm	25.4mm
Z	3-02-05	yellow tint	.187"	1"	4.7mm	25.4mm
AA	3-03-05	clear	.187"	1"	4.7mm	25.4mm
BB	3-04-05	clear with vents	.187"	1"	4.7mm	25.4mm
CC	3-08-05	yellow tint with vents	.187"	1"	4.7mm	25.4mm
DD	3-15-05	yellow with vents	.187"	1"	4.7mm	25.4mm
EE	3-01-06	purple	.406"	.75"	10.3mm	19.1mm
FF	3-02-06	purple tint	.406"	.75"	10.3mm	19.1mm
GG	3-03-06	clear	.406"	.75"	10.3mm	19.1mm
HH	3-04-06	clear with vents	.406"	.75"	10.3mm	19.1mm
II	3-08-06	purple tint with vents	.406"	.75"	10.3mm	19.1mm
JJ	3-15-06	purple with vents	.406"	.75"	10.3mm	19.1mm
KK	3-04-98	clear with vents	.155"	.71"	3.9mm	18.0mm
LL	3-04-99	clear with vents	.187"	1.50"	4.7mm	38.1mm
MM	3-05-27	red	.312"	1.50"	7.9mm	38.1mm
NN	3-05-28	navy	.437"	1"	11.1mm	25.4mm
OO	3-05-29	navy	1.25"	.25"	31.8mm	6.4mm
PP	3-05-30	navy	.156"	1"	4.0mm	25.4mm



FLAT CAPS

Radiopaque. Single use. Compatible with steam sterilization. 100 caps per package.

Number	Color	Details	Width inches	Length inches	Width mm	Length mm
A	orange		.375"	1"	9.5mm	25.4mm
B	orange tint		.375"	1"	9.5mm	25.4mm
C	clear		.375"	1"	9.5mm	25.4mm
D	clear	with vents	.375"	1"	9.5mm	25.4mm
E	orange tint	with vents	.375"	1"	9.5mm	25.4mm
F	orange	with vents	.375"	1"	9.5mm	25.4mm
G	brown		.625"	1"	15.9mm	25.4mm
H	brown tint		.625"	1"	15.9mm	25.4mm
I	clear		.625"	1"	15.9mm	25.4mm
J	clear	with vents	.625"	1"	15.9mm	25.4mm
K	brown tint	with vents	.625"	1"	15.9mm	25.4mm
L	brown	with vents	.625"	1"	15.9mm	25.4mm
M	black		1"	1"	25.4mm	25.4mm
N	black tint		1"	1"	25.4mm	25.4mm
O	clear		1"	1"	25.4mm	25.4mm
P	clear	with vents	1"	1"	25.4mm	25.4mm
Q	black tint	with vents	1"	1"	25.4mm	25.4mm
R	black	with vents	1"	1"	25.4mm	25.4mm
3-08-100	assortment (1 package of each)					
3-08-01, 3-08-02, 3-08-03, 3-08-04, 3-08-05, 3-08-06, 3-08-07, 3-08-08, 3-08-09, 3-08-11, 3-08-12, 3-08-13, 3-08-15, 3-08-17, 3-08-18, 3-08-19						

PRODUCTS SHOWN ACTUAL SIZE





OSTEOTOME CAPS

Radiopaque. Single use. Compatible with steam sterilization. 100 caps per package.

Number	Color	Details	Width inches	Length inches	Width mm	Length mm
A	clear		.25"	.75"	6.4mm	19.1mm
B	clear	with vents	.25"	.75"	6.4mm	19.1mm
C	white tint	with vents	.25"	.75"	6.4mm	19.1mm
D	clear		.375"	1"	9.5mm	25.4mm
E	clear	with vents	.375"	1"	9.5mm	25.4mm
F	blue tint	with vents	.375"	1"	9.5mm	25.4mm
G	clear		.50"	1"	12.7mm	25.4mm
H	clear	with vents	.50"	1"	12.7mm	25.4mm
I	green tint	with vents	.50"	1"	12.7mm	25.4mm
J	clear		.70"	1"	17.8mm	25.4mm
K	clear	with vents	.70"	1"	17.8mm	25.4mm
L	yellow tint	with vents	.70"	1"	17.8mm	25.4mm
M	clear		1"	1"	25.4mm	25.4mm
N	clear	with vents	1"	1"	25.4mm	25.4mm
O	black tint	with vents	1"	1"	25.4mm	25.4mm
P	clear	with vents	1.25"	1.25"	31.8mm	31.8mm
Q	purple tint	with vents	1.25"	1.25"	31.8mm	31.8mm
R	clear	with vents	1.50"	1.50"	38.1mm	38.1mm
S	red tint	with vents	1.50"	1.50"	38.1mm	38.1mm
T	clear	with vents	1.75"	1.25"	44.5mm	31.8mm
U	orange tint	with vents	1.75"	1.25"	44.5mm	31.8mm

PRODUCTS SHOWN ACTUAL SIZE



TAPERED CAPS

Radiopaque. Single use. Compatible with steam sterilization. 50 caps per package.

Number	Color	Details	Diameter inches	Length inches	Diameter mm	Length mm
A	white tint	with vents	.06" to .13"	1.50"	1.5mm to 3.3mm	38.1mm
B	blue tint	with vents	.09" to .17"	1.50"	2.5mm to 4.3mm	38.1mm
C	green tint	with vents	.14" to .21"	1.50"	3.5mm to 5.3mm	38.1mm
D	red tint	with vents	.18" to .25"	1.50"	4.5mm to 6.2mm	38.1mm
E	yellow tint	with vents	.22" to .29"	1.50"	5.5mm to 7.2mm	38.1mm
F	purple tint	with vents	.25" to .32"	1.50"	6.5mm to 8.2mm	38.1mm
G	orange tint	with vents	.29" to .36"	1.50"	7.5mm to 9.2mm	38.1mm
H	brown tint	with vents	.33" to .40"	1.50"	8.5mm to 10.2mm	38.1mm
I	black tint	with vents	.37" to .44"	1.50"	9.4mm to 11.2mm	38.1mm
J	white tint	with vents	.08" to .12"	2"	2mm to 3.1mm	50.8mm
K	blue tint	with vents	.10" to .19"	2"	2.5mm to 4.9mm	50.8mm
L	green tint	with vents	.14" to .23"	2"	3.5mm to 5.8mm	50.8mm
M	red tint	with vents	.18" to .27"	2"	4.5mm to 6.9mm	50.8mm
N	yellow tint	with vents	.22" to .31"	2"	5.5mm to 7.8mm	50.8mm
O	purple tint	with vents	.25" to .35"	2"	6.5mm to 8.9mm	50.8mm
P	orange tint	with vents	.29" to .38"	2"	7.5mm to 9.8mm	50.8mm
Q	brown tint	with vents	.33" to .42"	2"	8.5mm to 10.8mm	50.8mm
R	black tint	with vents	.37" to .46"	2"	9.4mm to 11.8mm	50.8mm
S	green tint	with vents	.53" to .58"	1.50"	13.4mm to 14.7mm	38.1mm
T	yellow tint	with vents	.61" to .66"	1.50"	15.4mm to 16.8mm	38.1mm
U	orange tint	with vents	.68" to .74"	1.50"	17.4mm to 18.8mm	38.1mm
V	green tint	with vents	.53" to .60"	2"	13.4mm to 15.2mm	50.8mm
W	yellow tint	with vents	.61" to .68"	2"	15.4mm to 17.2mm	50.8mm
X	orange tint	with vents	.68" to .76"	2"	17.4mm to 19.3mm	50.8mm
Y	blue tint		.250"	1.25"	6.35mm	31.8mm
Z	purple tint		.400"	1.90"	10.16mm	48.3mm
AA	orange tint		.450"	1.75"	11.43mm	44.5mm

PRODUCTS SHOWN ACTUAL SIZE

Instrument Protection – Skin Hook and Airflow Caps



SKIN HOOK CAPS

Protect delicate single and double prong skin hook instrumentation. Radiopaque. Single use. Compatible with steam sterilization. 100 caps per package.

Number	Color	Details	Width inches	Length inches	Width mm	Length mm
A 3-10-01	clear	single prong, 2.5-3.5mm hook	.190"	.75"	4.8mm	19.1mm
B 3-10-02	clear	single prong, 5mm hook	.250"	.75"	6.3mm	19.1mm
C 3-10-03	clear	double prong, 2mm apart	.165"	.75"	4.1mm	19.1mm
D 3-10-04	clear	double prong, 5mm apart	.265"	.75"	6.7mm	19.1mm
E 3-10-05	clear	double prong, 7mm apart	.360"	.75"	9.1mm	19.1mm

AIRFLOW CAPS™

An alternative to vented caps, grips on the inside of the cap allow more airflow during sterilization process and assist in easy removal. Radiopaque. Single use. Compatible with steam sterilization. 100 caps per package.

Number	Color	Diameter mm	Length inches	Diameter mm	Length mm
F 3-16-05	clear	.187"	1"	4.7mm	25.4mm
G 3-16-06	clear	.406"	.75"	10.3mm	19.1mm

Number	Color	Width inches	Length inches	Width mm	Length mm
H 3-16-07	clear	.375"	1"	9.5mm	25.4mm
I 3-16-08	clear	.625"	1"	15.8mm	25.4mm
J 3-16-09	clear	1"	1"	25.4mm	25.4mm

PRODUCTS SHOWN ACTUAL SIZE

Instrument Protection – Self-Retaining Retractor Caps



SELF-RETAINING RETRACTOR CAPS

Protect a variety of self-retaining retractors including Weitlaner, D'Errico-Adson, Beckman, and Adson retractors. Radiopaque. Single use. Compatible with steam sterilization. 100 caps per package.


Number	Color	Width inches	Length inches	Width mm	Length mm
A	blue tint	.850"	.80"	21.6mm	20.3mm
B	purple tint	1"	1"	25.4mm	25.4mm
C	orange tint	1"	1.25"	25.4mm	31.8mm
D	green tint	1"	1.75"	25.4mm	44.5mm



GELPI RETRACTOR CAPS

Protect a variety of Gelpi retractors. Radiopaque. Single use. Compatible with steam sterilization. 100 caps per package.

Number	Color	Diameter inches	Length inches	Diameter mm	Length mm
E	clear	.312"	1.50"	7.9mm	38.1mm
F	clear	.312"	2"	7.9mm	50.8mm
G	clear	.312"	2.50"	7.9mm	63.5mm
H	blue tint	.312"	1.50"	7.9mm	38.1mm
I	purple tint	.312"	2"	7.9mm	50.8mm
J	orange tint	.312"	2.50"	7.9mm	63.5mm



PRODUCTS SHOWN ACTUAL SIZE



TWIN CAPS

Protect instruments while keeping them in an open position for cleaning. Perfect for scissors and forceps. Radiopaque. Single use. Compatible with steam sterilization. 50 caps per package.

Number	Color	Details	Width inches	Length inches	Width mm	Length mm
A	white tint	with vents	.062"	.75"	1.6mm	19.1mm
B	blue tint	with vents	.078"	.75"	2.0mm	19.1mm
C	green tint	with vents	.109"	.75"	2.8mm	19.1mm
D	red tint	with vents	.125"	1"	3.2mm	25.4mm
E	yellow tint	with vents	.187"	1"	4.7mm	25.4mm
F	orange tint	with vents	.375"	1"	9.5mm	25.4mm
3-06-50	assortment (1 package of each) 3-06-01, 3-06-02, 3-06-03, 3-06-04, 3-06-05, 3-06-07					

MINIMALLY INVASIVE SURGICAL INSTRUMENT (MIS) CAPS

MIS Instrument Caps are designed to protect endoscopes, probes, veress needles, and trocars during sterilization. Radiopaque. Single use. Compatible with steam sterilization. 50 caps per package.

Number	Color	Details	Diameter inches	Length inches	Diameter mm	Length mm
G	clear	with vents	.17"	1.13"	4.3mm	28.7mm
H	clear	with vents	.17"	2"	4.3mm	50.8mm
I	clear	with vents	.36"	2.60"	9.1mm	66mm

DID YOU KNOW?

We believe that in addition to having the right products at the right time in your department, a vital part of the journey to improving patient safety and positive surgical outcomes is continual education. Key Surgical delivers training and education for industries of sterile processing, operating room, and endoscopy. While we are known for providing education in-person through in-services or at regional/national conferences, our Education Program is also online-based and CEs (continuing education credits) can be earned for participation. Whether you're on the go or on the job, we'll help you customize a learning pathway that tracks your progress and helps you work smarter not harder. Registration is free and so are the credits. From cutting edge to basic education, you can access our ever-growing library through our website or the Key Surgical Education App.

PRODUCTS SHOWN ACTUAL SIZE



PRODUCTS SHOWN ACTUAL SIZE

INSTRUMENT SLEEVES


Easy-to-fold protective sleeves conveniently protect your various delicate instruments during sterilization. Available in three versatile sizes. Single use. Compatible with steam and gas plasma sterilization. 1000 sleeves per package.

Number	Details	Width inches	Length inches	Width cm	Length cm
A CS-1010W*	tear-resistant, fits instruments up to 10mm wide	1"	3.5"	3cm	9cm
B CS-2025W*	tear-resistant, fits instruments up to 25mm wide	2"	6.3"	5cm	16cm
C CS-3050W*	tear-resistant, fits instruments up to 50mm wide	3"	8.8"	8cm	22cm



⚠️ *WARNING: This product can expose you to chemicals including benzophenone, which is known in the State of California to cause cancer. For more information go to www.p65warnings.ca.gov.





INSTRUMENTS NOT INCLUDED

DURAHOLDER® INSTRUMENT PROTECTION SYSTEM (IPS)

Designed to help protect and organize your delicate surgical instruments within a wrapped set or container system. Single use. Compatible with steam sterilization. Case of 100.

Number	Details	Length inches	Width inches	Length cm	Width cm
A CS-DH100	2 rows of 6 pockets	17.5"	11.5"	44.45cm	29.21cm
B CS-DH102	2 rows of 6 pockets with slots in pocket	17.5"	11.5"	44.45cm	29.21cm
C CS-DH200	2 rows of 6 notched pockets with tape closure strips	18.5"	10.25"	47cm	26cm
D CS-DH300	2 rows of 6 pockets	18"	16"	46cm	41cm
E CS-DH505	single pocket	5.5"	5.25"	14cm	13.3cm
F CS-DH512	single pocket	5"	12"	12.7cm	30.5cm



TAGS

Locking Tags are designed to help you identify your instruments, trays, and carts as they move from sterilization to staging areas to procedures. Our tags allow you to indicate special handling instructions, secure trays, and carts from tampering and quickly identify priority instruments for repair and sterilization. Tags are available in a wide range of styles, colors, and locking configurations. Printed tags with special handling instructions, as well as uniquely numbered tags, are available as a standard product. Need a customized tag? We can do that for you!

Breakaway Tags are also available in a range of colors and can be written on with a permanent lab marker. They must be cut or broken to open, providing a secure lock.

Write-On Tags and Smart Tags allow you to write specific handling instructions on each tag when more detailed communication is necessary. These tags attach to instruments and trays with ties and are larger and more visible than our other tags, ensuring that detailed instructions will be seen by all. Smart Tags include detachable barcodes for even easier follow-up.

Scope Tags are designed specifically with the needs of scope processing and care in mind. Perfect for recording important details and available in seven bright colors for easy color-coding. Attached loops make for a convenient one-piece tag.

Metal Tray Tags easily clip onto the edge of a tray and provide ample space for identification, tracking and barcoding needs. Stainless steel for durability.

Locking Tags

Case Locking Tags

Breakaway Tags

Write-On Tags

Smart Tags

Scope Tags



Plastic Tray Tags

Metal Tray Tags

One-Loop Tags



Please note: Unless specified, all products are not made with natural rubber latex.



PRODUCTS SHOWN ACTUAL SIZE

LOCKING TAGS

Versatile, durable, and economical. Tamper-resistant. Compatible with steam and gas plasma sterilization. Single use. Customizable. 1" length.

Number	Details	Color	Quantity
A LT-150* LT-150-1000*	CLEAN	white	100/pkg 1000/pkg
B LT-210* LT-210-1000*	QUICK TURNOVER	blue	100/pkg 1000/pkg
C LT-310* LT-310-1000*	LOANER	green	100/pkg 1000/pkg
D LT-450* LT-450-1000*	REPAIR	red	100/pkg 1000/pkg
E LT-550* LT-550-1000*	SHARPEN	yellow	100/pkg 1000/pkg
F LT-720* LT-720-1000*	BROKEN	orange	100/pkg 1000/pkg
G LT-100N* LT-100N-1000*	uniquely numbered	white	100/pkg 1000/pkg
H LT-300N* LT-300N-1000*	uniquely numbered	green	100/pkg 1000/pkg
I LT-400N* LT-400N-1000*	uniquely numbered	red	100/pkg 1000/pkg
J LT-500N* LT-500N-1000*	uniquely numbered	yellow	100/pkg 1000/pkg
K LT-900N* LT-900N-1000*	uniquely numbered	black	100/pkg 1000/pkg

WARNING: This product can expose you to chemicals including titanium dioxide, which is known in the State of California to cause cancer. For more information go to www.p65warnings.ca.gov.



PRODUCTS SHOWN ACTUAL SIZE

LOCKING TAGS (CONTINUED)

Versatile, durable, and economical. Tamper-resistant. Compatible with steam and gas plasma sterilization. Single use. Customizable. 1" length.


Number	Color	Quantity
A LT-100*	white	100/pkg
LT-100-1000*		1000/pkg
B LT-110*	pink	100/pkg
LT-110-1000*		1000/pkg
C LT-112*	gold	100/pkg
LT-112-1000*		1000/pkg
D LT-116*	purple	100/pkg
LT-116-1000*		1000/pkg
E LT-200*	blue	100/pkg
LT-200-1000*		1000/pkg
F LT-300*	green	100/pkg
LT-300-1000*		1000/pkg
G LT-400*	red	100/pkg
LT-400-1000*		1000/pkg
H LT-500*	yellow	100/pkg
LT-500-1000*		1000/pkg
I LT-900*	black	100/pkg
LT-900-1000*		1000/pkg
J LT-600*	gray	100/pkg
LT-600-1000*		1000/pkg
K LT-700*	orange	100/pkg
LT-700-1000*		1000/pkg
L LT-800*	brown	100/pkg
LT-800-1000*		1000/pkg

INDICATOR DOT LOCKING TAGS

Indicator dots on tags change color to indicate exposure to steam sterilization (blue dot to dark). Versatile, durable, and economical. Tamper-resistant. Compatible with steam sterilization. Single use. 1" length.

Number	Color	Quantity
M LS-100*	white	100/pkg
LS-100-1000*		1000/pkg
N LS-400*	red	100/pkg
LS-400-1000*		1000/pkg

⚠️ *WARNING: This product can expose you to chemicals including titanium dioxide, which is known in the State of California to cause cancer. For more information go to www.p65warnings.ca.gov.



PRODUCTS SHOWN ACTUAL SIZE

CASE LOCKING TAGS

A versatile, durable, and economical option for identifying and securing cases and containers. Tamper-resistant. Compatible with steam and gas plasma sterilization. Single use. 7/8" operational length. Smallest size slot tag will fit is 1/16" x 3/16".

Number	Color	Quantity
A LT-001*	white	100/pkg
LT-001-1000*		1000/pkg
B LT-002*	blue	100/pkg
LT-002-1000*		1000/pkg
C LT-003*	green	100/pkg
LT-003-1000*		1000/pkg
D LT-004*	red	100/pkg
LT-004-1000*		1000/pkg
E LT-007*	orange	100/pkg
LT-007-1000*		1000/pkg




CASE LOCKING TAGS WITH INDICATOR DOT

Indicator dots on tags changes color to indicate exposure to steam sterilization (blue dot to dark). Fits well with rigid containers. Tamper-resistant. Compatible with steam sterilization. Single use. 7/8" operational length.

Number	Color	Quantity
F LS-007*	orange	100/pkg
LS-007-1000*		1000/pkg



***WARNING:** This product can expose you to chemicals including titanium dioxide, which is known in the State of California to cause cancer. For more information go to www.p65warnings.ca.gov.



PRODUCTS SHOWN 70% ACTUAL SIZE

BREAKAWAY TAGS

Versatile, durable, and economical. Tamper-resistant. Tags cinch closed until tight and must be cut or broken to open, providing a secure lock. Select styles are compatible with steam and gas plasma sterilization. Single use.

Number	Color	Details	Length	Quantity
A LR-102*	blue	sterilizable	4"	100/pkg
B LR-103*	green	sterilizable	4"	100/pkg
C LR-104*	red	sterilizable	4"	100/pkg
D LR-105*	yellow	sterilizable	4"	100/pkg
E LR-103N*	green	sterilizable, uniquely numbered	4"	100/pkg
F LR-105N*	yellow	sterilizable, uniquely numbered	4"	100/pkg
G LR-312*	blue	sterilizable	5.5"	100/pkg
H LR-313*	green	sterilizable	5.5"	100/pkg
I LR-314*	red	sterilizable	5.5"	100/pkg
J LR-315*	yellow	sterilizable	5.5"	100/pkg
K LR-322*	blue	non-sterilizable, uniquely numbered	6"	100/pkg
L LR-323*	green	non-sterilizable, uniquely numbered	6"	100/pkg
M LR-324*	red	non-sterilizable, uniquely numbered	6"	100/pkg
N LR-325*	yellow	non-sterilizable, uniquely numbered	6"	100/pkg
O LR-327*	orange	non-sterilizable, uniquely numbered	6"	100/pkg
P LR-329*	black	non-sterilizable, uniquely numbered	6"	100/pkg



⚠️ *WARNING: This product can expose you to chemicals including titanium dioxide, which is known in the State of California to cause cancer. For more information go to www.p65warnings.ca.gov.



Tags



DOUBLE-SIDED

WRITE-ON TAGS ARE 3" x 5"

WRITE-ON TAGS

Durable and economical. Identify and record specific instructions with these versatile tags. Tamper-resistant. Compatible with steam and gas plasma sterilization. Single use. Customizable. 3" x 5". Each package includes 100 cable ties and a lab marker.

Number	Details	Color	Quantity/pkg	
A	WT-0009	BLANK	black	100 tags / 100 ties / 1 lab marker
B	WT-0102	VENDOR	blue	100 tags / 100 ties / 1 lab marker
C	WT-0203	LOANER	green	100 tags / 100 ties / 1 lab marker
D	WT-0304	PRIORITY	red	100 tags / 100 ties / 1 lab marker
E	WT-0405	RE-RUN	yellow	100 tags / 100 ties / 1 lab marker
F	WT-0507	BROKEN (double-sided)	orange	100 tags / 100 ties / 1 lab marker
G	WT-1911	MISSING	pink	100 tags / 100 ties / 1 lab marker
H	RT-100	BROKEN DO NOT USE	white	100 tags / 100 ties / 1 lab marker
I	RT-200	BLANK	white	100 tags / 100 ties / 1 lab marker



Tags

Tags – Write-On Tags, Lab Markers, and Ties



A DOUBLE-SIDED



B DOUBLE-SIDED



C DOUBLE-SIDED

WRITE-ON TAGS ARE 3" x 5"

WRITE-ON TAGS (CONTINUED)

Durable and economical. Identify and record specific instructions with these versatile tags. Tamper-resistant. Compatible with steam and gas plasma sterilization. Single use. Customizable. 3" x 5". Each package includes 100 cable ties and a lab marker.

Number	Details	Color	Quantity/pkg
A WT-2614	LOANER CHECK IN	light green	100 tags / 100 ties / 1 lab marker
B WT-2715	VENDOR CHECK IN	light blue	100 tags / 100 ties / 1 lab marker
C WT-3217	QUICK TURNOVER	blue	100 tags / 100 ties / 1 lab marker

LAB MARKERS

Ink is waterproof, alcohol resistant, and compatible with steam and gas plasma sterilization.

Number	Details	Color	Quantity
D RT-209P	lab marker	black	12/pkg

CABLE TIES

Ties are compatible with steam and gas plasma sterilization.


Number	Details	Color	Quantity
E RT-100T	ties (6")	white	1000/pkg



D



E



SMART TAGS ARE 2" x 5.5"

SMART TAGS

Durable, lockable tags have tear-off barcode for additional tracking. Identify and record specific instructions on these one piece versatile tags. Tamper-resistant. Compatible with steam and gas plasma sterilization. Single use. Each package includes a lab marker.

	Number	Details	Color	Length	Width	Quantity/pkg
A	ST-0009*	BLANK	black	5.5"	2"	50 tags / 1 lab marker
B	ST-0102*	VENDOR	blue	5.5"	2"	50 tags / 1 lab marker
C	ST-0203*	LOANER	green	5.5"	2"	50 tags / 1 lab marker
D	ST-0304*	PRIORITY	red	5.5"	2"	50 tags / 1 lab marker
E	ST-0405*	RE-RUN	yellow	5.5"	2"	50 tags / 1 lab marker
F	ST-0507*	BROKEN	orange	5.5"	2"	50 tags / 1 lab marker



 ***WARNING:** This product can expose you to chemicals including titanium dioxide, which is known in the State of California to cause cancer. For more information go to www.p65warnings.ca.gov.





FRONT



BACK



SCOPE TAGS

Double-sided tags designed to loop around a scope to identify process dates. Tear-resistant tags available in seven colors for easy color coding and include the month, day, and date selection (1-31) on one side of the tag. Compatible with steam and gas plasma sterilization. Single use. Customizable. Each package includes a lab marker.

Number	Color	Length	Width	Quantity/pkg	
A	ET-101*	white	13.25"	1.5"	100 tags / 1 lab marker
B	ET-102*	blue	13.25"	1.5"	100 tags / 1 lab marker
C	ET-103*	green	13.25"	1.5"	100 tags / 1 lab marker
D	ET-104*	red	13.25"	1.5"	100 tags / 1 lab marker
E	ET-105*	yellow	13.25"	1.5"	100 tags / 1 lab marker
F	ET-106*	purple	13.25"	1.5"	100 tags / 1 lab marker
G	ET-107*	orange	13.25"	1.5"	100 tags / 1 lab marker
H	ET-104-1*	red - BIOHAZARD	13.25"	1.5"	100 tags / 1 lab marker

ET-104-1 has BIOHAZARD printed on both sides. Does not include month, day, and date selection.



***WARNING:** This product can expose you to chemicals including benzophenone, which is known in the State of California to cause cancer. For more information go to www.p65warnings.ca.gov.

Scope Model #: _____

Serial #: _____

Preclean by: _____ **Date:** _____ **Time:** _____

Leak Test: Pass _____ Fail _____

Sterilizer Type: _____ **Number:** _____

Cycle #: _____ **Date:** _____ **Time:** _____

Parameters Met: Yes _____ No _____

Alcohol Flush: Yes _____ No _____

Air Flush: Yes _____ No _____

Stored: Yes _____ No _____

Hang Date: _____ **Time:** _____

PLACE BARCODE HERE

A



SCOPE PROCESSING TAG

Double-sided tag designed to be looped around a scope to capture detailed information such as scope model #, leak test (pass/fail), hang date/time, or a barcode. Tear-resistant. Compatible with steam and gas plasma sterilization. Single use. Customizable. Each package includes a lab marker.

Number	Color	Length	Width	Quantity/pkg
A ET-200*	white	10.5"	3"	100 tags / 1 lab marker

SCOPE REPAIR TAG


Double-sided tag designed to be looped around a scope to identify repair needs. Tear-resistant. Compatible with steam and gas plasma sterilization. Single use. Customizable. Each package includes a lab marker.

Number	Color	Length	Width	Quantity/pkg
B ET-204*	red	10.5"	3"	100 tags / 1 lab marker

REPAIR

KEYSURGICAL TO REORDER: 800.541.7995

www.keysurgical.com ET-204



B

PULL SCOPE FOR SERVICE

SERIAL NUMBER/IDENTIFIER: _____

DATE & TIME OF ISSUE: _____


DEPT./PROCEDURE ROOM: _____

DESCRIPTION OF DEFECT/ISSUE: _____

INITIALS: _____



WARNING: This product can expose you to chemicals including benzophenone, which is known in the State of California to cause cancer. For more information go to www.p65warnings.ca.gov.



WHITE TAGS SHOWN ACTUAL SIZE

PLASTIC TRAY TAGS

Durable, plastic tags designed for use on instrument trays. Tray tag provides a flat surface for barcode placement. Great for color coding. Compatible with steam and gas plasma sterilization. Tags include a built-in clip on the back.

Number	Details	Length	Width	Quantity
A	PT-01-01* with clip, white	1.18"	1.18"	10/pkg
B	PT-01-02* with clip, blue	1.18"	1.18"	10/pkg
C	PT-01-03* with clip, green	1.18"	1.18"	10/pkg
D	PT-01-04* with clip, red	1.18"	1.18"	10/pkg
E	PT-01-05* with clip, yellow	1.18"	1.18"	10/pkg
F	PT-01-06* with clip, purple	1.18"	1.18"	10/pkg
G	PT-01-07* with clip, orange	1.18"	1.18"	10/pkg
H	PT-01-09* with clip, black	1.18"	1.18"	10/pkg
I	PT-01-14* with clip, light blue	1.18"	1.18"	10/pkg
J	PT-02-01* with clip, white	2.16"	1.18"	10/pkg
K	PT-02-02* with clip, blue	2.16"	1.18"	10/pkg
L	PT-02-03* with clip, green	2.16"	1.18"	10/pkg
M	PT-02-04* with clip, red	2.16"	1.18"	10/pkg
N	PT-02-05* with clip, yellow	2.16"	1.18"	10/pkg
O	PT-02-06* with clip, purple	2.16"	1.18"	10/pkg
P	PT-02-07* with clip, orange	2.16"	1.18"	10/pkg
Q	PT-02-09* with clip, black	2.16"	1.18"	10/pkg
R	PT-02-14* with clip, light blue	2.16"	1.18"	10/pkg
S	PT-04-01* with clip, white	2.76"	1.18"	10/pkg
T	PT-04-02* with clip, blue	2.76"	1.18"	10/pkg
U	PT-04-03* with clip, green	2.76"	1.18"	10/pkg
V	PT-04-04* with clip, red	2.76"	1.18"	10/pkg
W	PT-04-05* with clip, yellow	2.76"	1.18"	10/pkg
X	PT-04-06* with clip, purple	2.76"	1.18"	10/pkg
Y	PT-04-07* with clip, orange	2.76"	1.18"	10/pkg
Z	PT-04-09* with clip, black	2.76"	1.18"	10/pkg
AA	PT-04-14* with clip, light blue	2.76"	1.18"	10/pkg
BB	PT-05-01* with clip, white	3.54"	1.18"	10/pkg
CC	PT-05-02* with clip, blue	3.54"	1.18"	10/pkg
DD	PT-05-03* with clip, green	3.54"	1.18"	10/pkg
EE	PT-05-04* with clip, red	3.54"	1.18"	10/pkg
FF	PT-05-05* with clip, yellow	3.54"	1.18"	10/pkg
GG	PT-05-06* with clip, purple	3.54"	1.18"	10/pkg
HH	PT-05-07* with clip, orange	3.54"	1.18"	10/pkg
II	PT-05-09* with clip, black	3.54"	1.18"	10/pkg
JJ	PT-05-14* with clip, light blue	3.54"	1.18"	10/pkg



***WARNING:** This product can expose you to chemicals including titanium dioxide, which is known in the State of California to cause cancer. For more information go to www.p65warnings.ca.gov.





TAGS SHOWN ACTUAL SIZE

PLASTIC TRAY TAGS (CONTINUED)

Durable, plastic tags designed for use on instrument trays. Tray tag provides a flat surface for barcode placement. Great for color coding. Compatible with steam and gas plasma sterilization. Tags designed to be attached to a tray using a fastener such as a ring or clip.

Number	Details	Length	Width	Quantity	
A	PT-03-01*	white	2.36"	1.57"	10/pkg
B	PT-03-02*	blue	2.36"	1.57"	10/pkg
C	PT-03-03*	green	2.36"	1.57"	10/pkg
D	PT-03-04*	red	2.36"	1.57"	10/pkg
E	PT-03-05*	yellow	2.36"	1.57"	10/pkg
F	PT-03-06*	purple	2.36"	1.57"	10/pkg
G	PT-03-07*	orange	2.36"	1.57"	10/pkg
H	PT-03-09*	black	2.36"	1.57"	10/pkg
I	PT-03-14*	light blue	2.36"	1.57"	10/pkg
J	PT-06-01*	white	3.54"	1.30"	10 /pkg
K	PT-06-02*	blue	3.54"	1.30"	10/pkg
L	PT-06-03*	green	3.54"	1.30"	10/pkg
M	PT-06-04*	red	3.54"	1.30"	10/pkg
N	PT-06-05*	yellow	3.54"	1.30"	10/pkg
O	PT-06-06*	purple	3.54"	1.30"	10/pkg
P	PT-06-07*	orange	3.54"	1.30"	10/pkg
Q	PT-06-09*	black	3.54"	1.30"	10/pkg
R	PT-06-14*	light blue	3.54"	1.30"	10/pkg

TRAY TAG RING

Durable, stainless steel ring designed to be used with Plastic Tray Tags to attach the tag to the tray.

Number	Details	Quantity
S	TG-100 stainless steel ring, .79" (20mm) diameter	100/pkg



TRAY TAG CLIP

Durable, stainless steel clip designed to be used with Plastic Tray Tags to attach to a tray. Specifically fits Plastic Tray Tags PT-06-01 through PT-06-14.

Number	Details	Quantity
T	TG-200 stainless steel clip	100/pkg



⚠️ WARNING: This product can expose you to chemicals including titanium dioxide, which is known in the State of California to cause cancer. For more information go to www.p65warnings.ca.gov.




A



B




C



D



FRONT



E

BACK

METAL TRAY TAGS

Stainless steel tags easily clip onto the side of a tray or container. Flat surface allows for placement of barcode, label, or other identifying instructions. Available in five sizes. Compatible with steam sterilization.

Number	Length inches	Width inches	Quantity
A MT-0303	1.18"	1.18"	1/pkg
MT-0303-10			10/pkg
B MT-0553	2.17"	1.18"	1/pkg
MT-0553-10			10/pkg
C MT-0703	2.76"	1.18"	1/pkg
MT-0703-10			10/pkg
D MT-0903	3.55"	1.18"	1/pkg
MT-0903-10			10/pkg
E MT-0905	3.55"	1.97"	1/pkg
MT-0905-10			10/pkg



PRODUCT SHOWN ACTUAL SIZE

CAN FIT A LABEL UP TO 3.5" x 2" IN SIZE.



ONE-LOOP TAGS

Designed to loop around tray handles to identify through color-coding or providing a flat surface for barcode placement. Compatible with steam and EtO sterilization. Single use.

	Number	Color	Length	Width	Quantity
A	PL-201	white	9.25"	2"	250/pkg
B	PL-202	blue	9.25"	2"	250/pkg
C	PL-203	green	9.25"	2"	250/pkg
D	PL-204	red	9.25"	2"	250/pkg
E	PL-205	yellow	9.25"	2"	250/pkg
F	PL-207	orange	9.25"	2"	250/pkg
G	PL-211	pink	9.25"	2"	250/pkg
H	PL-101	white	8"	4.25"	250/pkg



CAN FIT A LABEL UP TO 4.25" x 2.375" IN SIZE.

IDENTIFICATION TAPE

Our Identification Tape is designed to help you organize and identify instruments as they move between sterile processing and various departments. Using tape can help save you time, effort, and money by offering effective organization and identification of your surgical instrumentation. It is easy to apply and is compatible with steam sterilization. We have the broadest selection of tape available and it is colorful, easy to use, and economical.

Roll Tape is available in nine solid colors and three widths, while our Sheet Tape is available in over 250 colors and patterns. We are happy to help you customize your own identification tape as well!

Adhesive Remover is designed to dissolve tape residue on instruments, trays, basins, and other surfaces. Adhesive Remover is CARB compliant and biodegradable. Available in either (1) 32oz bottle or a case of (6) 4oz bottles. Features a citrus-scent.

Roll Tape

Adhesive Remover

Solid Color Sheet Tape

Patterned Sheet Tape

PLEASE NOTE: Unless specified, all products are not made with natural rubber latex.





ROLL TAPE

Compatible with steam sterilization. 250" total length. Available in three widths.

Color	1/8" Width	1/4" Width	3/8" Width
A white	T-01-01	T-02-01	T-03-01
B blue	T-01-02	T-02-02	T-03-02
C green	T-01-03	T-02-03	T-03-03
D red	T-01-04	T-02-04	T-03-04
E yellow	T-01-05	T-02-05	T-03-05
F purple	T-01-06	T-02-06	T-03-06
G orange	T-01-07	T-02-07	T-03-07
H brown	T-01-08	T-02-08	T-03-08
I black	T-01-09	T-02-09	T-03-09

ADHESIVE REMOVER

Key Surgical Adhesive Remover is designed to dissolve tape residue on instruments, trays, basins, or other surfaces. Adhesive Remover is CARB compliant, and biodegradable. Available in either (1) 32oz bottle or a case of (6) 4oz bottles. Features a citrus-scent.

Number	Details	Quantity/pkg
J MJ-7004	4oz	6 bottles/box
K MJ-7032	32oz	1 bottle




SOLID COLOR SHEET TAPE – 1/8" WIDTH


Compatible with steam sterilization. Customizable.

Each 8.5" x 11" sheet includes:

- 462 pieces of 1/8" x 1"
- 88 pieces of 1/8" x 1.5"
- 22 pieces of 1/8" x 2"
- 22 pieces of 1/8" x 5"
- (total of 748" of tape)



Number	Color
A T-110-01	white
B T-110-02	blue
C T-110-03	green
D T-110-04	red
E T-110-05	yellow
F T-110-06	purple
G T-110-07	orange
H T-110-08	brown
I T-110-09	black
J T-110-10	bright orange
K T-110-11	pink
L T-110-12	gray
M T-110-13	lime green
N T-110-14	light blue
O T-110-15	lavender
P T-110-16	beige
Q T-110-17	plum




SOLID COLOR SHEET TAPE – 1/4" WIDTH


Compatible with steam sterilization. Customizable.

Each 8.5" x 11" sheet includes:

- 231 pieces of 1/4" x 1"
- 44 pieces of 1/4" x 1.5"
- 11 pieces of 1/4" x 2"
- 11 pieces of 1/4" x 5"
- (total of 374" of tape)



Number	Color
A T-120-01	white
B T-120-02	blue
C T-120-03	green
D T-120-04	red
E T-120-05	yellow
F T-120-06	purple
G T-120-07	orange
H T-120-08	brown
I T-120-09	black
J T-120-10	bright orange
K T-120-11	pink
L T-120-12	gray
M T-120-13	lime green
N T-120-14	light blue
O T-120-15	lavender
P T-120-16	beige
Q T-120-17	plum
R T-120-21	emerald green



Identification Tape – Patterned Sheet Tape


PATTERNED SHEET TAPE – 1/4" WIDTH

Compatible with steam sterilization. Customizable. One 8.5" x 11" sheet pre-cut into 1", 1.5", 2", and 5" lengths. Cut as illustrated on page 111.



PATTERNED SHEET TAPE – 1/4" WIDTH


Compatible with steam sterilization. Customizable. One 8.5" x 11" sheet pre-cut into 1", 1.5", 2", and 5" lengths. Cut as illustrated on page 111.



Identification Tape – Patterned Sheet Tape



PATTERNED SHEET TAPE – 1/4" WIDTH

Compatible with steam sterilization. Customizable. One 8.5" x 11" sheet pre-cut into 1", 1.5", 2", and 5" lengths. Cut as illustrated on page 111.

 CV CV CV CV T-501-02 blue/CV	 CV CV CV CV T-501-04 red/CV				
 ER ER ER ER T-502-02 blue/ER	 ER ER ER ER T-502-03 green/ER	 ER ER ER ER T-502-04 red/ER	 ER ER ER ER T-502-06 purple/ER	 ER ER ER ER T-502-07 orange/ER	
 GI GI GI GI T-503-1609 beige/GI	 OB OB OB OB T-504-11 pink/OB	 OB OB OB OB T-504-14 light blue/OB	 OR OR OR OR T-523-03 green/OR	 OR OR OR OR T-523-07 orange/OR	 OR OR OR OR T-523-20 medical blue/ OR
 ASC ASC ASC ASC T-601-02 blue/ASC	 ASC ASC ASC ASC T-601-04 red/ASC	 ASC ASC ASC ASC T-601-07 orange/ASC	 CCU CCU CCU CCU T-602-1009 bright orange/ CCU	 DEN DEN DEN DEN T-603-0109 white/DEN	 DEN DEN DEN DEN T-603-04 red/DEN
 ENT ENT ENT ENT T-604-02 blue/ENT	 ENT ENT ENT ENT T-604-03 green/ENT	 ENT ENT ENT ENT T-604-04 red/ENT	 ENT ENT ENT ENT T-604-07 orange/ENT	 ENT ENT ENT ENT T-604-15 lavender/ENT	 TX TX TX TX T-506-02 blue/TX
 EYE EYE EYE EYE T-605-02 blue/EYE	 EYE EYE EYE EYE T-605-03 green/EYE	 EYE EYE EYE EYE T-605-04 red/EYE	 EYE EYE EYE EYE T-605-06 purple/EYE		

PATTERNED SHEET TAPE – 1/4" WIDTH

Compatible with steam sterilization. Customizable. One 8.5" x 11" sheet pre-cut into 1", 1.5", 2", and 5" lengths. Cut as illustrated on page 111.

 GEN GEN GEN GEN T-606-02 blue/GEN	 GEN GEN GEN GEN T-606-03 green/GEN	 GEN GEN GEN GEN T-606-04 red/GEN	 GEN GEN GEN GEN T-606-06 purple/GEN	 GEN GEN GEN GEN T-606-07 orange/GEN	 GEN GEN GEN GEN T-606-09 black/GEN
 GYN GYN GYN GYN T-607-03 green/GYN	 GYN GYN GYN GYN T-607-04 red/GYN	 GYN GYN GYN GYN T-607-06 purple/GYN	 GYN GYN GYN GYN T-607-07 orange/GYN	 GYN GYN GYN GYN T-607-19 fuchsia/GYN	 GYN GYN GYN GYN T-607-20 medical blue/GYN
 ICU ICU ICU ICU T-608-03 green/ICU	 ICU ICU ICU ICU T-608-04 red/ICU	 L&D L&D L&D L&D T-609-03 green/L&D	 L&D L&D L&D L&D T-609-04 red/L&D	 L&D L&D L&D L&D T-609-11 pink/L&D	 L&D L&D L&D L&D T-609-14 light blue/ L&D
 NEU NEU NEU NEU T-610-02 blue/NEU	 NEU NEU NEU NEU T-610-03 green/NEU	 NEU NEU NEU NEU T-610-04 red/NEU	 NEU NEU NEU NEU T-610-06 purple/NEU	 NEU NEU NEU NEU T-610-07 orange/NEU	 NEU NEU NEU NEU T-610-11 pink/NEU
 PLA PLA PLA PLA T-612-03 green/PLA	 PLA PLA PLA PLA T-612-04 red/PLA	 PLA PLA PLA PLA T-612-0509 yellow/PLA	 RAD RAD RAD RAD T-613-0709 orange/RAD	 SPN SPN SPN SPN T-614-02 blue/SPN	 SPN SPN SPN SPN T-614-04 red/SPN
 SPD SPD SPD SPD T-614-03 green/SPD	 SPD SPD SPD SPD T-614-07 orange/SPD	 SPD SPD SPD SPD T-614-17 plum/SPD	 PED PED PED PED T-611-14 light blue/PED	 PED PED PED PED T-611-18 baby pink/PED	 PED PED PED PED T-611-1309 lime green/ PED




Identification Tape

Identification Tape – Patterned Sheet Tape

PATTERNED SHEET TAPE – 1/4" WIDTH

Compatible with steam sterilization. Customizable. One 8.5" x 11" sheet pre-cut into 1", 1.5", 2", and 5" lengths. Cut as illustrated on page 111.

 T-616-03 green/URO	 T-616-04 red/URO	 T-616-06 purple/URO	 T-616-20 medical blue/URO	 T-617-04 red/VAS	 T-617-14 light blue/VAS
 T-622-03 green/H2O	 T-623-08 brown/BAR	 T-624-08 brown/LIV	 T-627-04 red/HIP		
 T-686-03 green/LOAN	 T-686-07 orange/LOAN	 T-686-20 medical blue/LOAN	 T-687-03 green/CLNC	 T-687-07 orange/CLNC	 T-687-20 medical blue/CLNC
 T-688-03 green/ANES	 T-688-07 orange/ANES	 T-688-20 medical blue/ANES	 T-689-03 green/NICU	 T-689-07 orange/NICU	 T-689-20 medical blue/NICU
 T-690-03 green/PACU	 T-690-07 orange/PACU	 T-690-20 medical blue/PACU			

PATTERNED SHEET TAPE – 1/4" WIDTH

Compatible with steam sterilization. Customizable. One 8.5" x 11" sheet pre-cut into 1", 1.5", 2", and 5" lengths. Cut as illustrated on page 111.

 T-701-0104 white/hearts	 T-701-02 blue/hearts	 T-701-03 green/hearts	 T-701-04 red/hearts	 T-701-06 purple/hearts	 T-701-07 orange/hearts
 T-701-09 black/hearts	 T-701-11 pink/hearts	 T-701-21 emerald green/hearts	 T-701-0509 yellow/hearts	 T-701-1309 lime green/hearts	 T-701-1509 lavender/hearts
 T-702-0104 white/bones	 T-702-02 blue/bones	 T-702-03 green/bones	 T-702-04 red/bones	 T-702-06 purple/bones	 T-702-07 orange/bones
 T-702-08 brown/bones	 T-702-09 black/bones	 T-702-11 pink/bones	 T-702-19 fuchsia/bones	 T-702-21 emerald green/bones	 T-702-0509 yellow/bones
 T-702-1309 lime green/bones	 T-702-1509 lavender/bones	 T-703-02 blue/stars	 T-703-03 green/stars	 T-703-04 red/stars	 T-703-06 purple/stars
 T-703-07 orange/stars	 T-703-19 fuchsia/stars	 T-703-20 medical blue/stars	 T-703-21 emerald green/stars		



Identification Tape

Identification Tape – Patterned Sheet Tape

PATTERNED SHEET TAPE – 1/4" WIDTH

Compatible with steam sterilization. Customizable. One 8.5" x 11" sheet pre-cut into 1", 1.5", 2", and 5" lengths. Cut as illustrated on page 111.

 T-704-06 purple/dots	 T-704-19 fuchsia/dots	 T-704-20 medical blue/dots	 T-704-21 emerald green/dots		
 T-705-06 purple/triangles	 T-705-19 fuchsia/triangles	 T-705-20 medical blue/ triangles	 T-705-21 emerald green/ triangles		
 T-706-02 blue/ smiley faces	 T-706-03 green/ smiley faces	 T-706-04 red/ smiley faces	 T-706-06 purple/ smiley faces	 T-706-0509 yellow/ smiley faces	 T-706-0709 orange/ smiley faces
 T-706-1109 pink/ smiley faces	 T-706-1309 lime green/ smiley faces	 T-706-1409 light blue/ smiley faces	 T-707-02 blue/diamonds	 T-707-04 red/diamonds	 T-707-17 plum/diamonds
 T-708-0104 white/no	 T-709-0103 white/clovers	 T-709-03 green/clovers	 T-710-0118 Key Surgical logo	 T-711-0109 white/yin yang	 T-712-0104 white/lips
 T-713-02 blue/footprints	 T-713-03 green/footprints	 T-713-04 red/footprints	 T-713-06 purple/footprints	 T-713-07 orange/footprints	 T-713-1109 pink/footprints

PATTERNED SHEET TAPE – 1/4" WIDTH

Compatible with steam sterilization. Customizable. One 8.5" x 11" sheet pre-cut into 1", 1.5", 2", and 5" lengths. Cut as illustrated on page 111.

 T-714-0504 yellow/nuclear	 T-715-0109 white/music notes	 T-716-03 green/crosses	 T-716-04 red/crosses	 T-720-02 blue/handprints	 T-720-04 red/handprints
 T-717-0103 white/pine trees	 T-718-0111 white/ribbon	 T-719-04 red/biohazard			
 T-401-0509 yellow/A	 T-402-0509 yellow/B	 T-403-0509 yellow/C	 T-404-0509 yellow/D	 T-405-0509 yellow/E	 T-407-0509 yellow/G
 T-412-0509 yellow/L	 T-414-0509 yellow/N	 T-416-0509 yellow/P	 T-416-1009 bright orange/P	 T-420-0509 yellow/T	 T-421-0509 yellow/U
 T-800-0509 yellow/O	 T-800-1009 bright orange/O	 T-801-0509 yellow/I	 T-802-0509 yellow/Z	 T-803-0509 yellow/3	 T-804-0509 yellow/4
 T-805-0509 yellow/5	 T-806-0509 yellow/6	 T-807-0509 yellow/7	 T-808-0509 yellow/8		



Identification Tape

PATTERNED SHEET TAPE – 1/4" WIDTH

Compatible with steam sterilization. Customizable. One 8.5" x 11" sheet pre-cut into 1", 1.5", 2", and 5" lengths. Each sheet contains numbers 1-44. Cut as illustrated on page 111.

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44
23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44
23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44
23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44

T-920-02
blue/1-44


1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44
23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44
23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44
23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44

T-920-04
red/1-44

PATTERNED SHEET TAPE – 1/4" WIDTH


Compatible with steam sterilization. Customizable. One 8.5" x 11" sheet pre-cut into 1" and 1.5" lengths. Cut as illustrated below.

Each 8.5" x 11" sheet includes:
308 pieces of 1/4" x 1"
44 pieces of 1/4" x 1.5"
(total of 374" of tape)




PATTERNED SHEET TAPE – 1/4" WIDTH

Compatible with steam sterilization. Customizable. One 8.5" x 11" sheet pre-cut into 1" and 1.5" lengths. Cut as illustrated on page 121.

 <p>T-401-19 red/A</p>	 <p>T-720-040109 black/red/card suits</p>
 <p>T-900-0109 white/0-8</p>	 <p>T-900-03 green/0-8</p>
 <p>T-900-06 purple/0-8</p>	 <p>T-900-07 orange/0-8</p>
 <p>T-900-19 fuchsia/0-8</p>	 <p>T-900-20 medical blue/0-8</p>
 <p>T-900-21 emerald green/0-8</p>	

PATTERNED SHEET TAPE – 1/4" WIDTH

Compatible with steam sterilization. Customizable. One 8.5" x 11" sheet pre-cut into 1" and 1.5" lengths. Cut as illustrated on page 121.

<table border="0"> <tr> <td>JAN</td><td>FEB</td><td>MAR</td><td>T-901-01</td> <td>APR</td><td>MAY</td><td>JUN</td><td>T-902-01</td> <td>JUL</td><td>AUG</td><td>SEP</td><td>T-903-01</td> <td>OCT</td><td>NOV</td><td>DEC</td><td>T-904-01</td> </tr> <tr> <td>JAN</td><td>FEB</td><td>MAR</td><td>white/Jan</td> <td>APR</td><td>MAY</td><td>JUN</td><td>green/Apr</td> <td>JUL</td><td>AUG</td><td>SEP</td><td>red/Jul</td> <td>OCT</td><td>NOV</td><td>DEC</td><td>lime green/Oct</td> </tr> <tr> <td>JAN</td><td>FEB</td><td>MAR</td><td>pink/Feb</td> <td>APR</td><td>MAY</td><td>JUN</td><td>orange/May</td> <td>JUL</td><td>AUG</td><td>SEP</td><td>brown/Aug</td> <td>OCT</td><td>NOV</td><td>DEC</td><td>yellow/Nov</td> </tr> <tr> <td>JAN</td><td>FEB</td><td>MAR</td><td>light blue/Mar</td> <td>APR</td><td>MAY</td><td>JUN</td><td>blue/Jun</td> <td>JUL</td><td>AUG</td><td>SEP</td><td>purple/Sep</td> <td>OCT</td><td>NOV</td><td>DEC</td><td>lavender/Dec</td> </tr> </table>	JAN	FEB	MAR	T-901-01	APR	MAY	JUN	T-902-01	JUL	AUG	SEP	T-903-01	OCT	NOV	DEC	T-904-01	JAN	FEB	MAR	white/Jan	APR	MAY	JUN	green/Apr	JUL	AUG	SEP	red/Jul	OCT	NOV	DEC	lime green/Oct	JAN	FEB	MAR	pink/Feb	APR	MAY	JUN	orange/May	JUL	AUG	SEP	brown/Aug	OCT	NOV	DEC	yellow/Nov	JAN	FEB	MAR	light blue/Mar	APR	MAY	JUN	blue/Jun	JUL	AUG	SEP	purple/Sep	OCT	NOV	DEC	lavender/Dec	 <p>T-910-02 blue/alphabet</p>
JAN	FEB	MAR	T-901-01	APR	MAY	JUN	T-902-01	JUL	AUG	SEP	T-903-01	OCT	NOV	DEC	T-904-01																																																		
JAN	FEB	MAR	white/Jan	APR	MAY	JUN	green/Apr	JUL	AUG	SEP	red/Jul	OCT	NOV	DEC	lime green/Oct																																																		
JAN	FEB	MAR	pink/Feb	APR	MAY	JUN	orange/May	JUL	AUG	SEP	brown/Aug	OCT	NOV	DEC	yellow/Nov																																																		
JAN	FEB	MAR	light blue/Mar	APR	MAY	JUN	blue/Jun	JUL	AUG	SEP	purple/Sep	OCT	NOV	DEC	lavender/Dec																																																		
	 <p>T-910-04 red/alphabet</p>																																																																
	 <p>T-910-19 fuchsia/alphabet</p>																																																																
	 <p>T-910-03 green/alphabet</p>																																																																
	 <p>T-910-0109 white/alphabet</p>																																																																



PERSONAL PROTECTIVE EQUIPMENT

Face Shields are available in several styles to fit the preferences of every individual. We offer full and short face shields, face shields with masks, with Tyvek™ expanders, with ties and with ear loops. All of our face shields offer different types of protection.

Barrier Gowns and Sleeves are lightweight, comfortable, and breathable.

Sterilizer Mitts are available in terry cloth or silicone. Cloth mitts provide you with a steam barrier for removing hot trays and racks from the sterilizer and are available in two lengths to protect most arm sizes. Silicone mitts are waterproof and heat resistant up to 480°F.

Nitrile gloves are powder-free, available in three lengths (9", 12", 24"). The 9" gloves have micro textured fingertips designed to increase dexterity and functionality and are incredibly durable. Our 24" nitrile gloves offer the most protection. They are treated with an antibacterial and anti-mildew agent and are also reusable.

Cleanroom Coveralls, Bouffant Caps, Masks, Beard Covers, and Shoe Covers are designed for use in cleanroom environments.

Adhesive Mats are designed to effectively capture dirt and dust particles from foot traffic.

Face Shields

Face Masks

Barrier Gowns

Protective Sleeves

Sterilizer Mitts

Nitrile Gloves

Cleanroom Apparel

Adhesive Mats

PLEASE NOTE: Unless specified, all products are not made with natural rubber latex.



Personal Protective Equipment – Face Shields



FULL FACE SHIELDS

Vented foam headband for increased airflow and comfort. Anti-fog, anti-glare mylar polyester film shields. Single use.

Number	Details	Quantity
A 4505	13" wide, 7.5" long, 1.625" wide foam headband	24/box
B 4506	13" wide, 9" long, 1.625" wide foam headband	24/box

FACE SHIELD WITH SPLASH GUARD


Vented foam headband for increased airflow and comfort. Anti-fog, anti-glare mylar polyester film shields with retractable Tyvek splash guard. Single use.

C 4509	13" wide, 5" long, 5" splash guard	24/box
--------	------------------------------------	--------

SELECT SHIELD

Disposable shields with reusable frames. Vented foam headband for increased airflow and comfort. Anti-fog, anti-glare mylar polyester film shields. Shields are single use.

D 4610	7.5" select shield starter kit	40 shields, 8 blue frames/box
E 4620	5.75" select shield starter kit	40 shields, 8 blue frames/box
F 4600	replacement frames	20/box
G 4615	7.5" replacement shields	40/box
H 4625	5.75" replacement shields	40/box





FOGFREE® MASKS WITH SHIELD

Fluid resistant outer layer with white facial tissue inner layer. Fog-free, non-woven barrier, distortion free wrap-around shield. Anti-glare mylar polyester film shields. Convenient ear loops. Single use.

Number	Details	BFE	PFE	Quantity
A 4508	blue mask with ear loops and shield	≥ 95% @ 3.0 micron	≥ 98% @ 0.5 micron	25/box
B 4517	blue mask with ear loops and shield	≥ 98% @ 3.0 micron	≥ 98% @ 0.1 micron	25/box

FACE MASK WITH SHIELD

Mask contains extra anti-fog foam strip and cellulose innerfacing with extra protection filter. Anti-fog, distortion free wrap-around shield. Anti-glare mylar polyester film shield. Form fitting ties. Single use.

C 4514	blue mask with fluid protection, anti-fog tie shield	> 99% @ 3.0 micron	99.4% @ 0.1 micron	25/box
--------	--	--------------------	--------------------	--------

1 - Bacterial Filtration Efficiency (BFE)

BFE is the effectiveness of a material to filter (prevent passage of) bacteria of a specified particle size, 3.0 µm (micron) ± µm. This is expressed as a percentage of a quantity that does not pass through the material.

2 - Particle Filtration Efficiency (PFE)

PFE is the effectiveness of a material to filter aerosol particles in the size range 0.1 µm to 5.0 µm. This is expressed as a percentage of a quantity that does not pass through the material.



Personal Protective Equipment – Barrier Gowns, Protective Sleeves, and Sterilizer Mitts



BARRIER GOWNS

Lightweight and comfortable protection. Fluid resistant. Polypropylene with a polyethylene coating. Fully closed back with closure ties at neck and waist with soft cloth cuffs for extra protection. Gowns wrapped individually in plastic. Single use.

	Number	Dimensions	Quantity
A	PS-110	40" (102cm) long	10/box
	PS-150		50/box
B	PS-160	57" (145cm) long	50/box
C	PS-170	63" (160cm) long	50/box

PROTECTIVE SLEEVES

Lightweight, strong, and durable. Sleeves are made from spun-bond polypropylene coated in polyethylene with elastic cuffs. Single use.

D	PS-52101	21" (53cm) long	200/box
	PS-51801	18" (46cm) long	200/box

STERILIZER MITTS - TERRY CLOTH

Built-in steam barrier for removing hot trays and racks from the sterilizer. Operates safely up to 500° F (260°C).

E	PS-3613	13" (33cm) long	1 pair
F	PS-3617	17" (43cm) long	1 pair

STERILIZER GLOVES - TERRY CLOTH

Built-in steam barrier for removing hot trays and racks from the sterilizer. Operates safely up to 250° F (121°C).

G	PS-3816	16" (39cm) long	1 pair
H	PS-3817	17" (43cm) long	1 pair

STERILIZER MITT - SILICONE

Waterproof mitt design allows firm grip for removing hot trays and racks from the sterilizer. Operates safely up to 480°F/ 248°C.

I	PS-3717	17" (43cm) long	1 mitt
---	---------	-----------------	--------





PPE SHOWN NOT INCLUDED

A



B



C



9" NITRILE GLOVES

Powder free, chemotherapy tested nitrile examination glove. They are ambidextrous and feature micro-textured fingertips for extra dexterity. Thermally activated by body heat for improved fit and comfort. Single use.

Number	Details	Quantity
A PS-3209	small	100/box
PS-3209-10		10 boxes/case
PS-3309	medium	100/box
PS-3309-10		10 boxes/case
PS-3409	large	100/box
PS-3409-10		10 boxes/case
PS-3509	x-large	100/box
PS-3509-10		10 boxes/case

12" NITRILE GLOVES

Powder free, nitrile examination glove. Gloves are made of synthetic polymers. Single use.

B PS-3212	small	100/box
PS-3212-10		10 boxes/case
PS-3312	medium	100/box
PS-3312-10		10 boxes/case
PS-3412	large	100/box
PS-3412-10		10 boxes/case
PS-3512	x-large	100/box
PS-3512-10		10 boxes/case
PS-3612	xx-large	100/box
PS-3612-10		10 boxes/case

24" NITRILE GLOVES

Extended sleeve provides arm and shoulder protection. Treated with an antibacterial and anti-mildew agent. Reusable.

C PS-3224	small	1 pair
PS-3324	medium	1 pair
PS-3424	large	1 pair
PS-3524	x-large	1 pair
PS-3624	xx-large	1 pair

GLOVE DISPENSER

Wall mount dispenser for personal protective equipment such as gloves. Durable ABS plastic construction. Easy to clean.

Number	Details	Width	Height	Depth	Quantity
D PS-7000	wall mount dispenser	12.5"	16"	4.5"	1/pkg



CLEANROOM COVERALLS

Designed for use in Class 7 cleanrooms and higher. These lightweight and breathable coveralls are designed for roominess and comfort. They feature zipper closures and elastic sleeve and ankle cuffs. Made from non-linting, non-woven 100% spun-bond polypropylene. Single use.

Number	Details	Quantity
A PS-4301	medium	25/box
PS-4401	large	25/box
PS-4501	x-large	25/box
PS-4601	xx-large	25/box
PS-4701	xxx-large	25/box
PS-4801	xxxx-large	25/box

BOUFFANT CAPS

Made from non-linting, non-woven 100% spun-bond polypropylene, these lightweight, breathable caps feature unexposed elastic bands for secure fit and comfort. Single use.

B PS-12101	21" (53cm)	1000/box
PS-12401	24" (61cm)	1000/box

BEARD COVERS

Designed for cool comfort while providing full coverage for cleanroom contamination control. non-woven 100% spun-bond polypropylene, these lightweight, breathable covers feature unexposed elastic bands for secure fit and comfort. Single use.

C PS-21801	one size fits all	1000/box
-------------------	-------------------	----------

CLEANROOM SHOE COVERS

Spun-bond polypropylene with stitched seams and anti-skid treads for safety. Single use.

D PS-6502	x-large	300/box
------------------	---------	---------

ADHESIVE MATS

Designed to effectively capture dirt and dust particles from foot-traffic and equipment wheels before they enter a controlled environment. Each numbered layer of the mat peels off easily to reveal a new clean surface. Acrylic adhesive layers in the mat do not support microbial growth.

E CR-32645-30	26" x 45"	120/box (30 mats/pad, 4 pads/box)
----------------------	-----------	-----------------------------------

PACKAGING AND PREPARATION

Mesh Trays hold smaller items that are difficult to keep track of during cleaning and sterilization as well as larger items with small parts. Our Mesh Trays, Instrument Trays, Pans, and Covers are heavy-duty stainless steel, for a tough, clean, and corrosion-resistant surface.

Surgical Tray Corner Protectors can be used with various trays and containers to help prevent wrap tears and wet packs.

Instrument Stringers are available in a wide selection of sizes and closure types, including Ball & Socket, Spring Lock, Pull-to-Open, Center Lock, U-Shaped, T-Bar, and Stringer Stand. Our Instrument Stringer Stand allows for vertical placement of instruments. No matter what your closure or style preference is, our stainless steel stringers are durable and versatile.

Sterilization Wrap is designed for use with EtO and steam sterilization. Highly permeable and durable, the wrap has a consistent texture giving you visual detection of small pinholes and tears.

Need Tray Liners, Peel Pouch Sterilization Racks, Band Clamps, or a Silicone Finger Mat (that you can cut to size)? We have those too.

Instrument Stringers

Band Clamps

Stringer Assistant

Tray Liners

Surgical Tray Corner Protectors

Wipers

Disinfectant

Lint Removers

Sterilization Pouches

Peel Pouch Sterilization Racks

Sterilization Wrap

Bowie-Dick Tests

Tape Dispensers

Silicone Finger Mat

Plastic Sterilization Trays

Mesh Trays

Scope Baskets

Perforated Mayo Trays

Oblong Instrument Trays




High-Sided Perforated Trays

Instrument Trays and Flat Covers

Catheter Trays and Dome Covers

Please note: Unless specified, all products are not made with natural rubber latex.





BALL & SOCKET


Highest quality German stainless steel stringers with stay-closed ball & socket design to prevent accidental opening. Expert hand-finishing eliminates sharp edges. To open or close, squeeze side bars together, then move cross bar.

Number	Length inches	Width inches	Length cm	Width cm	Quantity
A IS-12504	4"	2.5"	10cm	6cm	1/pkg
IS-12506	6"	2.5"	15cm	6cm	1/pkg
IS-12508	8"	2.5"	20cm	6cm	1/pkg
IS-12510	10"	2.5"	25cm	6cm	1/pkg

PULL-TO-OPEN

Stainless steel. Durable hinged design. Easy to pull open and close.

B IS-32504	4"	2.75"	10cm	7cm	1/pkg
IS-32506	6"	2.75"	15cm	7cm	1/pkg
IS-32508	8"	2.75"	20cm	7cm	1/pkg
IS-32510	10"	2.75"	25cm	7cm	1/pkg
IS-32512	12"	2.75"	30cm	7cm	1/pkg
IS-32514	14"	2.75"	36cm	7cm	1/pkg
IS-32516	16"	2.75"	41cm	7cm	1/pkg
IS-33018	18"	3"	46cm	8cm	1/pkg
IS-33020	20"	3"	51cm	8cm	1/pkg
C IS-35506	6"	5.5"	15cm	14cm	1/pkg
IS-35508	8"	5.5"	20cm	14cm	1/pkg
IS-35510	10"	5.5"	25cm	14cm	1/pkg
IS-35512	12"	5.5"	30cm	14cm	1/pkg



SPRING LOCK

Stainless steel. Expert hand-finishing eliminates sharp edges. Tapered spring design enables swift opening and closing.

Number	Length inches	Width inches	Length cm	Width cm	Quantity
A IS-23006	6"	2.75"	15cm	7cm	1/pkg
IS-23008	8"	2.75"	20cm	7cm	1/pkg
IS-23010	10"	2.75"	25cm	7cm	1/pkg
IS-23012	12"	2.75"	30cm	7cm	1/pkg
IS-23014	14"	2.75"	36cm	7cm	1/pkg
IS-23016	16"	2.75"	41cm	7cm	1/pkg

CENTER LOCK

Stainless steel. Expert hand-finishing eliminates sharp edges. To open or close, simply release center lock mechanism by unhooking center hinge.

B IS-42506	6"	3"	15cm	8cm	1/pkg
IS-42508	8"	3"	20cm	8cm	1/pkg
IS-42510	10"	3"	25cm	8cm	1/pkg
IS-42512	12"	3"	30cm	8cm	1/pkg
IS-42514	14"	3"	36cm	8cm	1/pkg
IS-42516	16"	3"	41cm	8cm	1/pkg

U-SHAPED

Stainless steel. Expert hand-finishing eliminates sharp edges. Designed for quick instrument handling. May be used individually or in pairs.

C IS-527506	6"	2.75"	15cm	7cm	1/pkg
IS-527508	8"	2.75"	20cm	7cm	1/pkg
IS-527510	10"	2.75"	25cm	7cm	1/pkg
IS-527512	12"	2.75"	30cm	7cm	1/pkg
IS-527514	14"	2.75"	36cm	7cm	1/pkg
IS-527516	16"	2.75"	41cm	7cm	1/pkg

U-SHAPED STRINGER END CAP

Stainless steel. Designed to fit over the end of U-Shaped stringers (2.75").


D IS-527-CAP	3.5"	3"	9cm	7cm	1/pkg
--------------	------	----	-----	-----	-------

U-SHAPED

Stainless steel. Expert hand-finishing eliminates sharp edges. Designed for quick instrument handling. Extended widths open instrument hinges further for cleaning. May be used individually or in pairs.

E IS-55506	6"	5.5"	15cm	14cm	1/pkg
IS-55508	8"	5.5"	20cm	14cm	1/pkg
IS-55510	10"	5.5"	25cm	14cm	1/pkg
IS-55512	12"	5.5"	30cm	14cm	1/pkg

Packaging and Preparation – Instrument Stringers, Stringer Assistant, and Band Clamps



T-BAR STRINGERS

Stainless steel. Expert hand-finishing eliminates sharp edges. Single bar stringer with a t-shaped stop at one end and a swiveling mechanism on the other that opens and locks shut. Designed for quick handling of instruments. May be used individually or in pairs.

Number	Length inches	Length cm	Quantity
A IS-7006	6"	15cm	1/pkg
B IS-7008	8"	20cm	1/pkg
IS-7012	12"	30cm	1/pkg

INSTRUMENT STRINGER STAND

Stainless steel. Designed for vertical placement of instruments within a tray.

Number	Base Measurements	Height	Quantity
C IS-62455*	2" x 4" x 5.25"	3"	1/pkg

STRINGER ASSISTANT

Designed to help with efficient instrument stringing. Features a tapered slot to firmly hold instrument stringers in a vertical position and an anti-skid bottom to keep it in place during use.

Number	Details	Dimensions	Quantity
D IS-675105*	white, acrylic, 7" tapered slot	10.25" x 7.5"	1/pkg

BAND CLAMPS

Instrument jaws need to be in an open position, and kept that way, during cleaning and sterilization. Band Clamps are a helpful tool for this, when an instrument stringer is not an option.

Number	Color	Width	Quantity
E IS-P0101	white	1.5"	1/pkg
F IS-P0105	yellow	2"	1/pkg
G IS-P0102	blue	2.25"	1/pkg
H IS-P0104	red	2.75"	1/pkg



WARNING: This product can expose you to chemicals including tetrafluoroethylene, which is known in the State of California to cause cancer. For more information go to www.p65warnings.ca.gov.



Packaging and Preparation – Corner Protectors and Tray Liners



SURGICAL TRAY CORNER PROTECTORS – Single Use

Designed to be used on the corners of sterilization trays to prevent rips and tears of sterilization wrap. Compatible with steam sterilization. Single use.

Number	Details	Dimensions	Quantity
A CS-9125C	foam, blue	6" x 6" x 1.25"	80 sheets/pkg (4 pre-cut protectors per sheet)

SURGICAL TRAY CORNER PROTECTORS – Reusable

Designed to be used on the corners of sterilization trays to prevent rips and tears of sterilization wrap. Compatible with steam and gas plasma sterilization. Reusable.

B CS-9113	medical-grade silicone, lime green	2.8" x 2.8" x 1.2"	80/pkg
-----------	------------------------------------	--------------------	--------



TRAY LINERS

Designed to help prevent residual moisture in your trays and reduce the incidents of wet packs.

Number	Dimensions	Details	Quantity
C TLP-012-018	12" x 18" (305mm x 457mm)	polyester/cellulose, 60gsm	600/box
TLP-020-025	20" x 25" (508mm x 635mm)	polyester/cellulose, 60gsm	400/box
D TL-202-500	13.78" x 22.83" (350mm x 580mm)	cellulose/crepe, folding sides, 107gsm	500/box
TL-204-500	13.78" x 13.78" (350mm x 350mm)	cellulose/crepe, folding sides, 107gsm	500/box
E TL-203-500	15.75" x 23.62" (400mm x 600mm)	cellulose/crepe, folding sides, 107gsm	500/box



Packaging and Preparation – Disinfectant, Lint Removers, and Wipers



A



B

SUPER SANI-CLOTH®

A nonwoven disposable cloth for use on hard, nonporous surfaces. Cleans, disinfects, and deodorizes. 55% alcohol content. Effective in 2 minutes or less against 22 microorganisms.

Number	Dimensions	Quantity
A DI-1160	6" x 6.75"	160 wipes/can (1 can/pkg)
DI-1160-12		12 cans/case

WALL MOUNT SANI-BRACKET

Adjustable Wall Mount Sani-Bracket for Large and Extra Large Sani-Cloth Wipes Canisters.

B DI-1000		1/pkg
-----------	--	-------

LINT REMOVER HANDLE AND LINT REMOVER REFILLS

Designed for heavy-duty hospital use.

C MJ-012	60 sheets/roll	12 rolls/box
D MJ-001	9" overall length	1 handle

POLYESTER WIPER

Designed for cleanroom duties. The wipes have ultrasonically sealed edges and are packaged in a controlled environment.

E CR-20909	9" x 9"	150 wipes/pkg
CR-20909-10		10 pkgs/case


KIMBERLY-CLARK WYPALL™ - L40 WIPERS

Kimberly-Clark Wypall - L40 Wipers are a general-purpose wiper featuring Kimberly-Clark's unique base sheet technology. Durable ¼ folded wipes easily absorb liquids yet are soft enough to use on delicate surfaces. Sold by the case.

F CR-21214	12.5" x 14.4"	56 wipes/pkg (18 pkgs/case)
------------	---------------	--------------------------------



C



E




D



F

Packaging and Preparation – Peel Pouch Sterilization Racks, Sterilization Pouches, and Sterilization Wrap



PEEL POUCH STERILIZATION RACKS

Designed to help keep peel pouches on edge while providing enough space for effective sterilant contact. These stainless steel racks feature slots for peel pouched items, aiding in efficient assembly and drying after steam sterilization.

Number	Dimensions	Quantity
A PR-090504-5	8.7" x 4.7" x 4", 5 slots, 1.5" apart	1/pkg
B PR-210806-4	20.5" x 7.5" x 6", 4 slots, 5" apart	1/pkg
C PR-150504	15.5" x 5" x 4", 10 slots, 1.5" apart	1/pkg
D PR-210806-8	20.5" x 8" x 6", 8 slots, 2.5" apart	1/pkg
E PR-160604-21	16.3" x 5.5" x 4", 21 slots, 0.75" apart	1/pkg

STERILIZATION POUCHES

Self-sealing sterilization pouches are made with both internal and external steam and EtO sterilization indicators. Sterilization indicators give the assurance of time, temperature, and pressure.

Number	Dimensions	Quantity
F CS-3553P	3.5" x 5.25"	200/box
CS-3590P	3.5" x 9"	200/box
CS-5310P	5.25" x 10"	200/box
CS-7513P	7.5" x 13"	200/box
CS-1218P	12" x 18"	100/box

STERILIZATION WRAP

Designed for use with ethylene oxide (EtO) and steam sterilization methods, it has the highest bacterial filtration efficiency available. Highly permeable to prevent wet packs. Bacterial barrier to both airborne and waterborne contamination.

Number	Dimensions	Quantity
G CS-1818	18" x 18"	500/pkg
CS-2020	20" x 20"	500/pkg
CS-2424	24" x 24"	500/pkg
CS-3030	30" x 30"	250/pkg
CS-3636	36" x 36"	250/pkg
CS-4545	45" x 45"	250/pkg
CS-5454	54" x 54"	100/pkg



Packaging and Preparation – Bowie-Dick Tests and Tape Dispensers



A



BOWIE-DICK TEST PACK – Reusable

Each test pack comes with 25 consecutively numbered test sheets, a Cube (top and bottom) and instructions for use. The metal clamp (sold separately) holds the indicator sheet in place and allows for the Lantor Cube to be placed directly over the drain. Certified to meet AAMI standards for daily use to evaluate the efficacy of air removal and steam penetration in pre-vacuum steam sterilizers. Reduces waste and inventory.

Number	Details	Quantity
A CS-BD001	Lantor Cube® Metal Clamp	1/pkg
CS-BD025	Lantor Cube® Bowie-Dick Test Pack	25 tests/box

BOWIE-DICK TEST PACK – Single Use

Bowie-Dick Test Pack has instructions for use pre-printed on the box for ease of use and staff training. Test sheet is easy to read and interpret. Certified to meet AAMI standards for daily air removal and steam penetration in pre-vacuum steam sterilizers.

Number	Details	Quantity
B CS-BD030	Airview™ Bowie-Dick Test Pack	30 packs/case



B



TAPE DISPENSER

Heavy-duty metal frame with slip-resistant base. 3" core. Designed to hold and dispense Process Indicator Tape. Cuts tape quickly and evenly. Safety blade cover for user safety.

Number	Details	Quantity
C CS-DISP-10	holds single roll of 1" wide tape	1/pkg
D CS-DISP-11	holds two rolls of 1" wide tape or one roll of 2" wide tape	1/pkg




C



D



Packaging and Preparation – Plastic Sterilization Trays and Silicone Finger Mat



PLASTIC STERILIZATION TRAYS


Designed for use with smaller or delicate instruments during sterilization and storage. Medical-grade silicone finger mats help protect and organize instruments within the tray. Perforated, slotted base and lid assist with circulation of steam during sterilization. Raised nodes in base assist with easy removal of silicone mat when needed. Transparent lid allows for easy viewing of contents. Compatible with steam and gas plasma sterilization.


Number	Details	Length	Width	Height	Quantity
A	PST-1100SB tray base, silicone mat, lid	7.5"	2.5"	.75"	1/pkg
B	PST-2100SB tray base, silicone mat, lid	6.5"	4"	.75"	1/pkg
C	PST-3100SB tray base, silicone mat, lid	10"	6"	.75"	1/pkg
D	PST-3200SB tray base, silicone mat, lid	10"	6"	1.5"	1/pkg
E	PST-3300SB tray base, insert tray, 2 silicone mats, lid	10"	6"	1.5"	1/pkg

SILICONE FINGER MAT

Medical-grade silicone mat with raised fingers is designed to hold, organize, and protect surgical instruments within a tray during sterilization. This versatile product can be cut to size to fit your specific needs. Can be used with Mesh Trays or as mat replacements for Plastic Sterilization Trays. Compatible with steam and gas plasma sterilization.

Number	Details	Length	Width	Height	Quantity
F	FM-1810 silicone finger mat	18.2"	10.2"	.75"	1/pkg






MICRO MESH TRAYS

Micro Mesh Trays are designed to hold smaller items during washing, sterilization, storage, and transportation. Micro mesh allows you to keep larger items that have small loose parts together. Micro mesh size is approximately 2mm x 2mm. Compatible with steam and gas plasma sterilization.


Number	Details	Length	Width	Height	Quantity
A	MT-5500 with lid	1.57"/40mm	1.57"/40mm	.79"/20mm	1/pkg
B	MT-5600 with lid	3.15"/80mm	1.57"/40mm	.79"/20mm	1/pkg
C	MT-5650 with lid	6.30"/160mm	1.57"/40mm	.79"/20mm	1/pkg
D	MT-5700 with lid	2.75"/70mm	2.75"/70mm	1.18"/30mm	1/pkg
E	MT-5800 with lid	3.15"/80mm	3.15"/80mm	1.57"/40mm	1/pkg
F	MT-5900 with lid	3.54"/90mm	3.54"/90mm	1.97"/50mm	1/pkg
G	MT-6020 with lid	7.09"/180mm	5.12"/130mm	.79"/20mm	1/pkg
H	MT-7100 with lid	8.89"/225mm	4.33"/110mm	1.97"/50mm	1/pkg
I	MT-7200 with lid	10.63"/270mm	5.31"/135mm	1.57"/40mm	1/pkg
J	MT-7300 with lid	15.75"/400mm	9.84"/250mm	1.97"/50mm	1/pkg
K	MT-9001 with drop handles	18.90"/480mm	9.84"/250mm	1.97"/50mm	1/pkg
L	MT-9101 with drop handles	9.84"/250mm	9.45"/240mm	1.97"/50mm	1/pkg
M	MT-9521 with drop handles	27.56"/700mm	8.27"/210mm	1.97"/50mm	1/pkg




MESH TRAYS

Mesh Trays are designed for a wide variety of instruments during washing, sterilization, storage, and transportation. Mesh size is approximately 5mm x 5mm. Compatible with steam and gas plasma sterilization.

Number	Details	Length	Width	Height	Quantity
A	MT-8005 with drop handles	10.24"/260mm	7.48"/190mm	1.97"/50mm	1/pkg
B	MT-8050 with drop handles	14.96"/380mm	9.84"/250mm	1.97"/50mm	1/pkg
C	MT-8500 with drop handles	19.69"/500mm	5.91"/150mm	3.94"/100mm	1/pkg
D	MT-9000 with drop handles	18.90"/480mm	9.84"/250mm	1.97"/50mm	1/pkg
E	MT-9010 with drop handles	18.90"/480mm	9.84"/250mm	2.75"/70mm	1/pkg
F	MT-9100 with drop handles	9.84"/250mm	9.84"/250mm	1.97"/50mm	1/pkg
G	MT-9110 with drop handles	9.84"/250mm	9.45"/240mm	2.75"/70mm	1/pkg
H	MT-9210 with drop handles	17.32"/440mm	4.33"/110mm	2.56"/65mm	1/pkg





SCOPE BASKETS

Stainless steel wire scope baskets designed to contain and protect rigid scopes during the washing and sterilization process. Fully removable lid with easy-to-use locking mechanism. The fixed silicon padding in the base and cover gently but firmly hold scopes in place. Compatible with steam, gas plasma, and EtO sterilization.

Number	Details	Length	Width	Height	Quantity
A SB-1210	stainless steel with silicon inserts	7.08"/180mm	3.15"/80mm	2.05"/52mm	1/pkg
B SB-4210	stainless steel with silicon inserts	18.11"/460mm	3.15"/80mm	2.05"/52mm	1/pkg
C SB-4610	stainless steel with silicon inserts	18.11"/460mm	6.30"/160mm	2.05"/52mm	1/pkg
D SB-6520	stainless steel with silicon inserts -includes small lens basket with fixed bracket	25.19"/640mm 3.15"/80mm	5.90"/150mm 3.15"/80mm	3.03"/77mm 1.57"/50mm	1/pkg



Packaging and Preparation – Perforated Mayo Trays, Oblong Instrument Trays, and High-Sided Perforated Trays



A B C D E



M

L



F G H I J K



N O



P Q R

PERFORATED MAYO TRAYS

Coved corners and large, open bead allow for easy handling, stacking, and cleaning. Perforations allow for drainage. 304 stainless steel. Finished to provide a durable and corrosion resistant surface.

Number	Length	Width	Height	Quantity
A PW-10F-PMT	10"/25cm	6.5"/17cm	.75"/2cm	1/pkg
B PW-13F-PMT	13"/33cm	9"/23cm	.75"/2cm	1/pkg
C PW-15F-PMT	15"/38cm	10"/25cm	.75"/2cm	1/pkg
D PW-17F-PMT	17"/43cm	11"/28cm	.75"/2cm	1/pkg
E PW-19F-PMT	19"/48cm	12"/30cm	.75"/2cm	1/pkg

OBLONG INSTRUMENT TRAYS

Coved corners and large, open bead allow for easy handling, stacking, and cleaning. 304 stainless steel. Finished to provide a durable and corrosion resistant surface.

F PW-10F	10"/25cm	6.5"/17cm	.75"/2cm	1/pkg
PW-10F-10	10"/25cm	6.5"/17cm	.75"/2cm	10/pkg
G PW-13F	13"/33cm	9"/23cm	.75"/2cm	1/pkg
PW-13F-10	13"/33cm	9"/23cm	.75"/2cm	10/pkg
H PW-15F	15"/38cm	10"/25cm	.75"/2cm	1/pkg
I PW-17F	17"/43cm	11"/28cm	.75"/2cm	1/pkg
J PW-19F	19"/48cm	12"/30cm	.75"/2cm	1/pkg
K PW-1520	21"/53cm	16"/41cm	.75"/2cm	1/pkg

HIGH-SIDED PERFORATED TRAYS

High sides protect fragile instruments during set-up and procedures. Perforations allow for drainage. 304 stainless steel. Finished to provide a durable and corrosion resistant surface.

L PW-952-HPT	8"/20cm	5"/13cm	2"/5cm	1/pkg
M PW-952-2 (cover for PW-952-HPT)	9.16"/29.3cm	5"/13cm	.40"/1cm	1/pkg
N PW-1002-HPT	10"/25cm	6.50"/17cm	2"/5cm	1/pkg
O PW-1202-HPT	12.25"/31cm	7.75"/20cm	2.25"/6cm	1/pkg
P PW-S12102-HPT	12"/30cm	10"/25cm	2.50"/6cm	1/pkg
Q PW-E1650-HPT	16.50"/42cm	10"/25cm	2.50"/6cm	1/pkg
R PW-S20122-HPT	20.75"/53cm	12.75"/32cm	2.50"/6cm	1/pkg

*TRAYS AND COVERS SOLD SEPARATELY

Packaging and Preparation – Instrument Trays and Flat Covers



INSTRUMENT TRAYS

Seamless trays have covered corners for easy handling and transportation. Finished to provide a durable and corrosion resistant surface. 304 stainless steel.

Number	Length	Width	Height	Quantity
A PW-952	8.94"/22.7cm	5"/13cm	2"/5cm	1/pkg
B PW-953	8.94"/22.7cm	5"/13cm	3"/8cm	1/pkg
C PW-1002	10"/25cm	6.50"/17cm	2"/5cm	1/pkg
D PW-1202-10	12.25"/31cm	7.75"/20cm	2.25"/6cm	10/pkg
E PW-E20124	20.75"/53cm	12.75"/32cm	4"/10cm	1/pkg
F PW-E20126	20.75"/53cm	12.75"/32cm	6"/15cm	1/pkg

FLAT COVERS

Covers have recessed strap handles for easy handling. 304 stainless steel. Finished to provide a durable and corrosion resistant surface.

Number	Details	Quantity
G PW-952-2	flat cover for PW-952 and PW-953	1/pkg
H PW-1002-2	flat cover for PW-1002 and PW-1002-HPT	1/pkg
I PW-1202-2	flat cover for PW-1202 and PW-1202-HPT	1/pkg
J PW-2000-2	flat cover for PW-E20124 and PW-E20126	1/pkg

*TRAYS AND COVERS SOLD SEPARATELY

Packaging and Preparation – Catheter Trays and Dome Covers



CATHETER TRAYS

Seamless trays. 304 stainless steel. Finished to provide a durable and corrosion resistant surface.

Number	Length	Width	Height	Quantity
A PW-83	8.75"/22cm	3.25"/8cm	1.50"/4cm	1/pkg
B PW-170-0	17.25"/44cm	4.25"/10.79cm	2"/5cm	1/pkg
C PW-952	8.94"/22.7cm	5"/13cm	2"/5cm	1/pkg
D PW-953	8.94"/22.7cm	5"/13cm	3"/8cm	1/pkg

DOMES COVERS

Dome covers feature strap handles for easy handling. 304 stainless steel. Finished to provide a durable and corrosion resistant surface.

Number	Details	Quantity
E PW-83-2	dome cover for PW-83	1/pkg
F PW-170-2	dome cover for PW-170-0	1/pkg
G PW-952D-2	dome cover for PW-952 and PW-953	1/pkg

*TRAYS AND COVERS SOLD SEPARATELY

WIRE AND PIN MANAGEMENT

Key Surgical provides the highest quality and broadest selection of wires and pins for a variety of procedures. To support your management of this specialized inventory, we also have you covered with products to sterilize and dispense, cut, protect, and remove wires and pins.

K-Wires and Steinmann Pins are available in various diameters and lengths. Our high-quality wires and pins are made from certified implantable stainless steel and are precision ground in Germany.

Wire Dispensers and Racks are designed to be used during sterilization of wires and pins as well as dispensing them in the operating room during a procedure. Whether you prefer single or multiple chamber dispensers or racks that can fold open to a tabletop position in the O.R., we've got what you need.

Cutting, protecting, and then removing wires and pins requires a variety of tools. We've got you there too.

[K-Wires and Steinmann Pins](#)

[K-Wire and Pin Covers](#)

[K-Wire Dispensers](#)

[Sterilization Racks](#)

[K-Wire Ruler and Pin Gauge](#)

[Wire and Pin Cutters](#)

[Pliers](#)

[Twist Drill Bits and Drill Rack](#)

Please note: Unless specified, all products are not made with natural rubber latex.



K-WIRES – DOUBLE TROCAR

Available in various lengths and diameters. Made from certified medical-grade stainless steel. Sold non-sterile. 6 wires per tube, 1 tube per package.

Number	Length inches	Length cm	Diameter inches	Diameter mm
KI-71-100	4"	10cm	.028"	.7mm
KI-71-101	4"	10cm	.035"	.9mm
KI-71-102	4"	10cm	.045"	1.1mm
KI-71-104	4"	10cm	.054"	1.4mm
KI-71-103	4"	10cm	.062"	1.6mm
KI-71-105	5"	13cm	.028"	.7mm
KI-71-106	5"	13cm	.035"	.9mm
KI-71-107	5"	13cm	.045"	1.1mm
KI-71-108	5"	13cm	.054"	1.4mm
KI-71-109	5"	13cm	.062"	1.6mm
KI-71-021	6"	15cm	.028"	.7mm
KI-71-022	6"	15cm	.035"	.9mm
KI-71-023	6"	15cm	.045"	1.1mm
KI-71-024	6"	15cm	.054"	1.4mm
KI-71-025	6"	15cm	.062"	1.6mm
KI-71-111	9"	23cm	.035"	.9mm
KI-71-112	9"	23cm	.045"	1.1mm
KI-71-114	9"	23cm	.054"	1.4mm
KI-71-113	9"	23cm	.062"	1.6mm
KI-71-121	12"	30cm	.035"	.9mm
KI-71-122	12"	30cm	.045"	1.1mm
KI-71-124	12"	30cm	.054"	1.4mm
KI-71-123	12"	30cm	.062"	1.6mm

ALL .062" K-WIRES CAN BE THREADED. PLEASE CALL FOR MORE DETAILS ON THREADED K-WIRES.



ADDITIONAL SIZES AND STYLES AVAILABLE UPON REQUEST

K-WIRES – SINGLE TROCAR WITH ROUND END

Available in various lengths and diameters. Made from certified medical-grade stainless steel. Sold non-sterile. 6 wires per tube, 1 tube per package.

Number	Length inches	Length cm	Diameter inches	Diameter mm
KI-71-130	4"	10cm	.028"	.7mm
KI-71-131	4"	10cm	.035"	.9mm
KI-71-132	4"	10cm	.045"	1.1mm
KI-71-133	4"	10cm	.054"	1.4mm
KI-71-134	4"	10cm	.062"	1.6mm
KI-71-170	5"	13cm	.028"	.7mm
KI-71-171	5"	13cm	.035"	.9mm
KI-71-172	5"	13cm	.045"	1.1mm
KI-71-174	5"	13cm	.054"	1.4mm
KI-71-175	5"	13cm	.062"	1.6mm
KI-71-176	6"	15cm	.028"	.7mm
KI-71-177	6"	15cm	.035"	.9mm
KI-71-182	6"	15cm	.045"	1.1mm
KI-71-184	6"	15cm	.054"	1.4mm
KI-71-185	6"	15cm	.062"	1.6mm
KI-71-151	9"	23cm	.035"	.9mm
KI-71-152	9"	23cm	.045"	1.1mm
KI-71-154	9"	23cm	.054"	1.4mm
KI-71-153	9"	23cm	.062"	1.6mm
KI-71-161	12"	30cm	.035"	.9mm
KI-71-162	12"	30cm	.045"	1.1mm
KI-71-164	12"	30cm	.054"	1.4mm
KI-71-163	12"	30cm	.062"	1.6mm

ALL .062" K-WIRES CAN BE THREADED. PLEASE CALL FOR MORE DETAILS ON THREADED K-WIRES.



K-WIRES – DOUBLE DIAMOND

Available in various lengths and diameters. Made from certified medical-grade stainless steel. Sold non-sterile. 6 wires per tube, 1 tube per package.

Number	Length inches	Length cm	Diameter inches	Diameter mm
KI-71-200	4"	10cm	.028"	.7mm
KI-71-201	4"	10cm	.035"	.9mm
KI-71-202	4"	10cm	.045"	1.1mm
KI-71-203	4"	10cm	.054"	1.4mm
KI-71-204	4"	10cm	.062"	1.6mm
KI-71-205	5"	13cm	.028"	.7mm
KI-71-206	5"	13cm	.035"	.9mm
KI-71-207	5"	13cm	.045"	1.1mm
KI-71-208	5"	13cm	.054"	1.4mm
KI-71-209	5"	13cm	.062"	1.6mm
KI-71-030	6"	15cm	.028"	.7mm
KI-71-031	6"	15cm	.035"	.9mm
KI-71-032	6"	15cm	.045"	1.1mm
KI-71-034	6"	15cm	.054"	1.4mm
KI-71-033	6"	15cm	.062"	1.6mm
KI-71-211	9"	23cm	.035"	.9mm
KI-71-212	9"	23cm	.045"	1.1mm
KI-71-214	9"	23cm	.054"	1.4mm
KI-71-213	9"	23cm	.062"	1.6mm



K-WIRES – SINGLE DIAMOND

Available in various lengths and diameters. Made from certified medical-grade stainless steel. Sold non-sterile. 6 wires per tube, 1 tube per package.

KI-71-230	4"	10cm	.028"	.7mm
KI-71-241	9"	23cm	.035"	.9mm
KI-71-242	9"	23cm	.045"	1.1mm
KI-71-243	9"	23cm	.062"	1.6mm
KI-71-251	12"	30cm	.035"	.9mm
KI-71-252	12"	30cm	.045"	1.1mm
KI-71-253	12"	30cm	.062"	1.6mm

ALL .062" K-WIRES CAN BE THREADED. PLEASE CALL FOR MORE DETAILS ON THREADED K-WIRES.



STEINMANN PINS

Available in 9" (23cm) length with various diameters. Made from certified medical-grade stainless steel. Sold non-sterile. 6 pins per tube, 1 tube per package.

A Double Trocar

Number	Diameter inches	Diameter mm
KI-71-450	.078"	2.0mm
KI-71-453	.093"	2.4mm
KI-71-456	.109"	2.8mm
KI-71-459	.125"	3.2mm
KI-71-462	.141"	3.5mm
KI-71-465	.156"	4.0mm
KI-71-468	.177"	4.5mm

B Single Trocar with Round End

Number	Diameter inches	Diameter mm
KI-71-510	.078"	2.0mm
KI-71-513	.093"	2.4mm
KI-71-516	.109"	2.8mm
KI-71-519	.125"	3.2mm
KI-71-522	.141"	3.5mm
KI-71-525	.156"	4.0mm
KI-71-528	.177"	4.5mm

C Single Trocar with 3 Shank End

Number	Diameter inches	Diameter mm
KI-71-480	.078"	2.0mm
KI-71-483	.093"	2.4mm
KI-71-486	.109"	2.8mm
KI-71-489	.125"	3.2mm
KI-71-492	.141"	3.5mm
KI-71-495	.156"	4.0mm
KI-71-498	.177"	4.5mm



STEINMANN PINS

Available in 9" (23cm) length with various diameters. Made from certified medical-grade stainless steel. Sold non-sterile. 6 pins per tube, 1 tube per package.

A Double Diamond


Number	Diameter inches	Diameter mm
KI-71-600	.078"	2.0mm
KI-71-603	.093"	2.4mm
KI-71-606	.109"	2.8mm
KI-71-609	.125"	3.2mm
KI-71-612	.141"	3.5mm
KI-71-615	.156"	4.0mm
KI-71-618	.177"	4.5mm

B Single Diamond with Round End

Number	Diameter inches	Diameter mm
KI-71-660	.078"	2.0mm
KI-71-663	.093"	2.4mm
KI-71-666	.109"	2.8mm
KI-71-669	.125"	3.2mm
KI-71-672	.141"	3.5mm
KI-71-675	.156"	4.0mm
KI-71-678	.177"	4.5mm

C Single Diamond with 3 Shank End

Number	Diameter inches	Diameter mm
KI-71-639	.125"	3.2mm
KI-71-642	.141"	3.5mm
KI-71-645	.156"	4.0mm
KI-71-648	.177"	4.5mm



STEINMANN PINS – THREADED

Available in 9" (23cm) length with various diameters. Made from certified medical-grade stainless steel. Sold non-sterile. 6 pins per tube, 1 tube per package.

A Double Trocar – Threaded

Number	Diameter inches	Diameter mm
KI-71-750	.078"	2.0mm
KI-71-753	.093"	2.4mm
KI-71-756	.109"	2.8mm
KI-71-759	.125"	3.2mm
KI-71-762	.141"	3.5mm
KI-71-765	.156"	4.0mm
KI-71-768	.177"	4.5mm

B Single Trocar With Round End – Threaded

Number	Diameter inches	Diameter mm
KI-71-770	.078"	2.0mm
KI-71-773	.093"	2.4mm
KI-71-776	.109"	2.8mm
KI-71-779	.125"	3.2mm
KI-71-782	.141"	3.5mm
KI-71-785	.156"	4.0mm
KI-71-788	.177"	4.5mm

C Double Diamond – Threaded

Number	Diameter inches	Diameter mm
KI-71-730	.078"	2.0mm
KI-71-733	.093"	2.4mm
KI-71-736	.109"	2.8mm
KI-71-739	.125"	3.2mm
KI-71-742	.141"	3.5mm
KI-71-745	.156"	4.0mm
KI-71-748	.177"	4.5mm

D Single Diamond With Round End – Threaded

Number	Diameter inches	Diameter mm
KI-71-700	.078"	2.0mm
KI-71-703	.093"	2.4mm
KI-71-706	.109"	2.8mm
KI-71-709	.125"	3.2mm
KI-71-712	.141"	3.5mm
KI-71-715	.156"	4.0mm
KI-71-718	.177"	4.5mm



PINS SHOWN NON-THREADED



SINGLE-CHAMBER K-WIRE DISPENSERS: GERMAN STAINLESS STEEL

Store, sterilize, and dispense one size of K-wire. Dispenses one wire at a time. High quality German stainless steel. Can hold lengths up to 12". Perforated for proper sterilization. Clearly marked with diameter for quick identification. 1 dispenser/pkg.

	Number	Length inches	Length cm	Diameter inches	Diameter mm	Quantity Held
A	KI-71-001	up to 6"	15cm	.028"	.7mm	100
B	KI-71-002	up to 6"	15cm	.035"	.9mm	50
C	KI-71-003	up to 6"	15cm	.045"	1.1mm	40
D	KI-71-004	up to 6"	15cm	.054"	1.4mm	25
E	KI-71-006	up to 6"	15cm	.062"	1.6mm	20
F	KI-71-012	up to 12"	30cm	.035"	.9mm	50
G	KI-71-013	up to 12"	30cm	.045"	1.1mm	40
H	KI-71-014	up to 12"	30cm	.054"	1.4mm	25
I	KI-71-016	up to 12"	30cm	.062"	1.6mm	20
J	KI-71-017	up to 12"	30cm	.078"	2.0mm	10
K	KI-71-019	up to 12"	30cm	.093"	2.4mm	5

MULTI-CHAMBER K-WIRE DISPENSERS: GERMAN STAINLESS STEEL

Store, sterilize, and dispense four different sizes of K-wires. Dial on bottom allows for easy dispensing of one size at a time. Chambers are clearly marked for easy identification. High quality German stainless steel. 1 dispenser/pkg.

	Number	Length	Diameter	Quantity Held	Diameter	Quantity Held	Diameter	Quantity Held	Diameter	Quantity Held
L	KI-46-542	4" (10cm)	.035" (.9mm)	30	.045" (1.1mm)	24	.054" (1.4mm)	24	.062" (1.6mm)	12
M	KI-46-541	6" (15cm)	.035" (.9mm)	30	.045" (1.1mm)	24	.054" (1.4mm)	24	.062" (1.6mm)	12
N	KI-46-538	12" (30cm)	.035" (.9mm)	30	.045" (1.1mm)	24	.054" (1.4mm)	24	.062" (1.6mm)	12
O	KI-46-539	12" (30cm)	.070" (1.8mm)	18	.078" (2.0mm)	12	.086" (2.2mm)	12	.094" (2.5mm)	6


MULTI-CHAMBER K-WIRE DISPENSERS: O.R. QUALITY STAINLESS STEEL

Store, sterilize, and dispense four different sizes of K-wires. Dial on bottom allows for easy dispensing of one size at a time. Chambers are clearly marked for easy identification. O.R. Quality stainless steel. 1 dispenser/pkg.

	Number	Length	Diameter	Quantity Held	Diameter	Quantity Held	Diameter	Quantity Held	Diameter	Quantity Held
P	KE-1210	4" (10cm)	.035" (.9mm)	30	.045" (1.1mm)	24	.054" (1.4mm)	24	.062" (1.6mm)	12
Q	KE-1220	12" (30cm)	.035" (.9mm)	30	.045" (1.1mm)	24	.054" (1.4mm)	24	.062" (1.6mm)	12

NOTE: DISPENSERS NOT RECOMMENDED FOR USE WITH THREADED WIRES. STERILIZATION RACKS SUGGESTED.

Wire and Pin Management – Sterilization Racks & Pin Gauge



K-WIRES AND PINS SOLD SEPARATELY

STERILIZATION RACKS

Sterilization racks are designed for use during washing, sterilizing, storing/transporting, and dispensing K-wires and Steinmann Pins. The racks are lightweight but made to withstand repeated use. They fold closed for sterilization and easily open to convert to a table top stand for use in the O.R. Dimensions are clearly marked on the rack to assist with quick identification.

K-Wire and Pin Rack

Number	Details	Dimensions	Quantity
A KI-71-RACK	anodized aluminum, holds 9" K-wires and pins Holds 6 of each of the following diameters:	12.5" x 5.25" x 1.5" (32 x 13 x 3.8cm)	1/pkg
	.028" (.7mm)	.045" (1.1mm)	.062" (1.6mm)
	.093" (2.4mm)	.125" (3.2mm)	.156" (4.0mm)
	.035" (.9mm)	.054" (1.4mm)	.078" (2.0mm)
	109" (2.8mm)	.140" (3.5mm)	.177" (4.5mm)


K-Wire Rack

B KI-71-RACK-46	anodized aluminum, holds 4-6" K-wires Holds 6 of each of the following diameters:	6.25" x 4.125" x .75" (15.9 x 10.5 x 1.9cm)	1/pkg
	.028" (.7mm)	.045" (1.1mm)	.062" (1.6mm)
	.035" (.9mm)	.054" (1.4mm)	

K-WIRE RULER AND PIN GAUGE

Measure 16 sizes of K-wire and pin lengths in inches, millimeters, and centimeters as well as diameters. Made from 420-series stainless steel with .025" thickness for durability.

C KI-39-710	6.75" (17.1cm) ruler	1/pkg
-------------	----------------------	-------






K-WIRE AND PIN COVERS

Post-surgical covers for protruding wires and pins. Plastic sheaths act as a guide for easy use. Available in twelve different sizes. Sold sterile or non-sterile. NS designates non-sterile. Compatible with steam and gas plasma sterilization. Single use.

Number	Diameter inches	Diameter mm	Color	Quantity
A C-028 C-028NS	.028"	.7mm	yellow	2 covers/pkg, 12 pkgs/box 50/pkg
B C-035 C-035NS	.035"	.9mm	blue	2 covers/pkg, 12 pkgs/box 50/pkg
C C-045 C-045NS	.045"	1.1mm	white	2 covers/pkg, 12 pkgs/box 50/pkg
D C-054 C-054NS	.054"	1.4mm	red	2 covers/pkg, 12 pkgs/box 50/pkg
E C-062 C-062NS	.062"	1.6mm	green	2 covers/pkg, 12 pkgs/box 50/pkg
F C-078 C-078NS	.078"	2.0mm	pink	2 covers/pkg, 12 pkgs/box 50/pkg
G C-093 C-093NS	.093"	2.4mm	gray	2 covers/pkg, 12 pkgs/box 50/pkg
H C-109 C-109NS	.109"	2.8mm	tan	2 covers/pkg, 12 pkgs/box 50/pkg
I C-125 C-125NS	.125"	3.2mm	lavender	2 covers/pkg, 12 pkgs/box 50/pkg
J C-141 C-141NS	.141"	3.5mm	brown	2 covers/pkg, 12 pkgs/box 50/pkg
K C-156 C-156NS	.156"	4.0mm	navy	2 covers/pkg, 12 pkgs/box 50/pkg
L C-188 C-188NS	.188"	4.5mm	orange	2 covers/pkg, 12 pkgs/box 50/pkg

PRODUCTS LARGER THAN ACTUAL SIZE



WIRE AND PIN CUTTERS


Expertly cuts everything from delicate wire to .25" pins. High quality German stainless steel.

Number	Details	Overall Length	Max. Diameter Cut
A	KI-46-622 delicate wire cutter, side cutting with tungsten carbide inserts	5"	.028"
B	KI-49-734 K-wire side cutter with tungsten carbide inserts	6"	.045"
C	KI-46-621 side cutter with tungsten carbide inserts	7"	.062"
D	KI-46-623 double action wire cutter, angled side cutting with tungsten carbide inserts	7"	.062"
E	KI-46-625 front and side wire cutter with tungsten carbide inserts	7"	.062"
F	KI-46-631 flush wire and cable cutter, front and side cutting with tungsten carbide inserts	7"	.062"
G	KI-46-647 cannulated pin cutter, end cutting	7.50"	.062"
H	KI-46-651 diamond pin cutter, end cutting, leaves flat-cut pins	6.50"	.078"
I	KI-46-624 double action wire cutter, angled side cutting with tungsten carbide inserts	8.86"	.093"
J	KI-49-736 side cutting with tungsten carbide inserts	9"	.094"
K	KI-46-626 large pin cutter, side cutting	22"	.25"

ADDITIONAL INSTRUMENTS AVAILABLE ON REQUEST



PRODUCTS NOT SHOWN AT RELATIVE SCALE



PLIERS: GERMAN STAINLESS STEEL

Bend, pull, and remove pins with durable, finely finished pliers. High quality German stainless steel.

Number	Details	Overall Length
A	needle nose pliers, delicate	5.51"
B	flat nose pliers, 5mm tip	5.50"
C	pin puller, with groove, tungsten carbide	5.50"
D	flat nose pliers, delicate	5.91"
E	flat nose pliers, cross-serrated surface	7"
F	wire extraction pliers, 2mm delicate jaw, tungsten carbide	7"
G	wire extraction pliers, 3mm delicate jaw, tungsten carbide	7"
H	narrow nose pliers, tapers to 3mm	7.50"
I	slip joint pliers, 10mm tip	7.87"
J	parallel pliers, 10mm tip	7"
K	parallel pliers, with cutter	7.50"
L	small locking pliers, tip width 7mm	6.75"
M	medium locking pliers, tip width 9mm	7.68"
N	large locking pliers, tip width 13mm	10"
O	needle nose locking pliers, jaw length 2.5", tip width 5mm	9.45"
P	large needle nose locking pliers, jaw length 3", tip width 7mm	12"



PRODUCTS NOT SHOWN AT RELATIVE SCALE
ADDITIONAL INSTRUMENTS AVAILABLE UPON REQUEST



O.R. Products – Twist Drill Rack and Twist Drill Bits



TWIST DRILL BITS

Jacobs Twist Drill Bits. German surgical stainless steel. Used to drill holes in bone for placement of bone screws during orthopedic surgery involving fixation of internal fractures. Available in a variety of diameters.

Number	Details	Quantity
A	KI-86-8415 5" str shank double fluted, 1.5mm diameter, 17mm, round end	1/pkg
B	KI-86-8420 5" str shank double fluted, 2.0mm diameter, 23mm, round end	1/pkg
C	KI-86-8425 5" str shank double fluted, 2.5mm diameter, 22mm, round end	1/pkg
D	KI-86-8427 5" str shank double fluted, 2.7mm diameter, 30mm, round end	1/pkg
E	KI-86-8432 5" str shank double fluted, 3.2mm diameter, 42mm, round end	1/pkg
F	KI-86-8435 5" str shank double fluted, 3.5mm diameter, 42mm, round end	1/pkg
G	KI-86-8440 5" str shank double fluted, 4.0mm diameter, 45mm, round end	1/pkg
H	KI-86-8445 5" str shank double fluted, 4.5mm diameter, 34mm, round end	1/pkg
I	KI-86-8448 5" str shank double fluted, 4.7mm diameter, 34mm, round end	1/pkg
J	KI-86-8450 5" str shank double fluted, 5.0mm diameter, 42mm, round end	1/pkg

TWIST DRILL RACK

O.R. Quality stainless steel sterilization rack designed to hold 13 drill bits of 5" length and 1.0 – 5.0mm diameter.

Number	Details	Dimensions	Quantity
K	KI-05-RACK-TD Holds the following diameters and quantities:	6.75" x 4.25" x .71" (17.1 x 10.5 x 1.8cm)	1/pkg
	1.5mm (2)	2.0mm (2)	2.5mm (2)
	2.7mm (1)	3.2mm (1)	3.5mm (1)
	4.0mm (1)	4.5mm (1)	4.7mm (1)
	5.0mm (1)		



TWIST DRILL BITS SOLD SEPARATELY

INSTRUMENT TRACKING

Key Surgical's KeyDot® is a highly effective and affordable solution for tracking individual surgical instruments beyond the tray level. KeyDot is an easy, and economical way to get the most from an asset tracking system and to maximize your return on investment. It can be used to track specialty, high-value, or peel-pack items on an individual basis. Whatever the application, KeyDot helps provide increased efficiency, accuracy, and productivity in tray assembly, general instrument identification, instrument inventory status, and management.

The KeyDot is a small, laser-etched, 2D Data Matrix barcode label that can be applied directly to an instrument's surface. Tested for sterilization efficacy and biocompatibility for use on stainless steel surgical instruments, it has been conveniently incorporated into most of the instrument tracking software applications available today.

With the KeyDot there is no specialized training or intricate equipment needed so application can be done on-site, eliminating the hassle of sending your instrument sets out for this tracking solution.

PLEASE NOTE: Unless specified, all products are not made with natural rubber latex.





HOW DOES KEYDOT WORK?

The KeyDot works like a license plate. It links a specific instrument to a data file within a software application, allowing that instrument to be tracked throughout the complete instrument processing cycle. Instruments marked with the KeyDot are scanned using a two-dimensional barcode reader. Scanners and software applications not available through Key Surgical.

THE KEYDOT IS EFFICIENT

KeyDot provides you with the convenience to decide what to mark and when to mark it – you're in control. Since there are no equipment requirements or complicated procedures and training involved, instruments can be marked with a KeyDot on a continuous basis by anyone in your department, starting with critical instruments and expanding at an individual pace.

KEYDOT IS EASY

With the KeyDot there is no specialized training or intricate equipment required. Applying a KeyDot to a surgical instrument is quick and easy: simply clean, lift, place, and press. The instrument then needs to be sterilized to cure the adhesive and assure firm adhesion.

1 clean application area



2 lift KeyDot from sheet



3 place KeyDot on instrument



4 apply firm pressure and sterilize to secure adhesive



FREQUENTLY ASKED QUESTIONS

How long will our KeyDot last?

KeyDot labels will remain readable for at least 200 processing cycles or two years, whichever comes first. The radial design of the KeyDot means there are no corners to lift up. The Data Matrix barcode is also highly redundant – even if up to 60% of the barcode is destroyed, the KeyDot is still readable. The KeyDot is designed to withstand harsh environmental and temperature conditions. If a KeyDot needs to be replaced, anyone can do it at a convenient time in a matter of seconds. And at around fifty cents per dot, the KeyDot remains an easy, efficient and economical way to mark instruments for the life of the instrument.

Is the KeyDot safe for our surgical instruments?

Yes. A licensed, independent medical testing laboratory has scientifically evaluated the KeyDot and found that the instrument surfaces under the KeyDot are completely sterilized during processing.

Do I have to program the KeyDots with my own numbers?

No. All KeyDots are manufactured and shipped with a unique number. We guarantee that there will be no duplicate numbers in the labeling sequence. Custom sequences can also be requested. Simply apply the KeyDot and easily link it to the appropriate data file in your software application.

How many characters are allowed on each KeyDot?

Human readable equivalents appear on the 1/8" and 3/16" KeyDot (up to 15 characters) and on the 1/4" and 3/8" (up to 17 characters).

KEYDOT GUARANTEE

Key Surgical guarantees that our KeyDot labels will remain readable for at least 200 processing cycles or two years, whichever comes first. Simply send us the supporting documentation from your instrument tracking system and we will replace, free of charge, any KeyDot that does not meet our performance benchmark.

KEYDOT

Alphanumeric barcodes are laser engraved on .003" destructible acrylate. 2D Data Matrix. Black background with white copy. Pressure-sensitive acrylic adhesive. Label copy includes barcode and human readable equivalents on the 1/4" and 3/8" KeyDot. 100 KeyDot labels per sheet.

	Number	Dimensions	
A	ID-405	1/8" diameter	
B	ID-400	3/16" diameter	
C	ID-410	1/4" diameter	
D	ID-420	3/8" diameter	

CUSTOM NUMBERED KEYDOTS AVAILABLE ON REQUEST. KEYDOTS SHOWN ACTUAL SIZE.




- A** Adhesive Remover 108-109
Airflow Caps 068-069
- B** Band Clamps 142-143
Barrier Gowns 130-131
Bowie-Dick Test Pack 150-151
Breakaway Tags 086-087
Brush Assortments 052-053
Brush Organizer 056-057
Brush Organizer with Brush Assortment 056-057
French Catheter Scale 056-057
Pegged Bar For Brush Organizer 056-057
- C** Cable Ties 090-091
Case Locking Tags 084-085
Case Locking Tags With Indicator Dot 084-085
Catheter Trays 164-165
Cleanroom Apparel 134-135
Adhesive Mats 134-135
Beard Covers 134-135
Bouffant Caps 134-135
Coveralls 134-135
Shoe Covers 134-135
Channel Cleaning Brushes 004-017
Acrylic Tip Brushes 026-029
Bone Reamer Brushes 018-019
Fan Tip Brushes 024-025
Lumen Guard™ Cleaning Brushes 021-023
Polypropylene Brushes 030-033
Corner Protectors 144-145
- D** Dome Covers 164-165
Dr. Weigert Neodisher® MultiZym 002-003
Dr. Weigert Neodisher® PreStop 002-003
Duraholder® Instrument Protection System 076-077
- E** Endoscope Brushes
Albarran Brush 048-049


- Brush Assortments 050-051
Brush Kits 042-043
Brushes 044-045
Channel Brushes 040-041
Double Ended 046-047
Premium Brushes 048-049
Pull-Through 046-047
Tapered 046-047
Valve Brushes 048-049
EfferZyme™ 002-03
- F** Face Masks 128-129
Face Shields 126-127
Flat Caps 062-063
Flat Covers 162-163
French Catheter Scale 056-057
- G** Gelpi Retractor Caps 070-071
General Instrument Cleaning Brushes 038-039
Acetabular Reamer Brush 038-039
Double-Ended Brush 038-039
Finger Loop Brush 038-039
Stainless Steel Bristles 038-039
Glove Dispenser 132-133



- H** High-Sided Perforated Trays 160-161
- I** **Identification Tape 106-123**
Indicator Dot Locking Tags 082-083
Instrument Cleaning 000-057
Instrument Protection 058-077
Instrument Sleeves 074-075
Instrument Stringers 138-143
Ball & Socket 138-139
Center Lock 140-141
Pull-To-Open 138-139
Spring Lock 140-141
Stringer Assistant 142-143
Stringer Stand 142-143
T-Bar Stringers 142-143
U-Shaped 140-141
U-Shaped Stringer End Cap 140-141
Instrument Tracking 186-191
Instrument Trays 162-163
- K** K-Wires 168-170
K-Wire Dispensers 174-175
KeyDot 188-191
- L** Lab Markers 090-091
Lint Remover Handle And Refills 146-147
Locking Tags 080-083
- M** Mesh Trays 156-157
Metal Tray Tags 102-103
Micro Mesh Trays 154-155
Minimally Invasive Surgical Instrument (MIS) Caps 072-073
- N** Nitrile Gloves 132-133



- O** Oblong Instrument Trays 160-161
One-Loop Tags 104-105
Osteotome Caps 064-065
- P** **Packaging And Preparation 136-165**
Peel Pouch Sterilization Racks 148-149
Pegged Bar For Brush Organizer 056-057
Perforated Mayo Trays 160-161
Personal Protective Equipment 124-135
Pin Covers 178-179
Pin Gauge 176-177
Pipe Cleaners 034-035
Double Cleaner 034-035
Fleece Cleaner 034-035
Flex Cleaner 034-035
Flex Ultra Fine Cleaner 034-035
Plastic Sterilization Trays 152-153
Plastic Tray Tags 098-101
Pliers 182-183
Polyester Wiper 146-147
Protective Sleeves 130-131



- R** Roll Tape 108-109
- Round Caps 060-061

- S** Scope Baskets 158-159
- Scope Processing Tags 096-097
- Scope Repair Tag 096-097
- Scope Tags 094-095
- Self-Retaining Retractor Caps 070-071
- Sheet Tape 110-123
 - Patterned Sheet Tape 112-123
 - Solid Color Sheet Tape 110-111
- Silicone Finger Mat 153
- Skin Hook Caps 068-069
- Smart Tags 092-093
- Steinmann Pins 171-173
- Sterilizer Cleaning Brushes 054-055
- Sterilizer Gloves 130-131
- Sterilizer Mitts 130-131
- Sterilization Pouches 148-149
- Sterilization Racks 176-177
- Sterilization Wrap 148-149
- Super Sani-Cloth® 146-147

- T** **Tags 078-105**
- Tape Dispenser 150-151
- Tapered Caps 066-067
- Toothbrush-Style Cleaning Brushes 036-037
 - Double Headed Toothbrush-Style 036-037
- Tray Liners 144-145
- Tray Tag Clip 100-101
- Tray Tag Ring 100-101
- Twin Caps 072-073
- Twist Drill Bits 184-185
- Twist Drill Rack 184-185

- W** Wall Mount Sani-Bracket 146-147
- Wire and Pin Cutters 180-181
- Wire And Pin Management 166-185**
- Write-On Tags 088-091
- Wypall™ - L40 Wipers 146-147



