



ENDOSCOPY SPONGES

The best offense is a good defense and for scope reprocessing pre-cleaning is a vital step. Endoscope cleaning sponges from Key Surgical help defend against biofilm; preparing the scope for manual cleaning and further high-level disinfection or sterilization.

“Manual cleaning... is the most important step in removing the microbial burden from an endoscope.”

-SGNA Standards for Reprocessing Endoscopes 2016, pg. 16

ABOUT:

- Convenient, pre-saturated sponge with MetriZyme™
- Designed for pre-cleaning
- Contour shape helps remove gross soil
- Individually wrapped
- 100/case (4 boxes of 25/box)
- Single use



PRODUCT NO. DESCRIPTION

ES-88 MetriSponge™



ABOUT:

- Dry sponge
- Absorbent, contour shape helps remove gross soil during pre-cleaning or manual cleaning
- 50/bag
- Single use

PRODUCT NO. DESCRIPTION

ES-50 Endo Sponge, dry

