



Face Shields

INTENDED USE:

Key Surgical® Face Shields are intended to protect personnel against splashing liquids/sprays.

INSTRUCTIONS FOR USE:

1. Remove face shield from packaging.
2. Place foam front against your forehead.
3. Secure elastic band around the back of your head.
4. Adjust shield for comfort.
5. Discard after single use.

 MANUFACTURED FOR KEY SURGICAL, LLC
8101 WALLACE ROAD / EDEN PRAIRIE, MN 55344 / USA
PH: 800.541.7995 KEYSURGICAL.COM

REF CATALOG
NUMBER

LOT LOT
NUMBER



DO NOT RE-USE

EC REP

MDSS GMBH
SCHIFFGRABEN 41, D-30175 / HANNOVER, GERMANY



BROC-3124, Rev B, 04/2018