

**Dry Cleaning Wipe (WIPE8X8)**

**EN** - **INTENDED USE:** Key Surgical® Dry Cleaning Wipes are intended to remove soils and clean surfaces.

**CONTRAINDICATIONS:** No known contraindications and/or adverse effects.

**PRECAUTIONS:** Store away from a source of ignition. Single use only; reuse may lead to patient cross-contamination.

**PREPARATIONS FOR USE:** Follow the item manufacturer's instructions for use. If utilizing a cleaning agent, be sure to follow the cleaning agent manufacturer's IFU for proper preparation of the cleaning solution.

**INSTRUCTIONS FOR USE:**

1. Open the packaging
2. Take out the product
3. Wipe down item or surface to remove residual soils or moisture
4. Dispose after single use

**DISPOSAL:** Dispose after use according to facility protocol. Do not reuse.

 KEY SURGICAL LLC  
8101 WALLACE ROAD / EDEN PRAIRIE, MN 55344 USA  
(US) TEL: 800.541.7995 (UK) TEL: +44(0)1628 810626  
(EU) TEL: +49 4563 905900

REF

CATALOG  
NUMBER

LOT

LOT  
NUMBER

 MANUFACTURER  
FOR US: MANUFACTURED FOR

UDI

UNIQUE DEVICE IDENTIFIER



DO NOT RE-USE

 DATE OF MANUFACTURE

 **STERIS**®