

Silicone Clamp Covers

EN - INTENDED USE: Key Surgical® Silicone Clamp Covers are disposable, single-use accessories to medical devices intended for transient-term (less than 60 minutes) use in surgical procedures requiring grasping and clamping. The Silicone Clamp Covers are not intended to control, diagnose, monitor or correct a defect of the heart or of the central circulatory system through direct contact with these parts of the body and are not intended for use in direct contact with the central circulatory system or central nervous system.

PRECAUTIONS: This device is single use only; reuse may lead to patient cross contamination and/or device failure. Do not re-sterilize. Re-sterilization of the device may result in loss of product functionality.

CAUTION: Prior to use visually inspect packaging. Do not use if damaged, opened or breached as sterility cannot be guaranteed.

CONTRAINDICATIONS: There are no known contraindications associated with the Silicone Clamp Covers.

INSTRUCTIONS FOR USE:

1. Remove silicone clamp covers from pouch.
2. Clamp covers may be trimmed to desired length.
3. Use during surgical procedure.

DISPOSAL: Upon completion of procedure, dispose of both used and unused clamp covers according to facility protocol. Do not reuse.


KEY SURGICAL LLC
 8101 WALLACE ROAD / EDEN PRAIRIE, MN 55344 USA
 (US) TEL: 800.541.7995 (UK) TEL: +44(0)1628 810626
 (EU) TEL: +49 4363 905900


MANUFACTURER
 FOR US: MANUFACTURED FOR

REF

CATALOG NUMBER

LOT

LOT NUMBER



CONSULT INSTRUCTIONS FOR USE

STERILE | **R**
 STERILIZED USING IRRADIATION

MD

MEDICAL DEVICE



DO NOT USE IF PACKAGE IS DAMAGED



CAUTION



DO NOT RE-USE



DO NOT RESTERILIZE



USE-BY DATE

UDI

UNIQUE DEVICE IDENTIFIER

 SINGLE STERILE BARRIER SYSTEM

 SINGLE STERILE BARRIER SYSTEM WITH PROTECTIVE PACKAGING OUTSIDE

 **STERIS**®