

Instruction for use– Distal tip protector

1. Product Intended Use

The device is used for protecting the distal end of the endoscope during storage and transportation.


2. Application Contraindication

- a. The product should be used before losing efficacy
- b. The product is used for only once and should be destroyed afterwards.
- c. The product should be stored in the shady, cool, dry, well-ventilated and clean condition.

3. The Main Raw Material

Serial No.	Material
1	PU

4. Product drawing



Jiangsu Kangjin Medical Instrument Co., Ltd.

Zhenglu Town, Changzhou, Jiangsu Province, PRC

Distributed by: Key Surgical LLC

8101 Wallace Road, Eden Prairie, MN 55344 USA