




INSTRUMENT STRINGERS

Key Surgical® Instrument Stringers are ideal for organizing ring-handled instruments during washing and sterilization. Various widths are available to ensure that hinged instruments remain in an open position for proper washing and sterilization. Several styles are available to fit your personal preferences in the department.

“Ratcheted instruments should be unlatched. Racks, pins, stringers, or other specifically designed devices can be used to hold the instruments in the unlatched position.”

-ANSI/AAMI ST79:2017, Section 8.2

STRINGER STYLES:



ABOUT:

- Durable, corrosion-resistant 304 (18/8) grade stainless steel
- Holds ring-handled, hinged instruments open for optimal cleaning
- Compatible with steam sterilization
- Various lengths, widths, and closure styles available
- 1/pkg




WIDTHS:

- 2.5", 2.75", 3", 4", or 5.5"

YOU MIGHT ALSO LIKE:



Mesh Trays



Surgical Tray Corner Protectors



Sterilization Wrap

