




LUMEN GUARD™ CLEANING BRUSHES

Surgical instrument cleaning just got a little smoother. Key Surgical® Lumen Guard Cleaning Brushes feature a smooth, inseparable, protective tip. Traditional exposed wire at the end of the brush is transformed through the use of Z-Tip™ technology, helping protect against potential scratches or damage to the lumen. Rigid, antimicrobial bristles aid in effective removal of debris.

ABOUT:

- Rigid, antimicrobial nylon bristles
- Smooth, inseparable tip protects against potential scratches to inner lumen
- Twisted stainless steel handle with loop end
- Diameters ranging from .050" (1.27mm) to 1" (25.4mm)
- Lengths from 6" (15.24cm) to 24" (61cm)
- Tapered style available



YOU MIGHT ALSO LIKE:

Brush Organizer

- Stainless steel
- 16" x 35"
- 12 pegs
- 2 additional bars (adhesive backing) with 4 pegs each



FOR MORE INFORMATION CALL 800.541.7995 OR VISIT KEYSURGICAL.COM



PRODUCT NUMBER	DIAMETER INCHES	DIAMETER MM	FR	BRISTLE AREA LENGTH	OVERALL LENGTH	QUANTITY
BZ-10-050 BZ-10-050-50	.050"	1.27mm	4	1"	10"	2 50
BZ-12-062 BZ-12-062-50	.062"	1.58mm	5	1"	12"	2 50
BZ-18-062 BZ-18-062-50	.062"	1.58mm	5	2"	18"	2 50
BZ-18-079 BZ-18-079-50	.079"	2.00mm	6	2"	18"	2 50
BZ-12-087 BZ-12-087-50	.087"	2.20mm	7	1"	12"	2 50
BZ-12-098 BZ-12-098-50	.098"	2.50mm	8	1"	12"	2 50
BZ-18-105 BZ-18-105-50	.105"	2.67mm	8	2"	18"	2 50
BZ-08-125 BZ-08-125-50	.125"	3.18mm	10	.75"	8"	2 50
BZ-12-125 BZ-12-125-50	.125"	3.18mm	10	1"	12"	2 50
BZ-18-125 BZ-18-125-50	.125"	3.18mm	10	2"	18"	2 50
BZ-24-125 BZ-24-125-50	.125"	3.18mm	10	2"	24"	2 50
BZ-18-131 BZ-18-131-50	.131"	3.33mm	10	2"	18"	2 50
BZ-18-144 BZ-18-144-50	.144"	3.37mm	11	2"	10"	2 50
BZ-12-158 BZ-12-158-50	.158"	4.00mm	12	1"	12"	2 50
BZ-18-158 BZ-18-158-50	.158"	4.00mm	12	2"	18"	2 50
BZ-18-197 BZ-18-197-50	.197"	5.00mm	15	2"	18"	2 50
BZ-18-236 BZ-18-236-50	.236"	6.00mm	18	2"	18"	2 50
BZ-16-250 BZ-16-250-50	.250"	6.35mm	19	2"	16"	2 50
BZ-18-250 BZ-18-250-50	.250"	6.35mm	19	2"	18"	2 50
BZ-18-312 BZ-18-312-50	.312"	8.00mm	24	2"	18"	2 50
BZ-18-354 BZ-18-354-50	.354"	9.00mm	28	2"	18"	2 50
BZ-18-393 BZ-18-393-50	.393"	10.00mm	30	2"	18"	2 50
BZ-18-433 BZ-18-433-50	.433"	11.00mm	34	2"	18"	2 50
BZ-18-472 BZ-18-472-50	.472"	12.00mm	N/A	2"	18"	2 50
BZ-18-1000 BZ-18-1000-50	1"	25.40mm	N/A	2"	18"	2 50
BZ-18-1543 BZ-18-1543-50	.157" to .433"	4mm to 11mm	N/A	2"	18"	2 50
BZ-06-2550 BZ-06-2550-50	.250" to .500"	6.35mm to 12.7mm	N/A	.75"	6"	2 50