



# PEEL POUCH STERILIZATION RACKS

Key Surgical® Peel Pouch Sterilization Racks are designed to keep peel pouches on edge while giving enough space for effective sterilant contact. These stainless steel racks feature multiple slots for peel pouched items, aiding in efficient assembly and drying after sterilization.

“Paper-plastic pouches should stand on edge in relation to the cart or shelf, with the paper side of one pouch next to the plastic side of the next pouch; holding racks or baskets specifically designed for pouches can be used.”


-ANSI/AAMI ST79: 8.5.2

## ABOUT:


- Durable, 304-grade stainless steel
- Organizes peel pouches for the sterilization process
- Sterilize pouched items in accordance with ANSI/AAMI standards




PR-090504-5



PR-210806-8



PR-210806-4



PR-160604-21

## PRODUCT NO.

## DESCRIPTION

PR-090504-5	8.7" x 4.7" x 4", 5 slots, 1.5" apart
PR-210806-4	20.5" x 7.5" x 6", 4 slots, 5" apart
PR-150504	15.5" x 5" x 4", 10 slots, 1.5" apart
PR-210806-8	20.5" x 8" x 6", 8 slots, 2.5" apart
PR-160604-21	16.3" x 5.5" x 4", 21 slots, 0.7" apart



PR-150504



**KEYSURGICAL®**