



AAMI LEVEL 1 & 2 GOWNS

Safety and infection prevention begin with the first interaction between the clinician and patient. Key Surgical now provides AAMI Level isolation gowns designed to protect you based on your need. AAMI standards classify a gown's ability to act as a barrier to penetration by liquids or liquid-borne pathogens. Regardless of the name (procedure, isolation, or cover gown) choose your gowns based on intended use and your desired level of protection based on your risk levels.

PRODUCT NO. — DESCRIPTION

PS-200 AAMI Level 1 Gown

- For use during basic patient interaction (MINIMAL risk situations such as basic care, hospital medical unit, minimal occurrence of splash)
- Features elastic at wrists
- Full closure with ties
- Material: 70% PP + 30% PE
- Material weight: 35 GMS (per sq. m)
- Blue gown
- One size fits all
- Single use
- 50/box



PRODUCT NO. — DESCRIPTION

PS-210 AAMI Level 2 Gown

- For use during non-invasive procedures (LOW risk situations such as blood draw, suturing, low occurrence of splash)
- Features soft cloth cuffs at wrists
- Full closure with ties
- Material: 100% non-woven SMS fabric
- Material weight: 45 GMS (per sq. m)
- Blue gown
- One size fits all
- Single use
- 50/box



a STERIS company