


CHANNEL CLEANING BRUSHES

BONE REAMER BRUSHES




When you've got a challenging instrument to clean, like a flexible bone reamer, you need a tool that's right for the job. Key Surgical® Bone Reamer Brushes are designed with a rigid stainless steel shaft to help you efficiently push out debris during the cleaning process. At 31" long, these help clean the entire length of the instrument for effective manual cleaning. Available in 3 diameters.

- Rigid stainless steel shaft
- 31" length
- Nylon bristles
- Loop end for easy hanging & storage
- 3/pkg



(brush tip details)



Product No.	Description	Qty
BR31118	Bone Reamer Brush, 31" long, rigid shaft, .118" (3mm) diameter	3
BR31197	Bone Reamer Brush, 31" long, rigid shaft, .197" (5mm) diameter	3
BR31276	Bone Reamer Brush, 31" long, rigid shaft, .276" (7mm) diameter	3

VISIT [STERIS.COM/HEALTHCARE](https://www.steris.com/healthcare) OR CONTACT YOUR STERIS REPRESENTATIVE TODAY