




# BONE REAMER BRUSHES

When you've got a challenging instrument to clean, like a flexible bone reamer, you need a tool that's right for the job. Key Surgical® Bone Reamer Brushes are designed with a rigid stainless steel shaft to help you efficiently push out debris during the cleaning process. At 31" in length, these help clean the entire length of the instrument for effective manual cleaning. Available in 3 diameters.




## FEATURES:

- Rigid stainless steel shaft
- 31" length
- Medical-grade nylon bristles
- Reusable design
- Loop end for easy handling & storage
- 3/pkg

PRODUCT NO.	DIAMETER	QUANTITY
BR-31-118	.118" (3mm)	3/pkg
BR-31-197	.197" (5mm)	3/pkg
BR-31-276	.276" (7mm)	3/pkg

### BRUSH TIP DETAIL



### YOU MIGHT ALSO LIKE:



TOOTHBRUSH-STYLE CLEANING BRUSHES



a STERIS company

VISIT OUR WEBSITE [KEYSURGICAL.COM](http://KEYSURGICAL.COM) OR CALL 800.541.7995 FOR MORE INFORMATION

