

CHANNEL CLEANING BRUSHES

Fan Tip and Acrylic Tip




Different cleaning tasks require different cleaning tools. It's important to pick the right brush for the job at hand. Key Surgical® Fan Tip and Acrylic Tip Brushes provide an extra level of protection and cleaning with either a fanned out tip for close-ended instruments or an acrylic tip which is great for lumens or channels with hard to maneuver curves or bends.

Fan Tip Brushes

- Fan-shaped brush tips help clean the distal end of close-ended lumens
- Twisted stainless steel handles with loop end
- Antimicrobial nylon bristles
- Diameters ranging from .197" (5mm) to .500" (12.7mm)
- Bristle area lengths: 2" or 3"
- Overall lengths: 12", 18", or 24"
- Available 2/pkg or 50/pkg

Acrylic Tip Brushes

- Acrylic covered tips provide extra protection during lumen cleaning
- Flexible handles resist kinking or coiling
- Antimicrobial nylon bristles
- Diameters ranging from .062" (1.58mm) to .440" (11.18mm)
- Bristle area lengths: 1", 1.5", 2", or 3"
- Overall lengths: 12", 18", or 24"
- Available 2/pkg or 50/pkg



CONTACT YOUR STERIS REPRESENTATIVE TO LEARN MORE

[steris.com](https://www.steris.com)
Mentor, OH 44060 | USA

Key Surgical is a registered trademark of STERIS
MKTG4008 Rev G, 2024
M12622EN.2024-04, Rev A
©2024 STERIS. All rights reserved.

 STERIS