



CLAMP COVERS



Clamp Covers are designed to reduce vessel trauma during surgical procedures requiring grasping or clamping (particularly when a hemostatic clamp is utilized). Thread-like indents provide non-slip tension for an atraumatic grip.

- Available in 3 colors (red, blue, or white)
- Available in 3 widths (.1" , .2", or .3") to fit most clamp sizes
- 5" lengths
- Sterile and non-sterile options
- Medical-grade silicone, free from natural rubber latex
- Radiopaque



Product No.	Description	Qty
VS2101	Clamp Jaw Covers, .1" x 5", red, sterile	2 per sterile pouch (10 pouches/box)
VS2102	Clamp Jaw Covers, .2" x 5", white, sterile	2 per sterile pouch (10 pouches/box)
VS2103	Clamp Jaw Covers, .3" x 5", blue, sterile	2 per sterile pouch (10 pouches/box)
VS2201	Clamp Jaw Covers, .1" x 5", red, non-sterile	20/pk
VS2202	Clamp Jaw Covers, .2" x 5", white, non-sterile	20/pk
VS2203	Clamp Jaw Covers, .3" x 5", blue, non-sterile	20/pk
VS201100	Clamp Jaw Covers, .1" x 5", red, non-sterile	100/pk
VS202100	Clamp Jaw Covers, .2" x 5", white, non-sterile	100/pk
VS203100	Clamp Jaw Covers, .3" x 5", blue, non-sterile	100/pk

[CONTACT YOUR STERIS REPRESENTATIVE TO LEARN MORE](#)