


DOUBLE-HEADER CLEANING BRUSH

Flexible Endoscope Channel Cleaning Brush

Let this combination brush pull a double-header in endoscope reprocessing, not you. With a classic endoscope channel cleaning brush on one end and a valve/control head brush on the other, you have just what you need for removing debris within endoscope channel diameters ranging from 2.0mm - 7.0mm.

FOR CHANNEL
SIZES
2.0MM - 7.0MM



- 230cm long channel brush
- Double-ended with channel brush on one end, valve/port brush on the other
- Tapered, nylon bristle brush head (5-7mm) with a protective tip at the end to help prevent against potential channel damage
- Valve/port brush attached
 - Full head of soft nylon bristles
 - Barbell handle with grooves for a steady, comfortable grip

Product no.	Description	Bristle Diameter (mm)	Valve Control Brush (mm)	Sheath Diameter (mm)	Overall Length (cm)	Channel Diameter (mm)	Quantity
BX00711609	Double-Header Combination Cleaning Brush	5 - 7 (tapered)	10	1.7	230	2.0 - 7.0	50/box
BX00711629	Double-Header Combination Cleaning Brush	5 - 7 (tapered)	10	1.7	230	2.0 - 7.0	100/box

To learn more, visit [STERIS.COM/HEALTHCARE](https://www.steris.com/healthcare) or contact your STERIS representative today.