



DURAHOLDER™ IPS


The Duraholder Instrument Protection System (IPS) is designed to prevent costly damage to delicate instruments. Comprised of durable Kimguard™ material and available in several designs, you can choose your Duraholder according to your specific requirements. Duraholder IPS is not intended to provide a sterile barrier and should be used in conjunction with an outer wrap or container system.

PHYSICAL FEATURES:

- 5 different designs/pouch configurations
- Built-in pockets and tape closure strips
- Material compatible with steam sterilization

BENEFITS:

- Holds delicate instruments in place within a closed container
- Helps protect against damage to instrumentation
- Improves the organization of instruments
- Alternative solution to peel packs



OTHER RELATED PRODUCTS



a STERIS company

FOR MORE INFORMATION CALL 800.541.7995 OR VISIT KEYSURGICAL.COM

