


ENDOSCOPE CLEANING AND PROTECTION

Cleaning Brushes and Protection for Flexible Endoscopes




Flexible endoscope manual cleaning has many complex steps; choosing the right brush for the job doesn't have to be one of them. Choose from single-ended, double-ended, combination designs, tapered, non-tapered, and our newest technology... a brush with a squeegee that removes more bioburden in just one pass.*

- **Comprehensive range** of diameters, lengths and styles such as single or double-ended brushes, tapered, non-tapered, and **unique styles** such as **DuoSwift®**
- **Protective tips** at the end of the brush help prevent inadvertent damage
- **Valve/port brushes** for cleaning valve/port channels; including a specialized brush for the elevator channel of a duodenoscope
- **Single-use** designs help promote patient safety and support infection prevention guidelines



	PRODUCT NO.	DESCRIPTION	BRUSH DIAMETER (MM)	FITS SCOPE CHANNEL (MM)	LENGTH (CM)	QTY
a.	EB10106	Channel cleaning brush (single-ended)	2	1.2-1.5	120	50
b.	BX00711616	Small Wonder channel cleaning brush	2.25	1.0-2.25	160	15
	BX00711617	Small Wonder channel cleaning brush	2.25	1.0-2.25	230	15
c.	EB10104	Channel cleaning brush with valve brush	5	2.8-4.2	230	25
d.	BX00711609	Double-Header combination cleaning brush	5-7	2.0-7.0	230	50
e.	BX00711607	Gemini channel cleaning brush (double-ended)	5-7 (tapered)	2.0-7.0	230	50
f.	BX00711612	Wooly Brush channel cleaning brush (single-ended)	5-7 (tapered)	2.0-7.0	230	50
	BX00711601	Channel cleaning brush (single-ended)	5-7 (tapered)	2.0-7.0	230	50
g.	BX00711618	DuoSwift Combination Squeegee Brush	5-7 (brush), 2.5-5.2 (squeegee)	2.8-7.0	230	50
	BX00711619	DuoSwift Combination Squeegee Brush (w/valve brush)	5-7 (brush), 2.5-5.2 (squeegee)	2.8-7.0	230	50
h.	BX00711603	Jumbo Brush channel cleaning brush	7-9 (tapered)	2.8-9.0	230	50
	BX00711606	Jumbo Brush channel cleaning brush	7-9 (tapered)	2.8-9.0	350	20

* Test report on file with STERIS



Intended to clean Savary-Gilliard and American dilators, our dilator brush features a long, flexible stainless-steel shaft and a brush head with a protective tip to allow you to maneuver through the lumen as necessary to remove gross soil and debris.


DILATOR BRUSH

PRODUCT NO.	DESCRIPTION	BRUSH DIAMETER (MM)	BRISTLE AREA LENGTH (MM)	OVERALL LENGTH (CM)	QTY
a. EB10140	Dilator cleaning brush	3	12	90	1



VALVE BRUSH PRODUCT INFORMATION

PRODUCT NO.	DESCRIPTION	BRUSH DIAMETER (MM)	QTY
b. EB10101	Double-ended valve brush	5 & 10	100
c. BX00711610	Valve control head cleaning brush	10	50
d. EB10108	Albarra channel brush	2.3	100



ENDOSCOPE TIP PROTECTION

PRODUCT NO.	DESCRIPTION	COMPATIBLE SCOPE OD (MM)	LENGTH (CM)	QTY
e. BX00711623	ProAmor® endoscope tip protector - sterile	4.9-15	10	50
BX00711624	ProAmor endoscope tip protector - non-sterile	4.9-15	10	50
f. BX00711621	Endo-Boot® endoscope tip protector - small	3.0-8.8	14	50
g. BX00711622	Endo-Boot endoscope tip protector - large	8.8-15	14	50

TO LEARN MORE, PLEASE CONTACT YOUR STERIS SALES REPRESENTATIVE TODAY