




ENDOSCOPE BRUSHES

Endoscopes are complex and often difficult to reprocess, which makes choosing the right brush critical for effective manual cleaning. Whether you need channel, valve, or combination brushes, our updated line of brushes has what you need. Our single-use brushes prevent cross-contamination. Let's get to it.

ABOUT:

Various styles available, including kits. Channel brushes feature protective tip, single use to prevent cross-contamination from re-use.



a STERIS company

VISIT OUR WEBSITE [KEYSURGICAL.COM](https://www.keysurgical.com) OR CALL 800.541.7995 FOR MORE INFORMATION



MKTG-4024 Rev E, 2021



ENDOSCOPE BRUSHES

PRODUCT NO.	DESCRIPTION	BRUSH DIAMETER	BRISTLE AREA LENGTH	LENGTH	QTY
EB10102	Double-Ended	3mm	20mm	2300mm	50/pkg
EB10103	Double-Ended	6mm	20mm	2500mm	50/pkg
EB10103.3	Triple-Cleaning (3 bristle areas)	6mm	20mm	2500mm	50/pkg
EB10136	Double-Ended	3mm	20mm	1200mm	50/pkg
EB10104	Double-Ended with Valve Brush Attached	5mm & 10mm	20mm & 40mm	2300mm	25/pkg
EB10105	Double-Ended, Brush and Pipe Cleaner	5mm	20mm & 125mm	2300mm	50/pkg
EB10106	Single-Ended	2mm	10mm	1200mm	50/pkg
EB10103.3C	Brush Kit (Includes Triple-Cleaning EB10103.3 and Valve Brush EB10101)				50/pkg
EB10107C	Brush Kit (Includes Double-Ended Channel EB10103 and Valve Brush EB10101)				50/pkg
EB10107G	Brush Kit (Includes Double-Ended Channel EB10102 and Valve Brush EB10101)				50/pkg
EB10101	Valve Brush	5mm & 10mm	20mm & 30mm		100/pkg
EB10108	Albarran Brush	2.3mm	13mm		100/pkg



a STERIS company

VISIT OUR WEBSITE KEYSURGICAL.COM OR CALL 800.541.7995 FOR MORE INFORMATION



MKTG-4024, Rev D 2020