

INSTRUMENT ORGANIZER



Key Surgical® Instrument Organizers are designed to help with efficient organization and protection of surgical instruments during the sterilization process. Designed to be used in conjunction with a sterile barrier system, the organizers feature multiple slots and rows to hold and protect your instruments. A place for everything, and everything in its place.


- Designed to help you efficiently organize and protect surgical instrumentation during the sterilization process
- Multiple rows of slots and various pocket configurations to fit a wide variety of instruments
- Each organizer can be cut between double-seams to create individual pieces so you can customize to your specific facility needs (does not apply to CSDH505 and CSDH512)
- SMS material provides strength and durability for protection while allowing steam sterilization of contents to occur*
- Instrument Organizers are not intended to provide a sterile barrier and must be used with a sterile barrier system to maintain sterility (see IFU)
- 100/case
- Single use



CSDH505




CSDH512




CSDH200


CSDH200 features self-adhesive strip



CSDH100



CSDH102



CSDH300

CONTACT YOUR STERIS REPRESENTATIVE TO LEARN MORE

[steris.com](https://www.steris.com)
Mentor, OH 44060 | USA

Key Surgical is a registered trademark of STERIS
MKTG4103 Rev B, 2024
M12706EN.2024-04, Rev A
©2024 STERIS. All rights reserved.



*Allows sterilization of contents to occur in a steam sterilizer using a sterilization and drying cycle that has been validated by the user for their equipment and procedures