

INSTRUMENT SLEEVES



Instrument protection is a vital part of the life cycle of an instrument. Using Key Surgical® instrument protection, such as Instrument Sleeves, supports a proactive plan for protecting instruments from damage during sterilization, storage, and transportation. Additionally, they can help prevent sharps injuries during handling. Use with hinged or non-hinged instruments to protect instrumentation during sterilization to point of use in the operating room.

“Ratcheted instruments should be unlatched. Racks, pins, stringers, or other specifically designed devices can be used to hold the instruments in the unlatched position.”

- AAMI/ANSI ST79, Section 8.2

- Helps **protect distal ends** of instruments (hinged or non-hinged) from damage during prep & pack processes, sterilization, and transportation
- **Transparent pocket** in front for visualization of distal end
- Instrument Sleeve **01457** is designed specifically **for hinged** instruments, featuring a folder over tab to keep instruments open
- Compatible with steam sterilization (refer to IFU)
- Single use



Product No.	Description	Dimensions	Qty
01455	Instrument Sleeve, paper/plastic, fits instruments up to 32mm wide. Designed for non-hinged instruments	2" x 5"	100
01457	Instrument Sleeve, paper/plastic, fits instruments up to 40mm wide. Designed for hinged instruments	3.5" x 6.6"	100

VISIT [STERIS.COM/HEALTHCARE](https://www.steris.com/healthcare) OR CONTACT YOUR STERIS REPRESENTATIVE TODAY