



INSTRUMENT ORGANIZER


Key Surgical Instrument Organizers are designed to help with efficient organization and protection of surgical instruments during the sterilization process. Designed to be used in conjunction with a sterile barrier system, the organizers feature multiple slots and rows to hold and protect your instruments. A place for everything, and everything in its place.

ABOUT:


- Designed to help you efficiently organize and protect surgical instrumentation during the sterilization process
- Multiple rows of slots and various pocket configurations to fit a wide variety of instruments
- Each organizer can be cut between double-seams to create individual pieces so you can customize to your specific facility needs (does not apply to CS-DH505 and CS-DH512)
- SMS material provides strength and durability for protection while allowing steam sterilization of contents to occur*
- Instrument Organizers are not intended to provide a sterile barrier and must be used with a sterile barrier system to maintain sterility. (see IFU).
- 100/case
- Single use



CS-DH505




CS-DH512




CS-DH200


CS-DH200
features self-
adhesive strip



CS-DH100



CS-DH102



CS-DH300



a STERIS company

FOR MORE INFORMATION CALL 800.541.7995 OR VISIT [KEYSURGICAL.COM](https://www.keysurgical.com)



*Allows sterilization of contents to occur in a steam sterilizer using a sterilization and drying cycle that has been validated by the user for their equipment and procedures

MKTG-4103 Rev A, 2022