



INSUFFLATION TUBING

Barb. Quick. Luer. 22mm. ISO. Which one is your favorite? No matter which connector style you choose, our Insufflation Tubing is designed to work with all standard systems and is sterile and ready-for-use. Prevent cross-contamination during your next (and future) laparoscopic procedures.

ABOUT:

- Used for rapid insufflation of abdominal cavity for laparoscopic procedures
- Lightweight, 10ft. tubing with filter
- High flow insufflation tubing
- Prevents particulate and bacterial infiltration from ambient air being delivered through the tube to the patient
- Virus Filtration Efficiency (VFE) of 99.999% for viruses
- Meets all ERCI and OTT recommendations
- Not made with natural rubber latex
- Sterile
- Single use
- 20/box



FL-2010



FL-2107



FL-2110



FL-2407



FL-2507

PRODUCT NO.

DESCRIPTION

FL-2010	Insufflation Tubing, Barb Connector
FL-2107	Insufflation Tubing, Quick Connector
FL-2110	Insufflation Tubing, Luer Connector
FL-2407	Insufflation Tubing, 22mm Connector
FL-2507	Insufflation Tubing, ISO Connector

