


MESH TRAYS



Key Surgical® Mesh Trays are designed for use during washing, sterilization, storage, and transportation of surgical instruments and/or sets.

- Constructed of durable 304-grade stainless steel mesh
- Standard Mesh Trays feature a mesh size of 5mm x 5mm
- Micro Mesh Trays feature a mesh size of 2mm x 2mm
- Choose from a wide range of sizes to support various needs in sterile processing. Small trays can be used holding small items such as screws, clips, etc. while larger trays can be used for instrument sets.
- Styles: Micro Mesh Tray with attached lid or standard Mesh Tray with drop handles (no lid)



Product No.	Description	Quantity
MT5500	Micro Mesh Tray with lid, 1.57" x 1.57" x .79" (40mm x 40mm x 20mm)	1
MT5600	Micro Mesh Tray with lid, 3.15" x 1.57" x .79" (80mm x 40mm x 20mm)	1
MT5650	Micro Mesh Tray with lid, 6.30" x 1.57" x .79" (160mm x 40mm x 20mm)	1
MT5700	Micro Mesh Tray with lid, 2.75" x 2.75" x 1.18" (70mm x 70mm x 30mm)	1
MT5800	Micro Mesh Tray with lid, 3.15" x 3.15" x 1.57" (80mm x 80mm x 40mm)	1
MT5900	Micro Mesh Tray with lid, 3.54" x 3.54" x 1.97" (90mm x 90mm x 50mm)	1
MT6020	Micro Mesh Tray with lid, 7.09" x 5.12" x .79" (180mm x 130mm x 20mm)	1
MT7100	Micro Mesh Tray with lid, 8.89" x 4.33" x 1.97" (225mm x 110mm x 50mm)	1
MT7200	Micro Mesh Tray with lid, 10.63" x 5.31" x 1.57" (270mm x 135mm x 40mm)	1
MT7300	Micro Mesh Tray with lid, 15.75" x 9.84" x 1.97" (400mm x 250mm x 50mm)	1
MT9001	Micro Mesh Tray with drop handles, 18.90" x 9.84" x 1.97" (480mm x 250mm x 50mm)	1
MT9101	Micro Mesh Tray with drop handles, 9.84" x 9.45" x 1.97" (250mm x 240mm x 50mm)	1
MT9521	Micro Mesh Tray with drop handles, 27.56" x 8.27" x 1.97" (700mm x 210mm x 50mm)	1
MT8005	Mesh Tray with drop handles, 10.24" x 7.46" x 1.97" (260mm x 190mm x 50mm)	1
MT8050	Mesh Tray with drop handles, 14.96" x 9.84" x 1.97" (380mm x 250mm x 50mm)	1
MT8500	Mesh Tray with drop handles, 19.69" x 5.91" x 3.94" (500mm x 150mm x 100mm)	1
MT9000	Mesh Tray with drop handles, 18.90" x 9.84" x 1.97" (480mm x 250mm x 50mm)	1
MT9010	Mesh Tray with drop handles, 18.90" x 9.84" x 2.75" (480mm x 250mm x 70mm)	1
MT9100	Mesh Tray with drop handles, 9.84" x 9.84" x 1.97" (250mm x 250mm x 50mm)	1
MT9110	Mesh Tray with drop handles, 9.84" x 9.45" x 2.75" (250mm x 240mm x 70mm)	1
MT9210	Mesh Tray with drop handles, 17.32" x 4.33" x 2.56" (440mm x 110mm x 65mm)	1

[CONTACT YOUR STERIS REPRESENTATIVE TO LEARN MORE](#)