




MESH TRAYS

Looking for a container to hold small items such as clips, screws, and box locks? How about a durable tray to help with washing and rinsing in decontamination? Or just a really great tray for building instrument sets? Key Surgical® Mesh Trays will do the trick. Mesh Trays are durable and designed to allow adequate sterilant contact of items contained within.


“Organizing containers and other organizing accessories may be placed in the set if they are designed and intended for sterilization.”
-ANSI/AAMI ST79:2017, Section 8.2

ABOUT:


- Durable 304-grade stainless steel
- Great for protecting and holding small items
- Perfect for building various instrument sets
- Wide variety of sizes
- Mesh Trays: 5mm x 5mm mesh
- Micro Mesh Trays: 2mm x 2mm mesh
- Available with attached lids or drop handles



YOU MIGHT ALSO LIKE



Silicone Finger Mat



Surgical Tray Corner Protectors



Sterilization Wrap



a STERIS company

VISIT OUR WEBSITE [KEYSURGICAL.COM](https://www.keysurgical.com) OR CALL 800.541.7995 FOR MORE INFORMATION





MESH TRAYS DIMENSIONS

MICRO MESH TRAYS: 2mm x 2mm mesh

PRODUCT NUMBER	DESCRIPTION	LENGTH	WIDTH	HEIGHT
MT-5500	Micro Mesh with Lid	1.57" (40mm)	1.57" (40mm)	.79" (20mm)
MT-5600	Micro Mesh with Lid	3.15" (80mm)	1.57" (40mm)	.79" (20mm)
MT-5650	Micro Mesh with Lid	6.30" (160mm)	1.57" (40mm)	.79" (20mm)
MT-5700	Micro Mesh with Lid	2.75" (70mm)	2.75" (70mm)	1.18" (30mm)
MT-5800	Micro Mesh with Lid	3.15" (80mm)	3.15" (80mm)	1.57" (40mm)
MT-5900	Micro Mesh with Lid	3.54" (90mm)	3.54" (90mm)	1.97" (50mm)
MT-6020	Micro Mesh with Lid	7.09" (180mm)	5.12" (130mm)	.79" (20mm)
MT-7100	Micro Mesh with Lid	8.89" (225mm)	4.33" (110mm)	1.97" (50mm)
MT-7200	Micro Mesh with Lid	10.63" (270mm)	5.31" (135mm)	1.57" (40mm)
MT-7300	Micro Mesh with Lid	15.75" (400mm)	9.84" (250mm)	1.97" (50mm)
MT-9001	Micro Mesh with Drop Handles	18.90" (480mm)	9.84" (250mm)	1.97" (50mm)
MT-9101	Micro Mesh with Drop Handles	9.84" (250mm)	9.45" (240mm)	1.97" (50mm)
MT-9521	Micro Mesh with Drop Handles	27.56" (700mm)	8.27" (210mm)	1.97" (50mm)

MESH TRAYS: 5mm x 5mm mesh

PRODUCT NUMBER	DESCRIPTION	LENGTH	WIDTH	HEIGHT
MT-8005	Mesh Tray with Drop Handles	10.24" (260mm)	7.48" (190mm)	1.97" (50mm)
MT-8050	Mesh Tray with Drop Handles	14.96" (380mm)	9.84" (250mm)	1.97" (50mm)
MT-8500	Mesh Tray with Drop Handles	19.69" (500mm)	5.91" (150mm)	3.94" (100mm)
MT-9000	Mesh Tray with Drop Handles	18.90" (480mm)	9.84" (250mm)	1.97" (50mm)
MT-9010	Mesh Tray with Drop Handles	18.90" (480mm)	9.84" (250mm)	2.75" (70mm)
MT-9100	Mesh Tray with Drop Handles	9.84" (250mm)	9.84" (250mm)	1.97" (50mm)
MT-9110	Mesh Tray with Drop Handles	9.84" (250mm)	9.45" (240mm)	2.75" (70mm)
MT-9210	Mesh Tray with Drop Handles	17.32" (440mm)	4.33" (110mm)	2.56" (65mm)



a STERIS company

VISIT OUR WEBSITE [KEYSURGICAL.COM](https://www.keysurgical.com) OR CALL 800.541.7995 FOR MORE INFORMATION



MKTG-4004 Rev H, 2022