



# PLASTIC STERILIZATION TRAYS


Key Surgical® Plastic Sterilization Trays are designed to store and protect delicate instruments during sterilization, storage, and transportation. Constructed of a polymer resin to provide strength, durability and chemical resistance, these trays are still lightweight and incredibly easy to use.

## DESCRIPTION:


- Highly durable, heat-resistant plastic base and lid
- Silicone finger mat protects and helps organize instruments
- Slots in lid and base facilitate the circulation of steam and aid in the drying process
- Raised nodes keep silicone mat from sticking to the base
- Transparent lid for easy viewing of contents with snap-shut closure
- Sold as a full unit – base, lid, mat
- Compatible with steam and gas plasma sterilization

## DIMENSIONS:


- PST-1100SB: 7.5" x 2.5" x .75"
- PST-2100SB: 6.5" x 4" x .75"
- PST-3100SB: 10" x 6" x .75"
- PST-3200SB: 10" x 6" x 1.5"
- PST-3300SB: 10" x 6" x 1.5"  
(double layer with insert tray and mat)




PST-3200SB




PST-1100SB



PST-2100SB



PST-3100SB




PST-3300SB

## YOU MIGHT ALSO LIKE:

### Silicone Finger Mat

Individual mat that can be cut to size for specific size needs. Versatile product may be used for replacements within Key Surgical Plastic Sterilization Trays.

- FM-1810: 18.2" x 10.2"



FM-1810



a STERIS company

FOR MORE INFORMATION CALL 800.541.7995 OR VISIT [KEYSURGICAL.COM](http://KEYSURGICAL.COM)

