

CHANNEL CLEANING BRUSHES

POLYPROPYLENE BRUSHES





Looking for a channel cleaning brush with just a little extra 'strength' to tackle challenging cleaning jobs? Try the Key Surgical® Polypropylene Brushes. These channel brushes feature rigid bristles (stiffer than standard nylon bristles) that help effectively and efficiently remove bioburden from instrument lumens/channels. The same high-quality twisted stainless steel handle and loop end* as our standard channel cleaning brushes but with rigid 'tough' bristles.

(*pull-through styles excluded)

- Channel cleaning brushes featuring rigid bristles for use during manual cleaning/decontamination. **Rigid bristles** help with efficient removal of bioburden from instrument lumens/channels
- Twisted stainless steel handle to help provide control as the brush is guided through the instrument during cleaning
- Brushes feature a loop end (excluding the pull-through styles)



- Brush diameters range from .039" (1mm) to .75" (19.05mm)
- Overall lengths: 16", 18", 22", 24", and 28"
- Bristle area lengths: 1" , 2", and 6"
- Available either 2/pkg or 50/pkg



See reverse side for individual product details

Product No.	Diameter (IN)	Diameter (MM)	FR	Bristle Area Length	Overall Length	Quantity
PP16039 PP1603950	.039"	1.00mm	3	1"	16"	2 50
PP22039 PP2203950	.039"	1.00mm	3	1"	22"	2 50
PP16062 PP1606250	.062"	1.58mm	5	1"	16"	2 50
PP22062 PP2206250	.062"	1.58mm	5	1"	22"	2 50
PP16087 PP1608750	.087"	2.20mm	7	1"	16"	2 50
PP22087 PP2208750	.087"	2.20mm	7	1"	22"	2 50
PP16105 PP1610550	.105"	2.67mm	8	2"	16"	2 50
PP22105 PP2210550	.105"	2.67mm	8	2"	22"	2 50
PP16125 PP1612550	.125"	3.18mm	10	2"	16"	2 50
PP22125 PP2212550	.125"	3.18mm	10	2"	22"	2 50
PP16158 PP1615850	.158"	4.00mm	12	2"	16"	2 50
PP22158 PP2215850	.158"	4.00mm	12	2"	22"	2 50
PP18197 PP1819750	.197"	5.00mm	15	2"	18"	2 50
PP22197 PP2219750	.197"	5.00mm	15	2"	22"	2 50
PP18250 PP1825050	.250"	6.35mm	19	2"	18"	2 50
PP24250 PP2425050	.250"	6.35mm	19	2"	24"	2 50
PP18393 PP1839350	.393"	10.00mm	30	2"	18"	2 50
PP18591 PP1859150	.591"	15.00mm	N/A	2"	18"	2 50
PP18750 PP1875050	.750"	19.05mm	N/A	6"	18"	2 50

PULL-THROUGH STYLE

PT28118 PT2811850	.118"	3.00mm	9	2"	28"	2 50
PT28197 PT2819750	.197"	5.00mm	15	2"	28"	2 50

TO LEARN MORE VISIT [STERIS.COM/HEALTHCARE](https://www.steris.com/healthcare) OR CONTACT YOUR STERIS REPRESENTATIVE TODAY