




POLYPROPYLENE BRUSHES

Looking for a channel cleaning brush that will give you the extra strength needed for challenging jobs? Key Surgical® Polypropylene Brushes are designed to help effectively remove soil and bioburden within a lumen or channel on a surgical instrument. The brushes feature polypropylene bristles that are much stiffer than standard nylon yet remain safe to use on the instrument. The durable, twisted stainless steel handle helps provide control as the brush is guided through the instrument during cleaning.

FEATURES:

- Non-absorbent bristles; maintains rigidity when wet
- Naturally resistant to bacterial growth
- Twisted stainless steel handle
- Looped end or pull-through style (centered bristles)
- Available lengths: 16", 18", 22", 24", and 28"
- Bristle area lengths: 1", 2", and 6"
- Available 2/pkg or 50/pkg



RELATED PRODUCTS:



TOOTHBRUSH-STYLE
CLEANING BRUSHES



BRUSH
ORGANIZER



a STERIS company

FOR MORE INFORMATION CALL 800.541.7995 OR VISIT [KEYSURGICAL](https://www.keysurgical.com).



MKTG-4025 Rev F, 2021



POLYPROPYLENE BRUSHES SPECIFICATIONS

LOOPED-END

PRODUCT NUMBER	DIAMETER INCHES	DIAMETER MM	FR	BRISTLE AREA LENGTH	OVERALL LENGTH	QUANTITY
PP-16-039	.039"	1.00mm	3	1"	16"	2
PP-16-039-50						50
PP-22-039	.039"	1.00mm	3	1"	22"	2
PP-22-039-50						50
PP-16-062	.062"	1.58mm	5	1"	16"	2
PP-16-062-50						50
PP-22-062	.062"	1.58mm	5	1"	22"	2
PP-22-062-50						50
PP-16-087	.087"	2.20mm	7	1"	16"	2
PP-16-087-50						50
PP-22-087	.087"	2.20mm	7	1"	22"	2
PP-22-087-50						50
PP-16-105	.105"	2.67mm	8	2"	16"	2
PP-16-105-50						50
PP-22-105	.105"	2.67mm	8	2"	22"	2
PP-22-105-50						50
PP-16-125	.125"	3.18mm	10	2"	16"	2
PP-16-125-50						50
PP-22-125	.125"	3.18mm	10	2"	22"	2
PP-22-125-50						50
PP-16-158	.158"	4.00mm	12	2"	16"	2
PP-16-158-50						50
PP-22-158	.158"	4.00mm	12	2"	22"	2
PP-22-158-50						50
PP-18-197	.197"	5.00mm	15	2"	18"	2
PP-18-197-50						50
PP-22-197	.197"	5.00mm	15	2"	22"	2
PP-22-197-50						50
PP-18-250	.250"	6.35mm	19	2"	18"	2
PP-18-250-50						50
PP-24-250	.250"	6.35mm	19	2"	24"	2
PP-24-250-50						50
PP-18-393	.393"	10.00mm	30	2"	18"	2
PP-18-393-50						50
PP-18-591	.591"	15.00mm	N/A	2"	18"	2
PP-18-591-50						50
PP-18-750	.750"	19.05mm	N/A	6"	18"	2
PP-18-750-50						50

PULL-THROUGH

PT-28-118	.118"	3.00mm	9	2"	28"	2
PT-28-118-50						50
PT-28-197	.197"	5.00mm	15	2"	28"	2
PT-28-197-50						50