





SCOPE BASKETS

Key Surgical® Scope Baskets are designed to hold rigid scopes in place during the washing and sterilization processes. We know that your scopes are an investment which is why our baskets are perfect for holding your delicate instrumentation in place and protecting it from damage.

FEATURES:

- Durable, stainless steel mesh design
- Silicon brackets help hold scopes in place during sterilization, storage, and transportation
- Fully removable lid, easy-to-use lock mechanism
- Compatible with steam, EtO, and gas plasma sterilization
- Various sizes available

PRODUCT NO.	DESCRIPTION
SB-1210	7.08" (180mm) x 3.15" (80mm) x 2.05" (52mm)
SB-4210	18.11" (460mm) x 3.15" (80mm) x 2.05" (52mm)
SB-4610	18.11" (460mm) x 6.30" (160mm) x 2.05" (52mm)
SB-6520	25.19" (640mm) x 5.90" (150mm) x 3.03" (77mm) *includes small mesh tray for accessories



RELATED PRODUCTS:

