

# SOAKER SHEETS




Key Surgical® Soaker Sheets are specifically designed for heavy duty use, such as reprocessing orthopedic sets which can be the most challenging for **packaging integrity** and **wet packs** due to weight.

- Designed for **use with heavier trays** (such as orthopedic trays) that are prone to condensation challenges
- Provides **absorption of condensate** during sterilization to help prevent wet packs
- Available in various sizes
- Manufactured with low chloride and sulfate content. Free from optical dyes
- May be used as a shelf liner (specific item numbers only)
- Compatible with steam and EO sterilization (refer to IFU)
- Single use

**"Absorbent tray liners can absorb condensation and reduce the incidents of wet packs."**

**-ANSI/AAMI ST79:2017**



Product No.	Description	Dimensions	Quantity
SOK006A	Soaker Sheet, cellulose, 140gsm	250mm x 250mm (9.84" x 9.84")	800
SOK003A	Soaker Sheet, cellulose, 140gsm	305mm x 610mm (12" x 24")	400
SOK002A	Soaker Sheet, cellulose, 140gsm	500mm x 375mm (19.69" x 14.76")	200
SOK010A*	Soaker Sheet, cellulose, 140gsm	500mm x 914mm (19.69" x 36")	96
SOK008A*	Soaker Sheet, cellulose, 140gsm	560mm x 1219mm (22" x 48")	100
SOK009A*	Soaker Sheet, cellulose, 140gsm	560mm x 1524mm (22" x 60")	100
SOK001A	Soaker Sheet, cellulose, 140gsm	750mm x 500mm (29.53" x 19.69")	96

\*may be used as a shelf liner

**CONTACT YOUR STERIS REPRESENTATIVE TO LEARN MORE**

**steris.com**  
Mentor, OH 44060 | USA

MKTG4114 Rev B, 2024  
Key Surgical® is a registered trademark of STERIS  
M12941EN.2024-07, Rev.B  
©2024 STERIS. All rights reserved.

