

# SOAKER SHEETS, TRAY MATS & TRAY LINERS

Managing Moisture in Sterilization Trays



## SOAKER SHEETS | 140GSM

- Designed for use with heavier trays (such as orthopedic trays) that are more prone to condensation challenges
- Use under trays to protect wraps from sharp corners and uneven surfaces of trays and baskets
- Compatible with steam and EtO sterilization





Image represents: MAT001, MAT002, MAT003



Image represents: DRYMAT001, DRYMAT002, DRYMAT003, DRYMAT018

## TRAY MATS | 60GSM - 75GSM

- Choose from two styles depending on your sterilization requirements
- Available in multiple sizes to fit a variety of tray/basket types
- Compatible with steam, VHP, or gas plasma (style dependent)

## TRAY LINERS | 60GSM - 107GSM

- Use under tray or within tray
- Available in various sizes, materials, and styles
- Compatible with steam and EtO sterilization



	Soaker Sheets	Tray Liners	Tray Mats	Tray Mats (Dry)
	140 gsm (cellulose)	80 gsm (cellulose)	60 gsm (polypropylene)	75 gsm (viscose-polypropylene)
Can help with moisture management	•	•		•
Can help with reducing tears	•	•	•	•
Compatible with steam sterilization	•	•	•	•
Compatible with EtO sterilization	•	•		
Compatible with low temperature sterilization			•	
Best for trays ≥13 lbs	•			
Best for trays <13 lbs		•	•	•
Best for batteries or sharp instruments			•	
Best for general instrument sets, trays with sharp corners				•

## SOAKER SHEETS

Product No.	Description	Sterilization	Material	Dimensions	Quantity
SOK006A	Soaker Sheet, 140gsm	Steam, EtO	Cellulose	250mm x 250mm (9.84" x 9.84")	800
SOK003A	Soaker Sheet, 140gsm	Steam, EtO	Cellulose	305mm x 610mm (12" x 24")	400
SOK002A	Soaker Sheet, 140gsm	Steam, EtO	Cellulose	500mm x 375mm (19.69" x 15")	200
SOK010A	Soaker Sheet, 140gsm	Steam, EtO	Cellulose	500mm x 914mm (19.69" x 36")	96
SOK008A	Soaker Sheet, 140gsm	Steam, EtO	Cellulose	560mm x 1219mm (22" x 48")	100
SOK009A	Soaker Sheet, 140gsm	Steam, EtO	Cellulose	560mm x 1524mm (22" x 60")	100
SOK001A	Soaker Sheet, 140gsm	Steam, EtO	Cellulose	750mm x 500mm (29.50" x 19.69")	96

## TRAY MATS

Product No.	Description	Sterilization	Material	Dimensions	Quantity
DRYMAT018	Tray Mat, 75gsm	Steam	Polypropylene/viscose	305mm x 305mm (12" x 12")	100
DRYMAT003	Tray Mat, 75gsm	Steam	Polypropylene/viscose	305mm x 610mm (12" x 24")	100
DRYMAT002	Tray Mat, 75gsm	Steam	Polypropylene/viscose	500mm x 375mm (19.69" x 15")	100
DRYMAT001	Tray Mat, 75gsm	Steam	Polypropylene/viscose	750mm x 500mm (29.50" x 19.69")	50
MAT003	Tray Mat, 60gsm	Steam, VHP, and gas plasma	Polypropylene	305mm x 610mm (12" x 24")	100
MAT002	Tray Mat, 60gsm	Steam, VHP, and gas plasma	Polypropylene	500mm x 375mm (19.69" x 15")	100
MAT001	Tray Mat, 60gsm	Steam, VHP, and gas plasma	Polypropylene	750mm x 500mm (29.50" x 19.69")	50

## TRAY LINERS

Product No.	Description	Sterilization	Material	Dimensions	Quantity
TL004A	Tray Liner, 80gsm	Steam, EtO	Cellulose	200mm x 305mm (7.87" x 12")	3,024
TL005A	Tray Liner, 80gsm	Steam, EtO	Cellulose	250mm x 250mm (9.84" x 9.84")	2,016
TL006A	Tray Liner, 80gsm	Steam, EtO	Cellulose	250mm x 500mm (9.84" x 19.69")	1,008
TL007A	Tray Liner, 80gsm	Steam, EtO	Cellulose	305mm x 410mm (12" x 16.14")	1,008
TL008A	Tray Liner, 80gsm	Steam, EtO	Cellulose	457mm x 550mm (18" x 21.65")	504
TL009A	Tray Liner, 80gsm	Steam, EtO	Cellulose	305mm x 305mm (12" x 12")	2,016
TL010A	Tray Liner, 80gsm	Steam, EtO	Cellulose	305mm x 457mm (12" x 18")	1,008
TL011A	Tray Liner, 80gsm	Steam, EtO	Cellulose	305mm x 610mm (12" x 24")	1,008
TL012A	Tray Liner, 80gsm	Steam, EtO	Cellulose	457mm x 457mm (18" x 18")	504
TL013A	Tray Liner, 80gsm	Steam, EtO	Cellulose	457mm x 610mm (18" x 24")	504
TLP012018	Tray Liner, 60gsm	Steam, EtO	Poly/cellulose	305mm x 457mm (12" x 18")	600
TLP020025	Tray Liner, 60gsm	Steam, EtO	Poly/cellulose	508mm x 635mm (20" x 25")	400
51448	Tray Liner, 107gsm (folded sides)	Steam, EtO	Cellulose/crepe	350mm x 580mm (13.78" x 22.83")	500
50110	Tray Liner, 107gsm (folded sides)	Steam, EtO	Cellulose/crepe	350mm x 350mm (13.78" x 13.78")	500
49937	Tray Liner, 107gsm (folded sides)	Steam, EtO	Cellulose/crepe	400mm x 600mm (15.75" x 23.62")	500

Contact your STERIS representative to learn more