

STERILIZATION RACKS



FOR K-WIRES, STEINMANN PINS, OR TWIST DRILL BITS



A selection of sterilization racks designed to help you conveniently wash, sterilize, store and dispense K-wires, Steinmann Pins, or Twist Drill Bits. Dimensions clearly marked on the rack to assist with quick identification. Folds closed but easily converts to a tabletop stand for use in the O.R.

K-WIRE RACK | K-WIRE AND PIN RACK


- Sterilization racks constructed of anodized aluminum, durable yet lightweight
- Lids fold back to act as a support, holding the rack in an upright, optimized position for use in the O.R.
- Wire and pin diameters clearly marked on racks provide efficient identification
- Compatible with steam sterilization



Product No.	Description	Qty
a. KI71RACK46	K-wire Rack, 6.25" x 4.135" x .65", use with 4-6" wires. Holds 6 each of: .028" (.7mm), .035" (.9mm), .045" (1.1mm), .054" (1.4mm), and .062" (1.6mm)	1
b. KI71RACK	K-wire and Pin Rack, 12.5" x 5.40" x 1.9", use with 9" wires and pins. Holds 6 each of: .028" (.7mm), .035" (.9mm), .045" (1.1mm), .054" (1.4mm), .062" (1.6mm), .078" (2mm), .093" (2.4mm), .109" (2.8mm), .125" (3.2mm), .140" (3.5mm), .156" (4mm), and .177" (4.5mm)	1

TWIST DRILL BIT RACK

- **Holds 13 drill bits** ranging from 1.5mm to 5mm in diameter. Drill bit diameters are marked inside rack for to help with loading and identification
- Rack constructed of **U.S. surgical stainless steel**, durable and withstands multiple uses
- Opens to a tabletop stand/position for efficient use in the O.R.
- Compatible with steam sterilization



Rack holds the following size drill bits: (sold separately)

Diameter	Length	Max held
1.5mm	5"	2
2.0mm	5"	2
2.5mm	5"	2
2.7mm	5"	1
3.2mm	5"	1
3.5mm	5"	1
4.0mm	5"	1
4.5mm	5"	1
4.7mm	5"	1
5.0mm	5"	1

Product No.	Description	Qty
c. KI05RACKTD	Twist Drill Bit Rack, 6.75" x 4.25" x .75", holds 13 drill bits	1

VISIT [STERIS.COM/HEALTHCARE](https://www.steris.com/healthcare) OR CONTACT YOUR STERIS REPRESENTATIVE TODAY