


STERILIZATION WRAP



Key Surgical® Sterilization Wrap is uniquely strong *and* flexible so you don't have to compromise one benefit for the other. Featuring the strength and durability (double-bonded, SMS) needed to help protect against microorganisms while allowing sterilization, our sterilization wrap also has the flexibility that can make wrapping even the most challenging sizes/shapes a lot more enjoyable.

- Up to 39% stronger* material is paired with flexibility for a top-performing sterile barrier
- Each sheet is **tear-resistant**, double bonded (SMS)
- Easy to use, **flexible** wrap for securely wrapping trays, devices, and other surgical items
- **Unique** diamond pattern and dual-color (blue/ green) supports visualization of holes or tears during **aseptic presentation** in the operating room
- Validated for steam, VHP, EtO, and gas plasma
- Single use

FOR WRAPPING
TECHNIQUES, PLEASE
REFER TO AAMI ST79



Product No.	Weight Wrap	Description	Quantity
WR1002424	100	24" x 24"	240/CASE

Product No.	Weight Wrap	Description	Quantity
WR2001818	200	18" x 18"	480/CASE
WR2002424	200	24" x 24"	240/CASE
WR2003030	200	30" x 30"	144/CASE
WR2003636	200	36" x 36"	144/CASE

Product No.	Weight Wrap	Description	Quantity
WR3002424	300	24" x 24"	240/CASE
WR3003030	300	30" x 30"	144/CASE
WR3003636	300	36" x 36"	144/CASE
WR3005454	300	54" x 54"	48/CASE

Product No.	Weight Wrap	Description	Quantity
S40024	400	24" x 24"	120/CASE
S40036	400	36" x 36"	72/CASE
WR4004040	400	40" x 40"	72/CASE
WR4004545	400	45" x 45"	48/CASE
S40048	400	48" x 48"	48/CASE

Product No.	Weight Wrap	Description	Quantity
WR5001818	500	18" x 18"	288/CASE
WR5002424	500	24" x 24"	120/CASE
WR5003030	500	30" x 30"	120/CASE
WR5003636	500	36" x 36"	72/CASE
WR5004545	500	45" x 45"	48/CASE
WR5004848	500	48" x 48"	48/CASE
WR5005454	500	54" x 54"	48/CASE

Product No.	Weight Wrap	Description	Quantity
S60036	600	36" x 36"	72/CASE
WR6004040	600	40" x 40"	48/CASE
WR6004545	600	45" x 45"	48/CASE
S60048	600	48" x 48"	48/CASE
S60054	600	54" x 54"	48/CASE

*ISO 9073-4 test method compared 100 weight wrap with market leading product showed a 39% improvement in tensile strength

CONTACT YOUR STERIS REPRESENTATIVE TO LEARN MORE

steris.com
Mentor, OH 44060 | USA

MKTG4091 Rev I, 2023
M12700EN.2023-09, Rev.A
©2023 STERIS. All rights reserved.

 **STERIS**