



TRAY MATS

Key Surgical® Tray Mats are designed to be placed under trays/baskets during the inspection and assembly stage of instrument reprocessing. With two different styles to choose from, you can protect sterile barrier packaging from sharp edges of a tray or use for rapid moisture dispersion.

Image represents: MAT001, MAT002, MAT003



ABOUT:

- Choose from two styles depending on your need and sterilization requirements.
- Available in multiple sizes to fit a variety of tray/basket types
- Tray Mats (DRYMAT001, DRYMAT002, DRYMAT003, and DRYMAT018) feature an added absorbent layer designed to provide rapid moisture dispersion during steam sterilization
- Strong, durable material (polypropylene or poly/viscose blend)
- Single use




Image represents: DRYMAT001, DRYMAT002, DRYMAT003, DRYMAT018

PRODUCT NO.	DESCRIPTION	STERILIZATION*	MATERIAL	DIMENSIONS	QUANTITY
DRYMAT018	Tray Mat, 75gsm	Steam	Polypropylene/viscose	305mm x 305mm (12" x 12")	100
DRYMAT003	Tray Mat, 75gsm	Steam	Polypropylene/viscose	305mm x 610mm (12" x 24")	100
DRYMAT002	Tray Mat, 75gsm	Steam	Polypropylene/viscose	500mm x 375mm (19.70" x 15")	100
DRYMAT001	Tray Mat, 75gsm	Steam	Polypropylene/viscose	750mm x 500mm (29.50" x 19.70")	50
MAT003	Tray Mat, 60gsm	Steam, VHP, and gas plasma	Polypropylene	305mm x 610mm (12" x 24")	100
MAT002	Tray Mat, 60gsm	Steam, VHP, and gas plasma	Polypropylene	500mm x 375mm (19.70" x 14.76")	100
MAT001	Tray Mat, 60gsm	Steam, VHP, and gas plasma	Polypropylene	750mm x 500mm (29.50" x 19.70")	50

*See IFU for more information

VISIT KEYSURGICAL.COM OR CALL 800.541.7995 FOR MORE INFORMATION