



# TWIST DRILL RACK

Twist Drill Rack from Key Surgical® is a perfect solution to managing your orthopedic twist drill bits. The rack is designed to hold 5" drill bits, up to 13 at a time, ranging from 1.5 to 5.0mm in diameter. Folds closed for convenient reprocessing but opens up to a tabletop position for easy use in the O.R.

## ABOUT:

- Durable, OR-grade stainless steel
- Holds 13 drill bits ranging from 1.5 - 5.0mm diameter
- Drill bit diameters are marked inside rack for easy identification
- Opens to tabletop stand for the O.R. and folds closed for sterilization, storage, and transportation
- Material compatible with steam sterilization
- Reusable

Rack holds the following size drill bits:  
(sold separately)

DIAMETER	LENGTH	MAX HELD
1.5mm	5"	2
2.0mm	5"	2
2.5mm	5"	2
2.7mm	5"	1
3.2mm	5"	1
3.5mm	5"	1
4.0mm	5"	1
4.5mm	5"	1
4.7mm	5"	1
5.0mm	5"	1



OPEN POSITION AS A TABLE TOP STAND



CLOSED POSITION FOR STORAGE AND TRANSPORTATION

PRODUCT NO.	DESCRIPTION
KI-05-RACK-TD	Twist Drill Rack 6.75" x 4.25" x .75"



a STERIS company