


VASCULAR LOOPS



Vascular Loops are used to support smooth and gentle handling, identifying, retracting, and occluding arteries and veins.

Unlike clamps and retractors, which can pinch vessels flat or even crack them, Vascular Loops feature an elliptical shape. Color-coded for identification purposes.

- Available in 4 colors (red, blue, yellow, or white)
- Available in 3 sizes (mini, maxi, extra maxi)
- Sterile, ready-to-use
- 2 loops per sterile pouch (10 pouches/box)
- Medical-grade silicone
- Not made with natural rubber latex
- Radiopaque



	Product No.	Description	Qty
a.	VL2201	Vascular Tie, red, .075" x 18", mini	10/box
b.	VL2202	Vascular Tie, blue, .075" x 18", mini	10/box
c.	VL2203	Vascular Tie, yellow, .075" x 18", mini	10/box
d.	VL2204	Vascular Tie, white, .075" x 18", mini	10/box
e.	VL2205	Vascular Tie, red, .12" x 18", maxi	10/box
f.	VL2206	Vascular Tie, blue, .12" x 18", maxi	10/box
g.	VL2207	Vascular Tie, yellow, .12" x 18", maxi	10/box
h.	VL2208	Vascular Tie, white, .12" x 18", maxi	10/box
i.	VL2215	Vascular Tie, red, .12" x 32", maxi	10/box
j.	VL2214	Vascular Tie, blue, .20" x 18", extra maxi	10/box
k.	VL2218	Vascular Tie, white, .20" x 22", extra maxi	10/box

CONTACT YOUR STERIS REPRESENTATIVE TO LEARN MORE