

WRITE-ON TAGS



Key Surgical® Write-On Tags offer a simple and colorful answer for for communicating between departments. Write-On Tags allow you to capture specific handling instructions when more detailed information is necessary. Pre-printed tags can help communicate the message clearly and efficiently.

- 3" x 5"
- Durable and economical
- Styles available:
 - pre-printed messages (one-sided or two-sided)
 - blank
- Compatible with steam and VHP (low-temp) sterilization methods. (refer to IFU)
- Not made with natural rubber latex
- Single use
- 100/pkg



WT0009



WT0203



WT0405




WT0304




WT1911




RT100




WT2614



WT2715



WT3217



CONTACT YOUR STERIS REPRESENTATIVE TO LEARN MORE

steris.com
Mentor, OH 44060 | USA

Key Surgical is a registered trademark of STERIS
MKTG4010 Rev D, 2024
M12624EN.2024-04, Rev B
©2024 STERIS. All rights reserved.

