



STOPCON DISPOSABLE ENDOSCOPE VALVE: AIR/WATER CLEANING ADAPTER

STOPCON disposable air/water cleaning adapter for flexible endoscope. Intended to be used during the pre-cleaning process, to help clear the air/water channel of Olympus® GI endoscopes post procedures.

ABOUT:

- Eliminates the traceability challenges commonly associated with reusable valves
- Oil-free valves (and no oil lubrication required, the use of silicone-based oil can create an environment for biofilm formation*)
- Single use
- 100/pkg



EN10920

*Residue from silicone-based lubricant products promotes an environment that protects microorganisms and promotes biofilm formation.

- SGNA, Standards of Infection Prevention in Reprocessing Flexible Gastrointestinal Endoscopes, 2018

