



# FACE MASKS: ASTM LEVEL PPE

Not all personal protective equipment is created equal. This is especially true for nose and mouth protection against the varied environments in which you work. Select the type of mask you need and have peace of mind with ASTM Level 2 and Level 3 face masks from Key Surgical. All masks are not made with natural rubber latex.

## ASTM LEVEL 2 FACE MASKS



**PRODUCT NO.**  
4521

**QUANTITY**  
50/box



**PRODUCT NO.**  
GCPBLSF

**QUANTITY**  
50/box

## ASTM LEVEL 3 FACE MASK



**PRODUCT NO.**  
4522

**QUANTITY**  
50/box

| PRODUCT NUMBER    | ASTM F2100 Level 2 Compliant | ASTM F2100 Level 3 Compliant | Meets Class 1 flammability requirement ( $\geq 3.5$ seconds) per 16 CFR Part 1610 | Bacterial Filtration Efficiency (BFE) $\geq 98\%$ | Particulate Filtration Efficiency (PFE) $\geq 98\%$ | FDA-cleared for use in surgical setting <sup>1</sup> | Physical features of mask   |
|-------------------|------------------------------|------------------------------|---|---|---|--|---|
| 4521 <sup>1</sup> | ●                            |                              | ●   | ●   | ●   | ●  | <ul style="list-style-type: none"> <li>Comfortable, elastic earloops</li> <li>Heat-sealed, 3-layer design with skin-friendly fabric</li> </ul>  |
| GCPBLSF           | ●                            |                              | ●   | ●   | ●   | ●  | <ul style="list-style-type: none"> <li>Comfortable, elastic earloops</li> <li>SecureFit mask technology for a custom fit for any face shape. Reduces venting at the top, bottom, and sides</li> <li>Aluminum strips above the nose and below the chin (extra-long aluminum nose piece)</li> <li>Fluid resistant layers (inner and outer)</li> </ul> |
| 4522              |                              | ●                            | ●   | ●   | ●   | ●  | <ul style="list-style-type: none"> <li>Ties at back of head and neck</li> <li>Ties provide a secure facial fit, reduces venting at the sides</li> <li>Heat-sealed, 3-layer design with skin-friendly fabric</li> </ul>  |



a STERIS company

FOR MORE INFORMATION CALL 800.541.7995 OR VISIT [KEYSURGICAL.COM](https://www.keysurgical.com)    

<sup>1</sup> Applicable to product number 4521, ASTM Level 2 Face Mask. Enforcement Policy for Face Masks and Respirators During the Coronavirus Disease (COVID-19) Public Health Emergency (Revised). Food and Drug Administration. May 2020.