



BRUSH ORGANIZER

Hang your instrument cleaning brushes out to dry and keep them organized with the Key Surgical® Brush Organizer. The wall mount design allows you to maximize usable space in the sterile processing department. The organizer features up to 20 pegs for hanging cleaning brushes.

ABOUT:

- Durable stainless steel
- Flat back for easy flush-mount to wall
- 16" x 35" (1.33ft x 2.91ft)
- 12 permanent pegs (3/4" long)
- Threaded, angled pegs help keep brushes in place
- Two additional pegged bars included (adhesive for universal placement)
- Drip tray catches excess water

PRODUCT NO.	DESCRIPTION
BR-7000	Brush Organizer, 16" x 35"
BR-7100	Brush Organizer with Brush Assortment (BR-4199)

YOU MIGHT ALSO LIKE:

Cleaning Brush Assortments



Note: Brushes shown here for suggested use. Specific brushes shown here not included.



a STERIS company

FOR MORE INFORMATION CALL 800.541.7995 OR VISIT KEYSURGICAL.com

