




ENDOSCOPE MANUAL CLEANING: DOUBLE-HEADER CLEANING BRUSH

Let this combination brush pull a double-header in endoscope reprocessing, not you. With a classic endoscope channel cleaning brush on one end and a valve/control head brush on the other, you have just what you need for removing debris within endoscope channel diameters ranging from 2.0mm – 7.0mm.

FOR CHANNEL
SIZES
2.0MM – 7.0MM



ABOUT:

- 230cm long channel brush
- Double-ended with channel brush on one end, valve/port brush on the other
- Tapered, nylon bristle brush head (5-7mm) with a protective tip at the end to help prevent against potential channel damage
- Valve/port brush attached
 - Full head of soft nylon bristles
 - Barbell handle with grooves for a steady, comfortable grip

PRODUCT NO.	DESCRIPTION	BRISTLE DIAMETER (MM)	VALVE CONTROL BRUSH (MM)	SHEATH DIAMETER (MM)	OVERALL LENGTH (CM)	CHANNEL DIAMETER (MM)	QUANTITY
BX00711609	Double-Header Combination Cleaning Brush	5 - 7 (tapered)	10	1.7	230	2.0 - 7.0	50/box
BX00711629	Double-Header Combination Cleaning Brush	5 - 7 (tapered)	10	1.7	230	2.0 - 7.0	100/box



a STERIS company

VISIT OUR WEBSITE [KEYSURGICAL.COM](https://www.keysurgical.com) OR CALL 800.541.7995 FOR MORE INFORMATION

