



ENDO-LINERS

Endo-Liners are color-coded protective films designed to clearly identify endoscopes that are either clean or contaminated. During temporary storage or transport to and from the procedure room, we've got your scope identification needs covered - literally.

ABOUT:

- Multiple sizes to fit any endoscope reprocessing tray.
- Endo-Liner (103580L and 103580XXL) consists of two protective films that are color-coded (red and green), one thin, opaque plastic tray liner, and one clear plastic self-adhesive tray bag.
- Green side printed with text clearly stating "CLEAN" and red side print with text clearly stating "BIOHAZARD/CONTAMINATED" to comply with applicable OSHA requirements.
- Sterile.
- Single use.



103580L / 103580XXL



PRODUCT NO.	DESCRIPTION	QUANTITY
103580L	Endo-Liner, Large	200/box
103580XXL	Endo-Liner, X-Large	200/box
103580E	Eco Endo-Liner, Large	200/box

YOU MIGHT ALSO LIKE:

Endoscope Cleaning Brushes



a STERIS company

VISIT OUR WEBSITE KEYSURGICAL.COM OR CALL 800.541.7995 FOR MORE INFORMATION

