



# ENDOSCOPE PROTECTION AND TRANSPORT

Successful endoscopy procedures are only possible with clean, working equipment. Managing the transportation of endoscopes at the start or end of a procedure can be challenging. Key Surgical Endoscope Tip Protector and Endo-Liner help with identification and protecting the distal tip from damage.

## ENDOSCOPE TIP PROTECTOR

- Protects the endoscope's distal end during transportation
- Features multiple size openings, designed to fit various scope diameters
- May be used to help remove gross soil during pre-cleaning
- 5.5" L X 2" D
- Porous foam, white
- Single use
- NOTE: if used for cleaning, discard after use and apply a new, unused tip protector for transportation and storage



PRODUCT NO.	QUANTITY
EN10111	50/box



## ENDO-LINER

- Protect reprocessed endoscopes from airborne contamination (particles and microorganisms) during transport and/or storage
- Green side with text stating "CLEAN" and red side with text stating "BIOHAZARD/CONTAMINATED" to comply with OSHA requirements
- Compatible with most cart & tray transportation systems
- Single Use
- Sterile

PRODUCT NO.	DESCRIPTION	QUANTITY
103580L	Large	200/box
103580XXL	Extra large	200/box



a STERIS company