

PERSONAL PROTECTIVE EQUIPMENT



Face Shields

FACE SHIELDS FOR USE DURING MANUAL CLEANING AND DECONTAMINATION

Proper protection from splashes and sprays in the CSSD is a non-negotiable. Our Key Surgical® Face Shields are designed to protect your face during the various tasks at the sink such as brushing and flushing.

- **13" wide, anti-glare, anti-fog, mylar polyester film** for clear visibility without distractions during manual cleaning
- **Vented foam headband and wide elastic strap** for a secure fit while remaining comfortable during long periods of use
- Not made with natural rubber latex
- Single use

Depending on your preference, choose from three (3) styles:

Product No.	Description	Quantity
4509	Full face shield, with splash guard	24/box
4506	Full face shield, 9" length	24/box
4505	Full face shield, 7.5" length	24/box



4509: Full face shield with splash guard for face, chin area, and upper neck. Features pull-down tab for easy expansion of splash guard



4506: Full face shield extended length of 9" to protect the face, chin, and upper neck area



4505: Full face shield length of 7.5" to protect the face

Visit KEYSURGICAL.COM or contact your sales representative for more information