



# ASTM LEVEL 2 FACE MASK

## CDC MEDICAL MASK GUIDELINES

The CDC indicates that EN 14683 Type II, IIR, and ASTM F2100 minimum Level 1 or equivalent are indicated for use for direct clinical care of patients with COVID-19.<sup>1</sup>

### ABOUT:


- ASTM F2100-19 Level 2 compliant
- Approved for use as a surgical mask<sup>2</sup>
- BFE ≥ 98%
- PFE ≥ 98%
- Heat sealed 3-layer design with skin-friendly fabric
- Comfortable, elastic high ear loop (reduces pressure on the ears)
- Meets Class 1 Flammability requirement (≥ 3.5 seconds) per 16 CFR Part 1610
- Not made with natural rubber latex
- Single use

### FLUID RESISTANCE:

- Meets fluid resistance requirement (120 mmHg) equivalent to ISO 22609 Standard Test Method for Fluid Resistance of Medical Face Masks



PRODUCT NO.	DESCRIPTION	QUANTITY
4521	Face Mask, ASTM Level 2	50/box



FOR MORE INFORMATION CALL 800.541.7995 OR VISIT KEYSURGICAL.COM



<sup>1</sup><https://www.cdc.gov/coronavirus/2019-ncov/hcp/non-us-settings/emergency-considerations-ppe.html>

<sup>2</sup> Enforcement Policy for Face Masks and Respirators During the Coronavirus Disease (COVID-19) Public Health Emergency (Revised). Food and Drug Administration. May 2020.