



FULL FACE SHIELD

ABOUT:

- Vented foam headband for airflow and comfort
- 9" long to cover face, chin, and upper neck area
- Protects against splashes and sprays
- Flexible elastic band for easy fit
- Face shield extends around to the side of the face to help provide some peripheral protection



MATERIAL:

4519 Face Shield is composed of the following materials

SHIELD: Polyester-based film

FOAM: Flexible polyurethane film

ADHESIVE: Technomelt 246A MX

ELASTIC BAND: Polyester textured yarn (polyisoprene thread embedded)

PRODUCT NO.	DESCRIPTION	QUANTITY
4519	Full face shield, 9" long, covers chin and upper neck area	24/box



FOR MORE INFORMATION CALL 800.541.7995 OR VISIT KEYSURGICAL.COM

