



# RELIANCE HEMOSTASIS CLIPS

Reliance Hemostasis Clips are used to prevent bleeding in the GI tract, which most often occurs after a polyp has been removed during a GI procedure. The clip is a small metal device that can be rotated and/or repositioned for accurate placement to join the surrounding tissue together.

## ABOUT:

- Designed to clip tissue together during a GI procedure
- Repositionable and rotatable for improved accuracy, workflow, and performance
- Available in 11mm and 16mm opening widths
- Radiopaque
- MR Conditional
- Sterile
- Single-use




EN10400A.16




EN10400A

PRODUCT NO.	DESCRIPTION	QUANTITY
EN10400A.16	Reliance Hemostasis Clip, Coated, 230cm length, 16mm (open width)	5/box
EN10400A	Reliance Hemostasis Clip, Coated, 230cm length, 11mm (open width)	5/box

## YOU MIGHT ALSO LIKE:



Biopsy Forceps



Polypectomy Snares



VISIT OUR WEBSITE [KEYSURGICAL.COM](http://KEYSURGICAL.COM) OR CALL 800.541.7995 FOR MORE INFORMATION

