



IDENTIFICATION TAPE

Key Surgical® Identification Tape is designed to help you identify and organize instruments as they move between sterile processing and various departments. It's easy to apply and withstands repeated sterilization cycles. With the broadest selection of tape available, and the ability to customize your identification tape, you will have a colorful, economical, easy-to-use solution for instrument identification.

“Instrument tape and plastic dipping material, when used properly, are ways of identifying specific instruments.”

-ANSI/AAMI ST79:2017, Section 7.4.1

ABOUT:

- Over 250 colors and patterns available
- Sold by the roll or by the sheet
- Pre-cut sheet tape reduces product waste
- Available widths: 1/8", 1/4", or 3/8"
- Not made with natural rubber latex
- Material compatible with steam sterilization
- CUSTOMIZABLE (Sheet Tape Only)



SHEET TAPE



MORE PATTERNS AND COLORS AVAILABLE

ROLL TAPE



YOU MIGHT ALSO LIKE:

Safely remove identification tape from surgical instrumentation with Key Surgical adhesive remover. The formula is plant-based and citrus-based which provides a pleasant scent during use. Packaged 6 bottles/box, each bottle is 250ml (8.45oz).



BC-200-6



a STERIS company

VISIT OUR WEBSITE KEYSURGICAL.COM OR CALL 800.541.7995 FOR MORE INFORMATION

